

# VLP-RNA Extraction Control

## Product Handling Guide

Shipping:	On Dry or Blue Ice
Catalog number:	MDX068 MDX069
Concentration	10 <sup>4</sup> copies/μL
Batch/Lot No.:	See vial

Store at -20 °C



### Storage and stability:

VLP-RNA Extraction Control is shipped on dry or blue ice. On arrival store at -20 °C for optimum stability. Repeated freeze/thaw cycles should be avoided. Solutions should be mixed/equilibrated after each thawing to avoid phasing.

### Expiry:

When stored under the recommended conditions and handled correctly, full activity of the kit is retained for one year upon arrival.

### Safety precautions:

Read and understand the SDS (Safety Data Sheets) before handling the reagents. Hardcopies of the SDS will be provided with the first shipment, thereafter they will be available upon request.

### Quality control:

Meridian operates under ISO 13485 Quality Management System. VLP-RNA Extraction Control and its components are extensively tested for functionality.

### Notes:

For research and further manufacturing use only

## Description

VLP-RNA Extraction Control contains an internal control RNA sequence, with no known homology to any organism, encapsidated in a virus-like particle. The VLP-RNA Extraction Control is spiked in the sample prior to RNA extraction. Following RNA extraction, VLP-RNA Extraction Control can be detected in RT-qPCR adding the VLP Detection mix to the reaction mix. The detection of VLP-RNA Extraction Control confirms the success of the extraction and amplification steps and reduces the chance of obtaining a false negative result.

## Kit components

Table 1

Component
VLP-RNA Extraction Control
VLP Detection Mix (Red or Orange)

## Users Guidelines:

### Notes:

- Optimal volumes of VLP-RNA Extraction Control may vary depending on sample type and RNA extraction technique. Protocol optimization may be needed.
- After first use, VLP-RNA Extraction control can be stored at 2 to 8 °C up to 3 months.
- Validation of multiplex PCR should be performed prior to high-throughput processes.
- A control reaction should contain all components required for amplification of sample RNA, including Internal Control RNA, to ensure the amplification of the VLP-RNA Extraction Control
- A control to verify the absence of cross-reactivity between the user-assay and Internal Control RNA should be carried out.

## Extraction step

- Brief spin down all tubes before opening.
- Standard Protocol:
  - Spike 4 μL\* of VLP-RNA Extraction Control into each sample.
  - Follow the manufacturer's protocol for total RNA extraction.
  - Elute total RNA in a volume of 50 μL.
- Use 5 μL of the elution volume for a 20 μL RT-qPCR reaction.

\* This volume has to be considered as a reference. The amount of VLP-RNA Extraction Control spiked should be adjusted depending on sample and extraction method used.

## Post-extraction set up master mix preparation

Recommended reagent volumes per 20 μL RT-qPCR mix are given in Table 2

Table 2

Component	Supplied	Volume
2x RT-qPCR master mix	No	10 μL
Target probe/primer mix	No	X μL
Sample RNA from extraction step	No	Up to 5 μL
VLP Detection Mix (Red or Orange)**	Yes	0.8 μL
Reverse transcriptase 100x	No	0.2 μL
RNase inhibitor	No	0.4 μL
Total Volume (for 1 reaction)		20 μL

\* \*Vortex Control Mix before making up the master mix

## Assay setup

The RT-qPCR conditions in Table 4 are suitable for amplicons of up to 200 bp, however they can be varied to suit different commercial RT-qPCR mixes and machine-specific protocols.

Table 3

Step	Temperature	Time	Cycles
Reverse transcription	42 °C	10-20 min	1
Polymerase activation	95 °C	3 min	1
Denaturation	95 °C	10 s	35-40
Annealing/Extension	60 °C	30-45 s	

Acquire VLP-RNA Extraction Control fluorescence signal on the appropriate channel:

- VLP Detection Mix Red (Cy5 - emission wavelength = 670nm)
- VLP Detection Mix Orange (HEX - emission wavelength = 555nm)

Related Products	Cat. No.
VLP-RNA Extraction Control CUS	MDX071

## Technical Support

For any technical enquiries, please contact our Technical Support team via email at: [mbi.tech@meridianlifescience.com](mailto:mbi.tech@meridianlifescience.com)