MERIFLUOR®

GullSORB 5 mL

IgG Inactivation Reagent for the Inactivation of IgG Antibodies in Human Serum

REF XX715

IVD In vitro diagnostic medical device

Rx Only

INTENDED USE

GullSORB is an antihuman IgG reagent used to eliminate IgG interference in human serum containing both IgM and IgG. This reagent has been standardized for use with IgM IFA assays manufactured by Meridian Bioscience, Inc.

SUMMARY AND EXPLANATION OF THE TEST

The detection of immune IgM antibodies in patient serum is a standard serological procedure for the early diagnosis of infectious disease.¹ However, the validity of IgM testing can be compromised by the presence of pathogen-specific (immune) IgG, thereby requiring its removal before testing.²

There are two primary ways in which immune IgG can interfere with assays for IgM antibodies: (a) by competing with specific IgM for substrate binding sites, and (b) by forming immune complexes with rheumatoid factor which can mimic immune IgM.

Competitive Binding

Immune IgG more effectively competes for available binding sites on an antigen substrate than does immune IgM. In some instances, such as samples with low IgM titers, the competition by immune IgG can effectively block IgM binding and cause a false negative test result.

Rheumatoid Factor (RF)

When both immune IgG and RF are present in an IgM-negative sample, interference presents a different testing problem. RF is usually an IgM class, anti-IgG antibody that can bind with the immune IgG in the sample to form an immune complex. The complex becomes a reactive site for the μ -chain specific conjugate, producing an IgM false positive result.³

BIOLOGICAL PRINCIPLES

The GullSORB Reagent provides a means of eliminating both types of interference caused by immune IgG through an immunoprecipitation reaction with antihuman IgG antibodies. Patient serum is mixed with GullSORB prior to testing. Patient sera treated with GullSORB Reagent can then be tested for IgM with confidence in the accuracy of the assay.

REAGENTS/MATERIAL PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

GullSORB Reagent, 1 dropper vial, 5 mL, Anti-human IgG (caprine), containing 0.1% (w/v) sodium azide. Ready for use and stable until the date on the vial label when stored at 2-8 C.

MATERIALS NOT PROVIDED

- Small volume capacity tube to accommodate 50 μ L with mixing (500 μ L conical micro centrifuge tubes are recommended).
- Microliter pipettors, 5 μ L or equivalent.

PRECAUTIONS

All reagents are for in vitro diagnostic use only.

WARNINGS

- Contains SODIUM AZIDE, which may form explosive azide salts if allowed to build up in lead or copper plumbing. When disposing of the reagent in sinks, be sure to flush with a large volume of water.
- Sodium azide is very toxic in case of ingestion. In contact with acid it emits a very toxic gas. In case of contact with the skin wash immediately and abundantly with water.

HAZARD and PRECAUTIONARY STATEMENTS

Refer to the SDS, available at www.meridianbioscience.com for Hazards and Precautionary Statements.

SHELF LIFE AND STORAGE

Store the kit at 2-8 C.

SPECIMEN COLLECTION AND PREPARATION

Blood obtained by venipuncture should be allowed to clot at room temperature. Serum should be separated and stored at 2-8 C or frozen at -20 C if not processed within one week. Sera exhibiting hemolysis, lipemia or microbial growth may give erroneous results.

TEST PROCEDURE

- Holding the dropper vial in a strictly vertical position, dispense one drop of GullSORB Reagent into a suitable tube.
- Add 5 μ L of patient serum and mix gently. The mixture is a 1:10 dilution of the patient's serum. The mixture may be used immediately – no incubation is required.
- Use an appropriate volume of the mixture according to the instructions for the IgM assay in use. The flocculent material in the mixture will not affect Meridian IFA tests.

LIMITATIONS OF THE PROCEDURE

When used according to instructions, GullSORB Reagent will neutralize up to 15 mg/mL of IgG antibody in human serum. Sera with IgG levels greater than 15 mg/mL should be treated with a proportional increase in amount of GullSORB Reagent.

SPECIFIC PERFORMANCE CHARACTERISTICS

The table shows the results of a comparison between GullSORB Reagent and ion exchange chromatography for efficiency of IgG removal. Subclass-specific radial immunodiffusion confirmed the removal of all IgG subclasses, including IgG4, from serum immunoprecipitated with GullSORB Reagent and centrifuged.

GullSORB Reagent showed no reactivity with antigens of Epstein-Barr Virus, Herpes Simplex Virus, Cytomegalovirus, Varicella-Zoster Virus, or *Chlamydia trachomatis*.

Total IgG Concentration by Nephelometry			
Sample	Untreated Serum (mg/mL)	After Ion-Exchange Chromatography (mg/mL)	After GullSORB Reagent (mg/mL)
24	08,80	0,39	0,00
25	19,00	0,67	0,38
26	15,00	0,20	0,18
27	11,00	0,39	0,00
28	11,30	0,70	0,38
29	11,20	0,87	0,10
30	09,20	0,60	0,08
31	13,20	0,40	0,15
32	13,80	0,76	0,55
Mean	12,50	0,55	0,16
1SD	03,17	0,22	0,14
	Efficiency	95,6% ±1,8%	98,7% ±1,1%

ITALIANO

MERIFLUOR®

GullSORB 5 mL

IgG Inactivation Reagent for the Inactivation of IgG Antibodies in Human Serum

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REFERENCES

1. Gallo D, et al. Multiple-Antigen slide test for detection of immunoglobulin m antibodies in newborn and infant sera by immunofluorescence. J Clin Microbiol 1981 13:631-36.
2. Doerr HW, Rentschler M and Scheiffer G. Serologic detection of active infections with human herpes viruses (CMV, EBV, HSV, VZV): diagnostic potential of IgA class and IgG subclass – specific antibodies. infection 1987 14:93-98.
3. Smith SB, Cochran SA and Repetti CF. Evaluation of a fluorescence immunoassay for toxoplasma IgM antibodies. Abstract from the 87th Annual Meeting of the American Society for Microbiology; 1987; Atlanta, Georgia. p. 417.

SN11758

REV. 04/16

Manufactured By

Meridian Bioscience, Inc.
USA/Corporate Office
3471 River Hills Drive
Cincinnati, Ohio 45244
Telephone: 513.271.3700
Orders/Customer Service:
800.543.1980
Technical Support Center:
800.343.3858
Information Fax: 513.272.5432
Ordering Fax: 513.271.0124

Meridian Bioscience Europe S. r. l.
Via dell'Industria, 7
20020 Villa Cortese, Milano
ITALY
Tel: +39 0331 43 36 36
Fax: +39 0331 43 36 16
Email: info@meridianbioscience.eu
WEB: www.meridianbioscience.eu

EC REP
Authorized Representative

Meridian Bioscience Europe s.a./n.v.
2 Avenue du Japon - 1420 Braine l'Alleud
BELGIUM
Tel: +32 (0) 67 89 59 59
Fax: +32 (0) 67 89 59 58
Email: info.bn@meridianbioscience.eu


Meridian Bioscience Europe France
34 rue de Ponthieu - 75008 Paris
FRANCE
Tel: +33 (0) 1 42 56 04 40
Fax: +33 (0) 9 70 06 62 10
Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.
Postbus 301 - 5460 AH Veghel
NETHERLANDS
Tel: +31 (0) 411 62 11 66
Fax: +31 (0) 411 62 48 41
Email: info.bn@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simbolos, Zeichenerklärung).

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
LOT	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
IVD	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	EC REP	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostika 98/79/EG.	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdünnungspuffer befindet
REF	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	SOLN STOP	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung	CONJ ENZ	Enzyme Conjugate / Coniugato enzimatico / Conjugué enzymatique / Conjugado enzimático / Enzymkonjugat
SN	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
TEST	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät	REAG	Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	BUF WASH	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
BUF	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
CONJ	Conjugate / Coniugato / Conjugué / Conjugado / Konjugat	DIL SPE	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer
SUBS	Substrate / Substrato / Substrat / Substrato / Substrat	BUF WASH 20X	Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
Rx Only	Caution: US Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner / La legge Federale negli Stati Uniti prevede che questo test sia venduto solo a un operatore abilitato ad operare nella salute / Attention: La vente ou la prescription de ce dispositif par un médecin agréé est soumise à des restrictions par la loi fédérale des US / Cuidado: leyes federales de los EU limitan la venta de dispositivos por o por la orden de licenciados en la práctica del cuidado de la salud / Vorsicht: Der Verkauf oder die Verschreibung dieses Geräts durch einen zugelassenen Arzt unterliegt den Beschränkungen des US-Bundesgesetzes.	DET REAG	Detection Reagent / Reagente Diretto / Réactif de Detection / Reactivo de Detección / Nachweise Reagenz

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM. USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.