
## Urea Breath Test For Detection Of H. pylori




1. Instruct patient to hold their breath for 4-5 seconds, and then exhale into the blue BASELINE breath sample bag



2. Push the cap until it clicks



3. Prepare the test drink by dissolving the citrica powder and the 13C-enriched urea tablet in 150 to 200 ml of tap water in the drinking cup provided



- 4. A. Close the cup firmly using both hands
  - B. Mix until it is completely dissolved



 Open the lid and instruct patient to drink the entire solution using the straw within 2 minutes



6. Wait 15 minutes



7. Instruct patient to hold their breath for 4-5 seconds, and then exhale into the gray POST INGESTION breath sample bag



8 Push the cap until it clicks



Place both breath sample bags in the provided transport bag

For detailed information regarding the step-by-step procedure, refer to the IDkit Hp® Two Package Insert





