

# Taq DNA Polymerase Product Handling Guide

Shipping:	On Dry/Blue Ice
Catalog numbers:	MDX001
Batch No.:	See vial
Concentration:	5 U/ $\mu$ L

Store at  $-20^{\circ}\text{C}$



## Storage and stability:

Taq DNA Polymerase is shipped on dry/blue ice. On arrival store at  $-20^{\circ}\text{C}$  for optimum stability. Repeated freeze/thaw cycles should be avoided. Thawing during transportation does not affect the product performance. Solutions should be mixed/equilibrated after each thawing to avoid phasing.

## Expiry:

When stored under the recommended conditions and handled correctly, full activity of the kit is retained until the expiry date on the outer box label.

## Safety precautions:

Read and understand the SDS (Safety Data Sheets) before handling the reagents. Hardcopies of the SDSs will be provided with the first shipment, thereafter they will be available upon request.

## Quality control:

Bioline operates under ISO 13485 Management System. The Taq DNA Polymerase and its components are extensively tested for activity, processivity, efficiency, heat activation, sensitivity, absence of nuclease contamination and absence of nucleic acid contamination.

## Notes:

This reagent has been manufactured under 13485 Quality Management System and is suitable for further manufacturing use as an IVD component.

## Description

Taq DNA Polymerase is a high performance PCR product containing a separate 5x Taq Reaction Buffer containing dNTPs,  $\text{MgCl}_2$  and enhancers at optimal concentrations, removing the need for optimization and providing superior fast amplification across a range of templates.

## Kit components

Table 1

Component
Taq DNA Polymerase
Taq Reaction Buffer, 5x

## Users Guidelines

The optimal conditions will vary from reaction to reaction and are dependent on the template/primers used.

The Taq Reaction Buffer, 5x comprises of 5 mM dNTPs, 15 mM  $\text{MgCl}_2$ , stabilizers and enhancers. The concentration of each component has been extensively optimized, additional PCR enhancers are not recommended.

Forward and reverse primers are generally used at the final concentration of 0.2-0.6  $\mu\text{M}$  each. As a starting point, we recommend using 0.4  $\mu\text{M}$  final concentration (*i.e.* 4 pmol of each primer per 20  $\mu\text{L}$  reaction volume).

For DNA templates with low structural complexity, such as plasmid DNA, we recommend using 50 pg - 10 ng DNA per 50  $\mu\text{L}$  reaction volume. For eukaryotic genomic DNA, we recommend a starting amount of 200 ng DNA per 50  $\mu\text{L}$  reaction, this can be varied between 5 ng - 500 ng.

## PCR reaction setup

Prepare a master mix of Taq DNA Polymerase and assay-specific primers (see recommended composition in Table 2).

Table 2

Reagent	Volume	Final Concentration
Taq Reaction Buffer, 5x	4 $\mu\text{L}$	1x
Template	As required	As required
20 $\mu\text{M}$ Forward Primer	0.4 $\mu\text{L}$	400 nM
20 $\mu\text{M}$ Reverse Primer	0.4 $\mu\text{L}$	400 nM
Taq DNA Polymerase	0.4 $\mu\text{L}$	0.1 U/ $\mu\text{L}$
Water (ddH <sub>2</sub> O)	$\leq 20 \mu\text{L}$	

## PCR amplification

The PCR conditions in Table 3 are suitable for amplicons of up to 1 kb.

Table 3

Step	Temperature	Time	Cycles
Initial denaturation	95 $^{\circ}\text{C}$	1 min	1
Denaturation	95 $^{\circ}\text{C}$	15 s	25-35
Annealing	User determined	15 s	
Extension	72 $^{\circ}\text{C}$	10 s	

## Technical Support

For any technical enquiries, please contact our Technical Support team via email at: [mbi.tech@meridianlifescience.com](mailto:mbi.tech@meridianlifescience.com)

Bioline Reagents Ltd  
UNITED KINGDOM

Tel: +44 (0)20 8830 5300  
Fax: +44 (0)20 8452 2822

Meridian Life Science Inc.  
USA

Tel: +1 901 382 8716  
Fax: +1 901 382 0027

Bioline GmbH  
GERMANY

Tel: +49 (0)337 168 1229  
Fax: +49 (0)3371 68 1244

Bioline (Aust) Pty. Ltd  
AUSTRALIA

Tel: +61 (0)2 9209 4180  
Fax: +61 (0)2 9209 4763

Bioline France  
FRANCE

Tel: +33 (0)1 42 56 04 40  
Fax: +33 (0)9 70 06 62 10