

# High Conc. Glycerol-Free Bst Product Handling Guide

Shipping:	On dry/blue ice
Catalog number:	MDX018
Batch No.:	See vial
Concentration:	100 units/ $\mu$ L

Store at  $-20\text{ }^{\circ}\text{C}$



## Storage and stability:

High Conc. Glycerol-Free Bst is shipped on dry or blue ice. On arrival store at  $-20\text{ }^{\circ}\text{C}$  for optimum stability. Repeated freeze/thaw cycles should be avoided. Solutions should be mixed/equilibrated after each thawing to avoid phasing.

## Expiry:

When stored under the recommended conditions and handled correctly, full activity of the kit is retained until the expiry date on the outer box label.

## Safety precautions:

Read and understand the SDS (Safety Data Sheets) before handling the reagents. Hardcopies of the SDS will be provided with the first shipment, thereafter they will be available upon request.

## Quality control:

Meridian operates under ISO 13485 Quality Management System. High Conc. Glycerol-Free Bst and its components are extensively tested for activity, processivity, efficiency, heat activation, sensitivity, absence of nuclease contamination and absence of nucleic acid contamination.

## Notes:

For research or further manufactured use only.

## Description

High Conc. Glycerol-Free Bst is a DNA polymerase (exonuclease minus), with strand-displacement properties. High Conc. Glycerol-Free Bst is used for Isothermal DNA amplification and LAMP (Loop-mediated Isothermal Amplification).

## Kit components

Table 1

Component
High Conc. Glycerol-Free Bst (100 U/ $\mu$ L)
Bst Reaction Buffer, 10x
Enzyme Dilution Buffer, 10x

## Users Guidelines

Thawing during transportation does not affect the product performance. Prior to use or storing at  $-20\text{ }^{\circ}\text{C}$ , the thawed reagents must be thoroughly mixed by 10 inversions.

## LAMP Optimisation

It is recommended that  $\text{MgSO}_4$  (not supplied) is supplemented to a total concentration of 8 mM as indicated in the example LAMP protocol, but this can be optimised according to individual assay requirements.

dNTPs can be optimised between 1.6 mM (0.4 mM each) and 6 mM (1.5 mM each) final concentration.

High Conc. Glycerol-Free Bst can be diluted using the Enzyme Dilution Buffer provided to between 0.12 and 0.32 U/ $\mu$ L final concentration, depending on individual assay requirements.

## Working Concentration of Bst

High Conc. Glycerol-Free Bst is diluted in 1x Enzyme Dilution Buffer, which is glycerol-free, to a working concentration of 8 U/ $\mu$ L as described in Table 2.

Table 2

Reagent	Ratio	Volume for 48 + 2 Reactions
Bst DNA Polymerase, 100 U/ $\mu$ L	0.08	4 $\mu$ L
Bst Reaction Buffer, 10x	0.1	5 $\mu$ L
Water (sterile, distilled)	0.82	41 $\mu$ L

## Typical LAMP reaction conditions:

Incubate at  $60\text{ }^{\circ}\text{C}$  for 60 minutes.

The following protocol is for a standard 25  $\mu$ L LAMP reaction to be used as a starting point for optimization.

Reagent	Volume	Final Concentration
Bst Reaction Buffer, 10x	2.5 $\mu$ L	1x
dNTP Mix (100 mM - 25 mM each)	0.4 - 1.5 $\mu$ L	1.6 - 6 mM
$\text{MgSO}_4$ (100 mM)	1.5 $\mu$ L	6 mM (8 mM total)
FIP/BIP Primers (25X)	1 $\mu$ L	1.6 $\mu$ M
F3/B3 Primers (25X)	1 $\mu$ L	0.2 $\mu$ M
Loop F/B Primers (25X)	1 $\mu$ L	0.4 $\mu$ M
High Conc. Glycerol-Free Bst (3-8 U/ $\mu$ L)	1 $\mu$ L	0.12 - 0.32 U/ $\mu$ L
Sample DNA	variable	> 10 copies
Water (ddH <sub>2</sub> O)	to 25 $\mu$ L	

Related Products	Cat. No.
Bst DNA Polymerase, 8 U/ $\mu$ L	MDX012
Bst Reaction Buffer, 10x	MDX076
Enzyme Dilution Buffer, 10x	MDX080
Enzyme Dilution Buffer, 1x	MDX078
dNTP Mix, 100mM	MDX051

## Technical Support

For any technical enquiries, please contact our Technical Support team via email at: [mbi.tech@meridianlifescience.com](mailto:mbi.tech@meridianlifescience.com)