

Fast 1-Step RT-qPCR Mix

Product Handling Guide

Shipping:	On Dry/Blue Ice
Catalog number:	MDX032
Batch No.:	See vial
Concentration:	2x

Store at -20 °C



Storage and stability:

Fast 1-Step RT-qPCR Mix is shipped on dry/blue ice. On arrival store at -20 °C for optimum stability. Repeated freeze/thaw cycles should be avoided. Thawing during transportation does not affect the product performance. Solutions should be mixed/equilibrated after each thawing to avoid phasing.

Expiry:

When stored under the recommended conditions and handled correctly, full activity of the kit is retained until the expiry date on the outer box label.

Safety precautions:

Read and understand the SDS (Safety Data Sheets) before handling the reagents. Hardcopies of the SDS will be provided with the first shipment, thereafter they will be available upon request.

Quality control:

Bioline operates under ISO 13485 Quality Management System. Fast 1-Step RT-qPCR Mix and its components are extensively tested for activity, processivity, efficiency, heat activation sensitivity, absence of nuclease contamination and absence of nucleic acid contamination.

Notes:

This reagent has been manufactured under 13485 Quality Management System and is suitable for further manufacturing use as an IVD component.

Description

Fast 1-Step RT-qPCR Mix is a combination of the latest advances in buffer chemistry and PCR enhancers and stabilizers, together with an antibody-mediated hot-start polymerase, dNTPs and MgCl₂ and separate reverse transcriptase and RNase inhibitor. Fast 1-Step RT-qPCR Mix has been designed for highly reproducible, accurate assay results under fast thermal cycling conditions, delivering excellent results in multiplex assays, making it ideal for, automated, high-throughput systems.

Kit components

Table 1

Component
Fast 1-Step RT-qPCR Reaction Mix, 2x
MMLV-RT
RNase inhibitor

Users Guidelines

Master mix preparation

Recommended reagent volumes per 20 µL qPCR mix are given in Table 2.

Table 2

Reagent	Volume
Fast 1-Step RT-qPCR Reaction Mix, 2x	10 µL
Primer-Probe Mix, 20x	1 µL
Template*	Up to 9 µL
MMLV-RT	0.2 µL
RNase inhibitor	0.4 µL
Water	As required
Total volume	Up to 20 µL

*Use 1 pg to 1 µg of total RNA or 0.01 pg to 20 µg purified mRNA in a single PCR

Assay setup

The qPCR conditions in Table 3 are suitable for amplicons of up to 200 bp. These cycling parameters have been optimized for Fast 1-Step RT-qPCR Mix on a number of platforms, however they can be varied to suit different machine-specific protocols.

Table 3

Step	Temperature	Time	Cycles
Reverse transcription**	45 °C	10 min	1
Polymerase activation	95 °C	2 min	1
Denaturation	95 °C	5 s	45
Annealing/Extension**	60 °C	20 s	

**When multiplexing, the reverse transcription reaction time can be extended up to 20 minutes and/or the temperature can be increased up to 48° C and the annealing/extension time can be extended up to 60 seconds and/or the temperature can be increased up to 65° C.

Technical Support

For any technical enquiries, please contact our Technical Support team via email at: mbi.tech@meridianlifescience.com

Bioline Reagents Ltd
UNITED KINGDOM

Tel: +44 (0)20 8830 5300
Fax: +44 (0)20 8452 2822

Bioline USA Inc.
USA

Tel: +1 901 382 8716
Fax: +1 901 382 0027

Bioline GmbH
GERMANY

Tel: +49 (0)337 160222 00
Fax: +49 (0)3371 60222 01

Bioline (Aust) Pty. Ltd
AUSTRALIA

Tel: +61 (0)2 9209 4180
Fax: +61 (0)2 9209 4763

Bioline France
FRANCE

Tel: +33 (0)1 42 56 04 40
Fax: +33 (0)9 70 06 62 10