

# Lyo-Ready 1-Step RT-qPCR ROX Mix

## Product Handling Guide

Shipping:	Blue Ice
Catalog number:	MDX174
Batch No.:	See vial
Concentration:	4x

Store at **-20 °C**



### Storage and stability:

Lyo-Ready 1-Step RT-qPCR ROX Mix is shipped on blue ice. On arrival store at -20 °C for optimum stability. Repeated freeze/thaw cycles should be avoided. Thawing during transportation does not affect the product performance. Solutions should be mixed/equilibrated after each thawing to avoid phasing.

### Safety precautions:

Read and understand the SDS (Safety Data Sheets) before handling the reagents. The SDS (Safety Data Sheets) are available upon request.

### Quality control:

Meridian operates under ISO 13485 Quality Management System. Lyo-Ready 1-Step RT-qPCR ROX Mix is tested functionally before its release (see Test Release document)

### Notes:

For research and further manufacturing use only.

## Description

Lyo-Ready 1-Step RT-qPCR ROX Mix is a glycerol-free, one tube formulation combining the latest advances in buffer chemistry and PCR enhancers, together with an optimized concentration of reverse transcriptase, antibody-mediated hot-start polymerase, dNTPs and MgCl<sub>2</sub>. Lyo-Ready 1-Step RT-qPCR ROX Mix has been designed for highly reproducible, accurate RNA and DNA target amplification, delivering excellent results in multiplex assays. In order to produce room temperature lyophilized RT-qPCR reagents, assay specific primers and probes can be added to Lyo-Ready 1-Step RT-qPCR ROX Mix for subsequent lyophilization.

## Kit components

Table 1

Component
Lyo-Ready 1-Step RT-qPCR ROX Mix, 4x

## Users Guidelines

### Master mix preparation

Recommended reagent volumes per 20 µL RT-qPCR mix are given in Table 2.

Table 2

Reagent	Volume
Lyo-Ready 1-Step RT-qPCR ROX Mix, 4x	5 µL
Primer-Probe Mix, 20x	1 µL*
Water	x µL
Total volume	Up to 20 µL

\* Primer and probe concentration needs be optimized

Dispense into reaction vials and immediately transfer into a freeze-dryer and run a suitable drying cycle.

For long-term storage at ambient temperatures, the lyophilized product from the freeze-dryer should be packaged with a silica sachet in a heat sealed pouch at low relative humidity conditions.

## Assay setup

Rehydrate the lyophilized RT-qPCR master mix in the reaction vials with 20 µL template-containing solution and run RT-qPCR.

The RT-qPCR conditions in Table 3 are suitable for amplicons of up to 200 bp. These cycling parameters have been optimized for Lyo-Ready 1-Step RT-qPCR ROX Mix on a number of platforms, however they can be varied to suit different machine-specific protocols.

Table 3

Step	Temperature	Time	Cycles
Reverse transcription	50 °C	10 mins	1
Polymerase activation	95 °C	2 min	1
Denaturation	95 °C	10 s	45
Annealing/Extension*	60 °C	30 s	

\* When multiplexing, the annealing/extension time can be extended up to 60 seconds and/or the annealing/extension temperature can be increased up to 65°C.

**Optional ROX:** Lyo-Ready 1-Step RT-qPCR ROX Mix is premixed with ROX (5-carboxy-X-rhodamine, succinimidyl ester), so that where necessary, ROX fluorescence can be optionally detected on certain real-time instruments. If your real-time instrument has the capability of using ROX and you wish to use this option, then this option must be selected by the user in the software

## Technical Support

For any technical enquiries, please contact our Technical Support team via email at: [mbi.tech@meridianlifescience.com](mailto:mbi.tech@meridianlifescience.com)