**INTENDED USE**  
The Revogene® *C. difficile* assay performed on the Revogene instrument is a qualitative *in vitro* diagnostic test that utilizes automated sample processing and real-time polymerase chain reaction (PCR) to detect the toxin B (*tcdB*) gene of toxigenic *Clostridium difficile* (*C. difficile*) in unformed (liquid or soft) stool specimens obtained from patients suspected of having *C. difficile* infection (CDI). The Revogene *C. difficile* assay is intended to be used by trained laboratory personnel in clinical, hospital and reference laboratory settings.

**SUMMARY AND EXPLANATION**  
*C. difficile* is a gram-positive bacillus, strict anaerobic spore forming, which is part of the normal intestinal flora in 1-3% of adults and up to 45% of infants<sup>1,2</sup>. In patients taking antibiotics in high doses or for long periods, the normal intestinal flora (bacteria) can be destroyed, allowing *C. difficile* to overgrow. In such cases, *C. difficile* produces toxins that can damage the intestines and cause mild to severe diarrhea and potentially fatal intestinal disorders such as pseudomembranous colitis (inflammation of the large intestine), toxic megacolon and sepsis<sup>3,4,5,6</sup>.

*C. difficile* infection (CDI) is the leading cause of infectious diarrhea in hospitals and long-term care facilities in industrialized countries.  
*C. difficile* is shed in feces. People can become infected if they touch items or surfaces that are contaminated with feces and then touch their mouth or nose. Healthcare workers can spread the bacteria to other patients or contaminate surfaces through hand contact.

The epidemiology of CDI has changed over the last decade<sup>7,8</sup>. Epidemic strains (BI / NAP1 / 027) of virulent *C. difficile* cause serious diseases that often require colectomy resulting in an increased mortality<sup>4,5</sup>. These strains are more likely to spread in hospitals due of antimicrobial resistance and spore formation.


The guidelines of the Centers for Disease Control and Prevention (CDC) recommend the use of a *C. difficile* diagnostic for patients presenting clinical (e.g. unformed stool) and epidemiological risk factors since the rapid identification of infections with *C. difficile* is of major importance on the management of the disease<sup>9</sup>.

Rapid, sensitive and specific diagnostic tests are necessary tools for the early detection of the disease and the quick implementation of effective infection control measures in the fight against CDI. Empiric treatment without a precise CDI diagnostic is inappropriate as even in an epidemic environment, only about 30% of hospitalized patients will acquire a CDI<sup>10</sup>. The efficiency and effectiveness of CDI diagnosis remains a challenge for clinicians and microbiologists.

The Revogene *C. difficile* assay can provide results from up to eight specimens in approximately 70 minutes. The assay minimizes operator intervention from the time the single-use disposable microfluidic cartridge (named PIE hereinafter) containing the sample is placed into the Revogene carousel until results are available.

**PRINCIPLE OF THE PROCEDURE**  
The Revogene automates sample homogenization, sample dilution, cell lysis, DNA amplification and detection of the amplified PCR products. User intervention is only required for discharging the patient specimen into the Sample Buffer Tube (SBT), transferring the sample from SBT into the PIE, and loading/unloading the PIEs into the Revogene carousel.

Each PIE is a completely integrated closed device in which a sample is dispensed and processed through different microfluidic chambers and channels which allow for the sample processing (i.e. sample homogenization, sample dilution, and cell lysis) and subsequent real-time PCR steps. The liquid from a single sample is transferred by centrifugation from one chamber to the next in sequence and all reagents specific for the PCR reaction are incorporated and dried within the PCR well(s) (Figure 1).

**Figure 1. Top View of a PIE.**  
1: Sample Loading chamber, 2: Homogenization chamber, 3: Overflow chamber, 4: Dilution/Lysis chamber, 5: Three PCR wells (#1 to #3 at the left) and one waste chamber (at the right).

A Process Control (PrC) is incorporated into each PIE to verify sample processing and amplification steps. The PrC allows for the verification of potential inhibitor substances as well as microfluidic, instrument or reagent failure. The amplified products are detected in real time using target-specific TaqMan® chemistry-based probes. No operator intervention is necessary once a PIE is loaded into the Revogene.

The Revogene can process from one up to eight samples simultaneously in the same run. The carousel must contain eight PIEs to maintain thermodynamic balance within the run. At run completion, the results are computed by the system from measured fluorescent signals and embedded calculation algorithms. Results that are displayed on the touchscreen may be printed, transferred and/or stored by the user using the USB port or the connectivity option.

**REAGENTS AND MATERIALS**

The *C. difficile* kit contains sufficient reagents and materials to process 24 specimens. The kit contains the following:

1. 24 Disposable Transfer Loop (DTL) on the side compartment: DTL consists of a single-use 5 microliters (µL) loop for transferring the unformed (liquid or soft) stool specimen to the SBT.
2. 24 individual pouches and each pouch contains the following materials:
  - a. 1 Disposable Transfer Tool (DTT): DTT consists of a single-use transfer pipette for transferring the sample from the SBT to the PIE.
  - b. 1 Sample Buffer Tube (SBT): Barcode-labeled tube containing TE 1X buffered solution (Tris-HCl pH 8.0/EDTA, Na<sub>2</sub>) as a dilution and preservation buffer for sample.
  - c. 1 *C. difficile* PIE: Barcode-labeled integrated device, composed of dried reagents allowing sample process and real-time PCR steps for PrC DNA and *C. difficile* toxin B gene DNA simultaneous amplification/detection. Each PIE contains PrC, PrC-specific primers and Taqman® chemistry-based probe, *C. difficile tcdB* gene-specific primers and Taqman® chemistry-based probe, dNTPs, buffer and DNA polymerase.

**MATERIALS/EQUIPMENTS REQUIRED BUT NOT PROVIDED**

1. Revogene® (cat# 610210)
2. Disposable gloves; powderless
3. Vortex mixer with a maximal speed of at least 3 200 rpm
4. Sample rack (cat# 132539; optional)
5. MOCK PIE (cat# 610210; optional)

## WARNING AND PRECAUTIONS

1. This product can only be used on the Revogene.
2. Do not use the kit if the label that seals the outer box is broken upon arrival.
3. Do not use *C. difficile* PIEs if the protective pouches are open or broken upon arrival.
4. Do not interchange DTT, SBT, and PIE between kit lots.
5. Each single-use DTL, DTT and *C. difficile* PIE are used to process one sample. Do not reuse DTL, DTT or PIE.
6. Always handle specimens as if they are infectious and in accordance with Good Laboratory Practices such as those described in Biosafety in Microbiological and Biomedical Laboratories<sup>11</sup> and in CLSI Document M29-A4<sup>12</sup>.
7. Wear disposable powderless gloves while handling specimens and thoroughly wash hands afterwards.
8. The *C. difficile* PIE contains dried reagents. The protective pouch should not be opened until ready to perform the test.
9. Dispose of unused materials and reagents and waste, including used PIEs, in accordance with country, federal, provincial, state and local regulations.
10. Do not open or break apart the PIE after use to avoid contamination with amplification products and/or infectious particles.
11. Do not use a PIE that has been dropped, shaken or inverted after the sample has been loaded as this may cause invalid results.
12. The *C. difficile* assay does not provide susceptibility results. Additional time is required to culture and perform susceptibility testing.
13. Do not use a kit that has passed its stated expiration date.
14. Do not refrigerate the loaded PIE.
15. An amount of stool exceeding the recommended amount may inhibit the *C. difficile* assay.
16. Each run must be performed with eight PIEs in the Revogene carousel to maintain thermodynamic and mechanical balance within the run.

## STORAGE AND STABILITY

1. Collected specimens should be stored between 2 C and 25 C during transport.
2. Stool specimens can be stored at 25 C for up to 2 days, or at 2-8 C for up to 4 days.
3. Inoculated SBT can be stored at 25 C for up to 2 days, or at 2-8 C for up to 3 days.
4. Store the *C. difficile* kit at 2-25 C. The expiration date is indicated on the box kit's label.
5. Do not open a pouch until ready to perform testing. Use the PIE within 1 hour after opening the pouch.

## INSTRUCTION FOR USE

### SPECIMEN COLLECTION AND TRANSPORT

Specimen type: Unformed (liquid or soft) stool specimens obtained from patients suspected of having CDI.

Collect the unformed stool in a dry, clean container according to the local guidelines or procedures.

1. Transfer liquid or soft stool (but not urine) into the container. Avoid mixing toilet paper, water or soap with the specimen.
2. Label the container with the specimen or patient identification (ID) and send to the laboratory for testing (refer to **Storage and Stability** section).

### SAMPLE PREPARATION AND HANDLING

**NOTE 1:** Start the test within 1 hour after opening the pouch containing the PIE.

**NOTE 2:** The content of one pouch is required for each specimen to be tested.


### SBT PREPARATION

1. For each specimen to be tested, unseal the right side of the pouch (when facing label) containing DTT and remove the SBT from the pouch.
2. Identify (or label) the SBT with the appropriate specimen identification without obscuring or writing over the barcodes. Place the SBT on the Sample rack, if used.
3. Vortex the specimen at maximal speed for 15 seconds. Dip the provided DTL into the stool. Remove any excess of soft stool present on the outside of the loop in order to take approximately 5 µL.
4. Remove the cap from the SBT, shake the DTL into the SBT for 2-3 seconds and/or swirl to dislodge sample from the loop. Make sure that only one SBT is open at once.
5. Replace the cap on the SBT, tightly close the SBT cap, and place it on the Sample rack, if used.
6. Prepare any additional specimens for testing by repeating steps 1 to 5 then proceed to step 7.
7. When all samples are prepared, proceed to *C. difficile* PIE preparation (next section).

### PIE PREPARATION

**NOTE 1:** Processing one sample at a time:

8. Vortex the SBT for 15 seconds at maximal speed.
9. Unseal the left side of the pouch (when facing label) containing the PIE, removing it from the pouch.
10. Place the *C. difficile* PIE on the Sample rack, if used.
11. Using the DTT, aspirate the sample buffer by squeezing the entire bulb. The liquid level into the DTT must be anywhere between the two marks (**Figure 2**). If the liquid level is not between the two marks, discharge the SB volume completely in the SBT by squeezing the entire bulb and repeat step 11.
12. Discharge completely the SB into the sample loading chamber of the PIE (**Figure 1**).
13. Close the cap of the PIE tightly, making sure the cap lock is well in place. Do not refrigerate the loaded PIE. Make sure that only one PIE is open at once.
14. Repeat steps 8 to 13 for any additional samples then proceed to step 15.
15. The stool specimens and the inoculated SBT can be stored at 2-8 C or at 25 C within the timeframe defined in the **Storage and Stability** section.

**Figure 2.**  
Representation of an appropriate sample buffer level using the Disposable Transfer Tool (DTT).

### REVOGENE OPERATION

**NOTE 1:** A maximum of eight samples can be processed simultaneously in a single run using the Revogene (including external controls).

**NOTE 2:** Each run must be performed with eight PIEs in the Revogene. When less than eight samples are processed, the empty places must be filled with MOCK PIEs\*.

**NOTE 3:** Refer to the Revogene Operator's Manual<sup>13</sup> for further information regarding Revogene set-up and operation.

1. Power on the Revogene (if not already done). The software will launch automatically.
2. Log in by entering the <user name> and <password> and tap <Login>. The main menu will appear automatically.
3. Tap <Setup Run>.
4. Enter the sample identification using either the barcode scanner or manual entry. Manual entry can be done by tapping the pencil icon of the <Scan or Enter Sample ID> line.
5. Enter the SBT and *C. difficile* PIE barcodes using the Revogene barcode scanner. Gently positioning the PIE vertically in front of the scanner. Alternatively, SBT and PIE barcodes may be entered manually (tap the pencil icon of their respective lines). Handle the PIE carefully without dropping, shaking or inverting it.
6. (Optional) Tap the pencil icon of the <Add Comments> line and type a comment.
7. Insert the *C. difficile* PIE into the Revogene, at any position of the carousel. The software will automatically associate sample and SBT to the correct *C. difficile* PIE.
8. Confirm that the PIE is inserted into instrument by tapping <OK> on the <insert PIE into instrument> line and repeat steps 4 to 8 for all samples.
9. When all *C. difficile* PIEs are inserted into instrument and MOCK PIEs when necessary, tap <Next>.
10. Scan the retention ring and place it on the carousel. Close the instrument lid.
11. Initiate the test run by tapping <Start>.

\*If MOCK PIEs are not available, use unused assay PIEs filled with un-inoculated SB (BLANK) or with external controls.

### VIEWING AND EXPORTING RESULTS

**NOTE 1:** Refer to the Revogene Operator's Manual<sup>13</sup> for further information regarding the acquisition of test results.

1. Once the run is completed, the lid opens automatically.
2. Tap the home icon.
3. If the Revogene has logged-out, re-enter <user name> and <password> and tap <Login>. The main menu will appear automatically.
4. Tap Results icon to access test results.
5. Tap <Last Run> to see the latest test results.
6. From the <Last Run>, select samples for which results report(s) has (have) to be exported. All samples can be selected in one time by clicking the first box to the left of the "Sample ID" column.
7. Tap <Export> and save where appropriate (e.g. USB key).
8. Remove the retention ring and *C. difficile* PIEs from the Revogene. Used *C. difficile* PIEs should be discarded in appropriate waste containers according to the institution's standard practices.

## REPEAT TESTING PROCEDURE

### UNRESOLVED OR INDETERMINATE RESULT FOR A SPECIMEN

When an Unresolved (UNR) or an Indeterminate (IND) result is obtained for a specimen, a repeat test from the corresponding inoculated SBT must be performed within the specified timeframe described in the **Storage and Stability** section. Only one repeat testing from the SBT is allowed.

Vortex the SBT for a minimum of 15 seconds at maximal speed using a vortex mixer. Using a new pouch, follow steps 9 to 13 of the **Sample Preparation and Handling / PIE Preparation** section, then follow the **Revogene Operation** section.

### UNRESOLVED, INDETERMINATE, FALSE NEGATIVE OR FALSE POSITIVE RESULT FOR AN EXTERNAL CONTROL

When an unresolved, an indeterminate, a false negative or a false positive result is obtained for an external control, the run is invalid. Specimens included in the run should be repeated using the corresponding inoculated SBT, along with freshly prepared external controls, within the specified timeframe described in the **Storage and Stability** section. Refer to the next **Quality Control** section for preparation of fresh external controls.

For the repeat testing using the corresponding inoculated SBT, vortex the SBT for a minimum of 15 seconds at maximal speed using a vortex mixer. Using a new pouch, follow steps 9 to 13 of the **Sample Preparation and Handling / PIE Preparation** section, then follow the **Revogene Operation** section.

## QUALITY CONTROL

Quality control procedures monitor the accuracy and precision of the analytical process. Each laboratory must establish the number, type and frequency of testing control materials per applicable regulations or accrediting agencies. The procedure described below may be employed, if appropriate, based on local policies and procedures.

**NOTE 1:** Separate DTL, DTT, SBT and PIE must be used for each external control preparation.

- Each *C. difficile* PIE contains a Process Control (PrC) that verifies for sample homogenization, sample dilution, cell lysis, inhibition of DNA amplification and assay reagents failure.
- Good laboratory practice recommends the use of control materials. User should follow the appropriate guidelines concerning the running of external controls. It is recommended that one positive external control and one negative external control should be run at least on a daily basis until adequate process validation is achieved with the *C. difficile* assay on the Revogene in each laboratory setting.
- External control materials are not provided by Meridian Bioscience, Inc. External controls are not used by the Revogene software for the purpose of sample test result interpretation. External controls are treated as if they are specimens.
- Various types of external controls are recommended to allow the user to select the most appropriate for its laboratory quality control program. Process and test external control preparations according to the **Sample Preparation and Handling** section.

### Positive External Control:





- A freshly prepared cell suspension of a toxigenic *C. difficile* strain, bearing the *tcdB* gene, from commercially available control material (e.g., ATCC® 43255™) prepared at 0.5 ± 0.05 McFarland and diluted 1/2 in saline (e.g., BD BBL™ Prepared Saline Solution, cat# 221819) is recommended for use as a positive external control.
- Alternatively, a previously characterized stool specimen positive for toxigenic *C. difficile* could also be used as a positive external control.

### Negative External Control:

- A freshly prepared cell suspension of a non-toxicogenic *C. difficile* strain from commercially available control material (e.g., ATCC® 43593™) prepared at 0.5 ± 0.05 McFarland in saline (e.g., BD BBL™ Prepared Saline Solution, cat# 221819) is recommended for use as a negative external control.
- Alternatively, a previously characterized stool specimen negative for toxigenic *C. difficile* could also be used as a negative external control.

## RESULTS INTERPRETATION

The results are computed by the Revogene from measured fluorescent signals and embedded calculation algorithms and are available on the "Results" window. Possible results are:

Sample	Symbol	Result	Interpretation
Patient specimen		Positive	Toxigenic <i>C. difficile</i> target DNA detected.
		Negative	Toxigenic <i>C. difficile</i> target DNA not detected.
		Unresolved	Amplification/detection failure of the process control as well as for the toxigenic <i>C. difficile</i> target DNA. Could be caused by inhibitory specimens, microfluidic or reagent failure. Repeat testing must be performed (refer to the <b>Repeat Testing Procedure</b> section for further guidance).
		Indeterminate	No reportable result due to possible Revogene detection error during the assay processing, the data analysis, or if the run is interrupted by the user. Repeat testing must be performed (refer to the <b>Repeat Testing Procedure</b> section for further guidance).
Positive External Control		Positive	Valid positive external control result.
		Negative	A positive external control that yields a negative result is indicative of a specimen handling/preparation problem. <b>The run is invalid.</b> Review the specimen handling/preparation technique. Repeat testing must be performed (refer to the <b>Repeat Testing Procedure</b> section for further guidance).
		Unresolved	Incorrect positive external control result. <b>The run is invalid.</b> Repeat testing must be performed (refer to the <b>Repeat Testing Procedure</b> section for further guidance).
		Indeterminate	Incorrect positive external control result. <b>The run is invalid.</b> Repeat testing must be performed (refer to the <b>Repeat Testing Procedure</b> section for further guidance).
Negative External Control		Positive	A negative external control that yields a positive assay result is indicative of a specimen handling and/or contamination event. <b>The run is invalid.</b> Review the specimen handling technique. Repeat testing must be performed (refer to the <b>Repeat Testing Procedure</b> section for further guidance).
		Negative	Valid negative external control result.
		Unresolved	Incorrect negative external control result. <b>The run is invalid.</b> Repeat testing must be performed (refer to the <b>Repeat Testing Procedure</b> section for further guidance).
		Indeterminate	Incorrect negative external control result. <b>The run is invalid.</b> Repeat testing must be performed (refer to the <b>Repeat Testing Procedure</b> section for further guidance).

## LIMITATIONS OF THE PROCEDURE

- The *C. difficile* assay must only be used with the Revogene by trained personnel.
- The *C. difficile* assay is not intended to differentiate carriers of *C. difficile* from those with *C. difficile* infection.
- The *C. difficile* assay does not provide susceptibility results. Culture isolates are needed for performing susceptibility testing.
- Performance characteristics of the *C. difficile* assay were established with unformed (liquid or soft) stool specimens collected from patients suspected of having *C. difficile* infection. Use of the *C. difficile* assay for clinical specimen types other than those specified has not been evaluated and performance characteristics are not established.
- Results from the *C. difficile* assay should be used as an adjunct to clinical observations and other information available to the physician.
- Assay results may be affected by concurrent antimicrobial therapy as *C. difficile* DNA may continue to be detected.
- A positive assay result does not necessarily indicate the presence of viable organisms. It is, however, indicative for the presence of toxigenic *C. difficile* DNA.
- Erroneous assay results may occur from improper specimen collection, handling or storage, technical error or sample mix-up. Careful compliance with the instructions of this insert, the Revogene Operator's Manual<sup>13</sup> and to established guidelines is necessary to avoid erroneous results.
- A negative result does not rule out the possibility of *C. difficile* colonization. False negative results may occur when the *C. difficile* concentration is below the limit of detection of the assay. If the patient has signs or symptoms of infection, other laboratory tests and clinical information should be used to confirm the negative result.
- Contamination or false negative results may occur if a PIE cap is incorrectly closed.
- While there are no known strains/isolates of toxigenic *C. difficile* lacking the *tcdB* gene, the occurrence of such a strain could lead to an erroneous result using the *C. difficile* assay.
- Mutations or polymorphisms in primer- or probe-binding regions may affect detection of *C. difficile* *tcdB* gene variants, resulting in a false negative result with the *C. difficile* assay.
- Calcium Carbonate (e.g. Tums®) or Aluminum Hydroxide/Magnesium Hydroxide (e.g. Stomax®) may potentially inhibit the detection of toxigenic *C. difficile* when either of these substances is present in sample buffer at a concentration of > 0.5 mg/mL or > 0.5 µL/mL respectively.
- Presence of *Clostridium sordellii* at a load of > 10<sup>5</sup> CFU/mL of sample buffer may lead to detection of false positive results.
- A combinatory effect of *Enterobacter aerogenes*, *Enterobacter cloacae* subsp. *cloacae*, *Enterococcus faecalis*, *Escherichia coli* O157:H7 and *Helicobacter pylori*, at ≥ 10<sup>6</sup> CFU/mL in sample buffer, may have an inhibitory effect on the detection of toxigenic *C. difficile*.
- A combinatory effect of *Shigella boydii*, *Shigella dysenteriae*, *Shigella sonnei*, *Staphylococcus aureus* subsp. *aureus*, *Streptococcus agalactiae* and *Vibrio parahaemolyticus*, at ≥ 10<sup>6</sup> CFU/mL in sample buffer, may have an inhibitory effect on the detection of toxigenic *C. difficile*.

## EXPECTED VALUES

The prevalence of *C. difficile* infection (CDI) depends upon a variety of factors including predisposition for infection due to prior therapy with broad-spectrum antibiotics, the presence of symptoms and the standard of care test. In a combined prospective and retrospective study, specimens were collected at 8 geographically diverse clinical sites from 2 581 subjects in the range of age between 0–2 years old to more than 60 years old. Of the 2 461 specimens that met all inclusion criteria without meeting any of the exclusion criteria, 333 were positive based on the combined results of direct and enriched toxigenic culture for an observed prevalence of 13.5% [333/2461; 95%CI: 12.2 – 15%]. The percentage of positive results observed with Revogene *C. difficile* assay in the study population was 11.5% [283/2461; 95%CI: 10.3 – 12.8%].

## CLINICAL PERFORMANCE CHARACTERISTICS

The clinical performance of the Revogene *C. difficile* assay was established in an Independent Ethics Committee (IEC)-approved, prospective and retrospective, multi-site, investigation comparing the results with toxigenic culture using leftover, de-identified, unformed (liquid or soft) stool specimens from subjects suspected of having CDI. Two thousand five hundred and eighty-one (2581) specimens were prospectively collected at seven geographically diverse sites across the US and Canada between February 4, 2017 and July 15, 2017 from symptomatic eligible subjects. Only specimens that met the study inclusion criteria and did not meet any of the exclusion criteria were enrolled. Two thousand four hundred and sixty-three (2463) specimens were used to establish the performance of the Revogene *C. difficile* assay by comparison with combined direct and enriched culture method. All 2463 freshly collected specimens were tested in culture, 798 were tested with Revogene *C. difficile* assay as fresh specimens and a subset of 1665 specimens were frozen to be tested later with the Revogene *C. difficile* assay.

The toxigenic culture was performed at a central reference laboratory and the Revogene *C. difficile* assay was performed at the seven proficient designated sites. The toxigenic culture method included direct and enriched culture followed by cytotoxicity testing. The direct culture method consisted in the transfer of a swab of the stool specimen from anaerobic transport medium to pre-reduced selective anaerobic media, a standard Cycloserine Cefoxitin and Fructose Agar plate (CCFA), followed by cytotoxicity testing on *C. difficile* colonies isolated from stool. Briefly, colonies isolated in direct cultures that were morphologically resembling *C. difficile* (i.e. yellow ground-glass-like in appearance with barnyard-like odor) were confirmed for their identity by aero intolerance test on chocolate agar plates and a Blood Agar Plate (BAP) under anaerobic conditions.

Once colonies were identified, a second BAP plate containing a vancomycin (5 mcg) disc was plated, and an inoculum loop was used to inoculate an anaerobic chopped meat broth. The inoculated broth was incubated anaerobically at 35–37 C for 48h for cytotoxicity testing using Cell Cytotoxicity Neutralization Assay (CCNA; Cytotoxicity Assay for *Clostridium difficile* Toxin, Diagnostic Hybrids). For the enriched culture method, the same swab that was utilized to inoculate the CCFA plate was used to inoculate a Cycloserine Cefoxitin Mannitol Broth with Taurocholate And Lysozyme (CCMB-TAL) tube. The enrichment broth was sub-cultured on another CCFA plate and followed the same procedure used for the direct method. A specimen was considered positive for toxigenic *C. difficile* if *C. difficile* was recovered from stool either by direct or enriched culture and if bacterial isolates tested positive by CCNA. If *C. difficile* was isolated from the direct culture and the isolate tested positive by cytotoxicity testing, the enrichment culture was not further analyzed. Specimens were classified as negative for toxigenic *C. difficile* only if they tested negative by both direct, and combined culture i.e. direct and enriched culture.

The sensitivity, specificity, Positive Predictive Value (PPV) and Negative Predictive Value (NPV) values were calculated by comparing Revogene *C. difficile* assay results with the combined results of direct and enriched culture method (reference method). Discrepant analysis was performed on a portion of specimens with discordant results between the Revogene *C. difficile* assay and combined culture method, using four sites' routine PCR assay. Finally, the Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) were determined comparing the Revogene *C. difficile* assay with the direct culture results.

## RESULTS

The study population demographics are presented in Table 1.

Table 1. Study population demographics of all compliant specimens at all levels.

Subjects	All Subjects	Fresh	Frozen
	N=2461	N=797	N=1664
<b>Source of specimen</b>			
In-patient	1804 (73.3%)	617 (77.4%)	1187 (71.3%)
Out-patient	420 (17.1%)	123 (15.4%)	297 (17.8%)
Emergency room	234 (9.5%)	57 (7.2%)	177 (10.6%)
Missing	3 (0.1%)	0 (0.0%)	3 (0.2%)
<b>Age Class</b>			
< 2	9 (0.4%)	4 (0.5%)	5 (0.3%)
3-18	105 (4.3%)	30 (3.8%)	75 (4.5%)
19-60	1199 (48.7%)	399 (50.1%)	800 (48.1%)
> 60	1148 (46.6%)	364 (45.7%)	784 (47.1%)

Of the 798 fresh and 1665 frozen eligible specimens that were compliant at the specimen and PCR level, 9 and 13 were respectively reported unresolved at initial testing (1.1% for the fresh specimens and 0.8% for the frozen specimens) and only one fresh specimen remained unresolved following repeat testing. The unresolved rate after repeat testing was 0.1% (1/798) for the fresh specimens and 0.0% (0/1665) for the frozen specimens.

Of the 798 fresh and 1665 frozen eligible specimens that were compliant at the specimen and PCR level, 12 and 28 were respectively reported indeterminate at initial testing (1.5% for the fresh specimens and 1.7% for the frozen specimens) and only one frozen specimen remained indeterminate following repeat testing. The indeterminate rate after repeat testing was 0.0% (0/798) for the fresh specimen and 0.1% (1/1665) for the frozen specimens.

The overall initial non-reportable rate was 2.6% (21/798) and 0.1% (1/798) after repeat testing for the fresh specimens. The overall initial non-reportable rate was 2.5% (41/1665) and 0.1% (1/1665) after repeat testing for the frozen specimens.

## COMPARISON WITH COMBINED DIRECT AND ENRICHED CULTURE

The clinical performance of the Revogene *C. difficile* assay compared to the combined results of direct and enriched toxigenic culture are presented in Table 2 for the fresh specimens and in Table 3 for the frozen specimens.

Of the 2463 eligible specimens, 2461 had valid results for direct toxigenic culture, enriched toxigenic culture and the Revogene *C. difficile* assay. Of the 2461 specimens, 333 were positive based on the combined results of direct and enriched toxigenic culture for an observed prevalence of 13.5% [333/2461; 95%CI: 12.2 – 15.0%].

Performance characteristics of the *C. difficile* assay obtained with the fresh stool specimens demonstrated 80.5% sensitivity (91/113; 95%CI: 72.0 – 87.4) and 97.1% specificity (664/684; 95%CI: 95.5 – 98.2%) in comparison to the reference method (combined direct and enriched culture method) obtained with fresh stool specimens (Table 2). Performance characteristics of the *C. difficile* assay obtained with the frozen stool specimens demonstrated 87.3% sensitivity (192/220; 95%CI: 82.1 – 91.4%) and 97.3% specificity (1405/1444; 95%CI: 96.3 – 98.1%) in comparison to the reference method (Table 3).

Table 2. Performance characteristics of the Revogene *C. difficile* assay in comparison to the reference method (combined direct and enriched culture method) obtained with fresh stool specimens.

Specimen Types	Site	True Positive	False Positive	False Negative	True Negative	N	Sensitivity	Specificity
Fresh	#2	14	6	2	63	85	87.5% (14/16; 95% CI: 61.7 – 98.5%)	91.3% (63/69; 95% CI: 82.0 – 96.7%)
	#3	18	2	3	135	158	85.7% (18/21; 95% CI: 63.7 – 97.0%)	98.5% (135/137; 95% CI: 94.83 – 99.8%)
	#4	6	0	1	72	79	85.7% (6/7; 95% CI: 42.1 – 99.6%)	100.0% (72/72; 95% CI: 95.0 – 100.0%)
	#5	21	5	4	115	145	84.0% (21/25; 95% CI: 63.9 – 95.5%)	95.8% (115/120; 95% CI: 90.5 – 98.6%)
	#7	16	5	7	148	176	69.6% (16/23; 95% CI: 47.1 – 86.8%)	96.7% (148/153; 95% CI: 92.6 – 98.6%)
	#9	16	2	5	131	154	76.2% (16/21; 95% CI: 52.8 – 91.8%)	98.5% (131/133; 95% CI: 94.7 – 99.8%)
Total		91	20 <sup>A</sup>	22 <sup>B</sup>	664	797		

<sup>A</sup> Of the 20 specimens with false-positive Revogene *C. difficile* assay results relative to the combined direct and enrichment culture 8 were positive and 4 were negative by a second NAAT method (sites' routine PCR assay).

<sup>B</sup> Of the 22 specimens with false-negative Revogene *C. difficile* assay results relative to the combined direct and enrichment culture 13 were negative and 4 were positive by a second NAAT method (sites' routine PCR assay).

Table 3. Performance characteristics of the Revogene *C. difficile* assay in comparison to the reference method (combined direct and enriched culture method) obtained with frozen stool specimens.

Specimen Types	Site	True Positive	False Positive	False Negative	True Negative	N	Sensitivity	Specificity
Frozen	#1	58	14	6	473	551	90.6% (58/64; 95% CI: 80.7 – 96.5%)	97.1% (473/487; 95% CI: 95.2 – 98.4%)
	#3	67	15	11	464	557	85.9% (67/78; 95% CI: 76.2 – 92.7%)	96.9% (464/479; 95% CI: 94.9 – 98.2%)
	#9	67	10	11	468	556	85.9% (67/78; 95% CI: 76.2 – 92.7%)	97.9% (468/478; 95% CI: 96.2 – 99.0%)
	Total		192	39 <sup>C</sup>	28 <sup>D</sup>	1405	1664	

<sup>C</sup> Of the 39 specimens with false-positive Revogene *C. difficile* assay results relative to the combined direct and enrichment culture 17 were positive and 15 were negative by a second NAAT method (sites' routine PCR assay).

<sup>D</sup> Of the 28 specimens with false-negative Revogene *C. difficile* assay results relative to the combined direct and enrichment culture 14 were negative and 12 were positive by a second NAAT method (sites' routine PCR assay).

Statistical analyses (Chi-square test) demonstrated that site poolability was not achieved for specificity and prevalence but was achieved for sensitivity.

**COMPARISON WITH DIRECT CULTURE**

The clinical agreement of the Revogene *C. difficile* assay compared to the results of direct toxigenic culture are presented in **Table 4** for fresh specimens and **Table 5** for frozen specimens.

**Table 4.** Performance characteristics of the Revogene *C. difficile* assay in comparison to the direct culture method obtained with fresh stool specimens.

Specimen Types	Site	True Positive	False Positive	False Negative	True Negative	N	PPA	NPA
Fresh	#2	11	9	1	64	85	91.7% (11/12; 95% CI: 61.5 – 99.8%)	87.7% (64/73; 95% CI: 77.9 – 94.2%)
	#3	14	6	0	138	158	100.0% (14/14; 95% CI: 76.8– 100.0%)	95.8% (138/144; 95% CI: 91.2 – 98.5%)
	#4	4	2	0	73	79	100.0% (4/4; 95% CI: 39.8– 100.0%)	97.3% (73/75; 95% CI: 90.7 – 99.7%)
	#5	11	15	1	118	145	91.7% (11/12; 95% CI: 61.5 – 99.8%)	88.7% (118/133; 95% CI: 82.1 – 93.6%)
	#7	12	9	0	155	176	100.0% (12/12; 95% CI: 73.5– 100.0%)	94.5% (155/164; 95% CI: 89.8 – 97.5%)
	#9	11	7	1	135	154	91.7% (11/12; 95% CI: 61.5– 99.8%)	95.1% (135/142; 95% CI: 90.1 – 98.0%)
	Total	63	48	3	683	797		

**Table 5.** Performance characteristics of the Revogene *C. difficile* assay in comparison to the direct culture method obtained with frozen stool specimens.

Specimen Types	Site	True Positive	False Positive	False Negative	True Negative	N	PPA	NPA
Frozen	#1	44	28	3	476	551	93.6% (44/47; 95% CI: 82.5 – 98.7%)	94.4% (476/504; 95% CI: 92.1 – 96.3%)
	#3	58	24	3	472	557	95.1% (58/61; 95% CI: 86.3 – 99.0%)	95.2% (472/496; 95% CI: 92.9 – 96.9%)
	#9	58	19	2	477	556	96.7% (58/60; 95% CI: 88.5 – 99.6%)	96.2% (477/496; 95% CI: 94.1 – 97.7%)
		Total	160	71	8	1425	1664	

Statistical analyses (Chi-square test) demonstrated that site poolability was achieved for sensitivity, specificity and prevalence.

**ANALYTICAL PERFORMANCE CHARACTERISTICS**

**ANALYTICAL SENSITIVITY**

The analytical sensitivity (Limit of Detection or LoD) of the CDiff assay was determined using clinical liquid stool matrix previously tested negative for toxigenic *C. difficile* and spiked with various concentrations of toxigenic *C. difficile* bacterial suspensions. Two strains of toxigenic *C. difficile* (ATCC® 43255™, Ribotype 087, Toxinotype 0 and ATCC® BAA-1805™, NAP1, Ribotype 027, Toxinotype IIIb) were tested in 24 replicates per concentration. The LoD is defined as the lowest concentration at which 95% or more of all replicates tested positive. In regards to the two strains tested, the LoD of the *C. difficile* assay was 1 500 CFU/mL of Sample Buffer (SB). Results are summarized in **Table 6**.

**Table 6.** LoD of the *C. difficile* assay.

Toxigenic <i>C. difficile</i> strain ATCC® number	LoD (CFU/mL of SB)
ATCC® 43255™	1 500
ATCC® BAA-1805™	1 500

**INCLUSIVITY**

Inclusivity of the *C. difficile* assay was determined for 20 strains of toxigenic *C. difficile* representing 8 different toxinotypes from diverse geographic origins. Each strain was tested from a quantitated cell culture spiked in a toxigenic *C. difficile*-negative liquid stool matrix at a load of 3 750 CFU/mL of SB, corresponding to two to three times the LoD value of ATCC® 43255™ strain. Three replicates per strain were tested using 3 different *C. difficile* kit lots. All toxigenic *C. difficile* strains tested were detected at 3 750 CFU/mL SB. Strains tested are described in **Table 7**.

**Table 7.** Toxigenic *C. difficile* strains tested for inclusivity with the *C. difficile* assay.

Toxigenic <i>C. difficile</i> strains	Toxinotype
ATCC® 9689™	(Toxinotype 0, A+, B+)
ATCC® 700792™	(Toxinotype 0, A+, B+)
ATCC® 17858™	(Toxinotype 0, A+, B+)
ATCC® BAA-1382™	(Toxinotype 0, A+, B+)
ATCC® 51695™	(Toxinotype 0, A+, B+)
ATCC® 43600™	(Toxinotype 0, A+, B+)
ATCC® 43599™	(Toxinotype 0, A+, B+)
ATCC® 43596™	(Toxinotype 0, A+, B+)
ATCC® 43594™	(Toxinotype 0, A+, B+)
ATCC® BAA-1804™	(Toxinotype 0, A+, B+)
ATCC® 17857™	(Toxinotype 0, A+, B+)
ATCC® BAA-1873™	(Toxinotype 0, A+, B+)
ATCC® BAA-1870™	(Toxinotype IIIb, NAP1, A+, B+)
ATCC® BAA-1803™	(Toxinotype IIIc, NAP1, A+, B+)
ATCC® BAA-1875™	(Toxinotype V, A+, B+)
ATCC® 43598™	(Toxinotype VIII, A-, B+)
CCUG 8864	(Toxinotype X, A-, B+)
ATCC® BAA-1812™	(Toxinotype XII, A+, B+)
ATCC® BAA-1814™	(Toxinotype XXII, A+, B+)
ATCC® BAA-2155™	(Toxinotype XXII, A+, B+)

In addition, an *in silico* analysis was performed on July 20<sup>th</sup>, 2017 to assess inclusivity of the primers and probe of the *C. difficile* assay target regarding 52 toxigenic *C. difficile* strains listed in the National Center for Biotechnology Information (NCBI) database. The alignment results showed no mismatch with the selected 52 sequences. The analysis predicted the detection of all these toxigenic *C. difficile* strains.

**CROSS-REACTIVITY**

The cross-reactivity of the *C. difficile* assay was assessed with high loads of organisms that are not targeted by the assay, or phylogenetically related to *C. difficile*, or non-toxigenic *C. difficile* strains, or present in the normal intestinal flora. The study included 50 bacteria, 1 yeast, 7 viruses and human DNA (**Table 8**). Bacteria and yeast were tested at a load of  $\geq 10^5$  CFU/mL of SB. Nucleic acids from 6 viruses and human DNA were tested at a load of  $\geq 10^5$  DNA or RNA copies/mL of SB. These organisms were tested using quantitated cell cultures or nucleic acid solutions spiked in a *C. difficile*-negative liquid stool matrix. Each organism was tested in PCR triplicates.

Under the conditions of the study, *Clostridium sordellii* was detected by the *C. difficile* assay at a load of approximately  $10^5$  CFU/mL of SB for one replicate out of 3, but was found non-reactive at a load of approximately  $10^5$  CFU/mL of SB. *Clostridium novyi* and *Clostridium scindens* strains produced false positive reactions in one replicate out of six tested at approximately  $10^5$  CFU/mL of SB. No reactivity was observed for three replicates tested at  $10^5$  CFU/mL of SB. *Enterococcus faecalis* strain produced false positive reactions in one replicate out of three tested at approximately  $10^7$  CFU/mL of SB. No reactivity was observed for three replicates tested at  $10^5$  CFU/mL of SB. The other organisms and nucleic acids tested were found non-reactive with the *C. difficile* assay.

For the *Coxsackievirus* only, the cross-reactivity with primers and probes of the *C. difficile* assay was evaluated by an *in silico* analysis performed on all strain sequences of *Coxsackievirus* contained in the National Center for Biotechnology Information (NCBI) database between February 7<sup>th</sup>, 2017 and August 11<sup>th</sup>, 2017. The analysis suggested that the *Coxsackievirus* strains should not be reactive with the *C. difficile* assay.

**Table 8.** List of organisms tested for cross-reactivity with the *C. difficile* assay.

Bacteria	
<i>Abiotrophia defectiva</i>	<i>Acinetobacter baumannii</i>
<i>Alcaligenes faecalis</i> subsp. <i>faecalis</i>	<i>Aeromonas hydrophila</i>
<i>Campylobacter jejuni</i> ( <i>Campylobacter coli</i> )	<i>Bacillus cereus</i>
<i>Clostridium bifementans</i>	<i>Bacteroides fragilis</i>
<i>Clostridium butyricum</i>	<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>
<i>Clostridium haemolyticum</i>	<i>Citrobacter freundii</i>
<i>Clostridium novyi</i>	<i>Clostridium perfringens</i>
<i>Flavonifractor plautii</i> ( <i>Clostridium orbiscindens</i> )	<i>Enterobacter aerogenes</i>
<i>Clostridium scindens</i>	<i>Enterobacter cloacae</i> subsp. <i>cloacae</i>
<i>Clostridium septicum</i>	<i>Enterococcus faecalis</i>
<i>Clostridium sordellii</i>	<i>Escherichia coli</i>
<i>Clostridium difficile</i> (non-toxicogenic) – ATCC® 43593™	<i>Klebsiella oxytoca</i>
<i>Clostridium difficile</i> (non-toxicogenic) – ATCC® 43601™	<i>Lactobacillus acidophilus</i>
<i>Clostridium sporogenes</i>	<i>Peptostreptococcus anaerobius</i>
<i>Edwardsiella tarda</i>	<i>Porphyromonas asaccharolytica</i>
<i>Escherichia coli</i> O157:H7	<i>Prevotella melaninogenica</i>
<i>Helicobacter pylori</i>	<i>Proteus mirabilis</i>
<i>Listeria monocytogenes</i>	<i>Pseudomonas aeruginosa</i>
<i>Plesiomonas shigelloides</i>	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i>
<i>Providencia alcalifaciens</i>	<i>Serratia marcescens</i> subsp. <i>marcescens</i>
<i>Salmonella enterica</i> subsp. <i>arizonae</i>	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Choleraesuis</i>	<i>Staphylococcus epidermidis</i>
<i>Serratia liquefaciens</i>	<i>Streptococcus agalactiae</i>
<i>Shigella boydii</i>	<i>Vibrio parahaemolyticus</i>
<i>Shigella dysenteriae</i>	<i>Shigella sonnei</i>
Yeast	
<i>Candida albicans</i>	
Viruses	
Human Adenovirus 1 (DNA)	Rotavirus (RNA)
Enterovirus D68 (RNA)	Norovirus (RNA)
Echovirus 4 (RNA)	Human Herpes virus 5 (Cytomegalovirus) (DNA)
<i>Coxsackievirus (in silico)</i>	
Human DNA	
Human gDNA	

**INTERFERING ORGANISMS**

The potentially inhibitory effect of 30 organisms, that may be present in the normal intestinal flora and which are not targeted by the assay, was assessed using organisms selected from the cross-reactivity study (Table 8). Each organism category (i.e., bacteria, yeast, viruses) was represented with a special attention to include the most frequent causative agents of intestinal tract infections. Groups of 2 to 6 organisms were prepared in toxigenic *C. difficile*-negative liquid stool matrix, and tested in duplicate in presence of either 3 750 CFU/mL of SB of the toxigenic *C. difficile* ATCC® 43255™ strain or 4 500 CFU/mL of SB of the toxigenic *C. difficile* ATCC® BAA-1805™ strain, to assess their potential interference on detection of toxigenic *C. difficile* or PrC. Each organism within group was diluted to reach a load of  $\geq 10^5$  CFU/mL of SB for bacterium and yeast, and  $\geq 10^5$  copies/mL of SB for virus. The 30 organisms included in the study are presented in Table 9.

**Table 9.** List of organisms tested for interference with the *C. difficile* assay.

Group 1	
<i>Aeromonas hydrophila</i>	<i>Bacillus cereus</i>
<i>Bacteroides fragilis</i>	<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>
<i>Campylobacter jejuni</i> ( <i>Campylobacter coli</i> )	
Group 2	
<i>Candida albicans</i>	<i>Citrobacter freundii</i>
<i>Clostridium difficile</i> (non-toxicogenic) – ATCC® 43593™	<i>Clostridium difficile</i> (non-toxicogenic) – ATCC® 43601™
<i>Clostridium perfringens</i>	<i>Clostridium sordellii</i>
Group 3	
<i>Enterobacter aerogenes</i>	<i>Enterobacter cloacae</i> subsp. <i>cloacae</i>
<i>Enterococcus faecalis</i>	<i>Escherichia coli</i> O157:H7
<i>Helicobacter pylori</i>	
Group 4	
<i>Lactobacillus acidophilus</i>	<i>Peptostreptococcus anaerobius</i>
<i>Plesiomonas shigelloides</i>	<i>Proteus mirabilis</i>
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Choleraesuis</i>	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i>
Group 5	
<i>Shigella boydii</i>	<i>Shigella dysenteriae</i>
<i>Shigella sonnei</i>	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>
<i>Streptococcus agalactiae</i>	<i>Vibrio parahaemolyticus</i>
Group 6	
Rotavirus (RNA)	Norovirus (RNA)

None of the 30 organisms present at  $\geq 10^5$  CFU/mL of SB for bacteria and yeast and  $\geq 10^5$  copies/mL of SB for viruses interfered with detection of PrC and with the toxigenic *C. difficile* ATCC® BAA-1805™ strain.

Group 3 and group 5 showed a potentially inhibitory effect on detection of the toxigenic *C. difficile* strain ATCC® 43255™. Nevertheless, when each bacterium from these groups was tested individually at a load of  $\geq 10^6$  CFU/ mL of SB in presence of *C. difficile* strain ATCC® 43255™, none interfered.

## INTERFERING SUBSTANCES

The potentially inhibitory effect of 16 exogenous and 5 endogenous substances that may be present in the intestinal tract was assessed using toxigenic *C. difficile* negative samples and toxigenic *C. difficile* strains ATCC® 43255™ and ATCC® BAA-1805™ at two to three times the LoD (3 750 CFU/mL of SB and 4 500 CFU/mL of SB respectively) in presence of liquid stool matrix. Substances were tested at their potentially highest concentration that could be found in a stool specimen. The results for the 21 substances are presented in **Table 10**.

Results demonstrated no reportable interference on PrC. Calcium Carbonate (e.g. Tums®) and Aluminum Hydroxide/Magnesium Hydroxide (e.g. Stomax®) showed a potentially inhibitory effect on the detection of toxigenic *C. difficile* when either of these substances was present in SBT at a concentration of 5 mg/mL (0,5% W/V) and 5 µL/mL (0,5% V/V) respectively. When tested at 0,5 mg/mL (0,05% W/V) or 0,5 µL/mL (0,05% V/V) respectively, these substances showed no reportable interference with the *C. difficile* assay.

**Table 10.** List of exogenous and endogenous substances tested with the *C. difficile* assay.

Exogenous substances		
Substance (commercial name)	Concentration or amount in SBT <sup>1</sup>	Results <sup>2</sup>
Vaginal antifungal / anti-itch (Nystatin)	0,5% W/V	NI
Creams / ointments (Personnelle Hydrocortisone cream)	0,5% V/V	NI
Anti-hemorrhoidal creams / ointments (Preparation H®)	0,5% V/V	NI
Antacids (Tums®)	0,5% W/V	I <sup>3</sup>
Antacids (Stomax®)	0,5% V/V	I <sup>4</sup>
Enemas (Life BRAND™ Heavy Mineral Oil USP)	0,5% V/V	NI
Enemas (Mesalazine or 5-aminosalicylic Acid)	0,5% W/V	NI
Condom with spermicidal lubricant (Trojan® with spermicidal lubricant condom)	Square of 2 mm <sup>2</sup>	NI
Anti-diarrheal medication (Pepto Bismol™)	0,5% V/V	NI
Anti-diarrheal medication (Imodium®)	0,5% V/V	NI
Laxatives (Senokot®)	0,5% V/V	NI
Oral and topical antibiotics (Vancomycin)	0,5% V/V	NI
Oral and topical antibiotics (Metronidazole)	0,5% W/V	NI
Non-steroidal anti-inflammatory (Aleve®)	0,5% W/V	NI
Moist towelettes (Equate™ Flushable Moist Wipes)	Square of 2 mm <sup>2</sup>	NI
Moist towelettes (Wet Ones®)	Square of 2 mm <sup>2</sup>	NI
Endogenous substances		
Substance	Concentration or amount in SBT <sup>1</sup>	Results <sup>2</sup>
Fecal fat, triglycerides mix (C2-C10)	0,5% V/V	NI
Fecal fat, Palmitic acid	1,0% W/V	NI
Fecal fat, Stearic acid	0,5% W/V	NI
Whole blood	0,5% V/V	NI
Mucus	0,5% V/V	NI

<sup>1</sup>W/V: Weight/Volume; V/V: Volume/Volume

<sup>2</sup>I: Interference with the *C. difficile* assay; NI: No Interference with the *C. difficile* assay

<sup>3</sup>No Interference at 0,05% W/V

<sup>4</sup>No interference at 0,05% V/V

## CARRY-OVER AND CROSS-CONTAMINATION

The within-run and between-run carry-over and cross-contamination were assessed using positive samples prepared in a toxigenic *C. difficile*-negative liquid stool matrix to reach a final concentration of > 10<sup>7</sup> CFU/mL of SB of the toxigenic *C. difficile* ATCC® 43255™ strain. Negative samples, prepared with the toxigenic *C. difficile*-negative liquid stool matrix only, were also tested.

For the within-run study, a total of 10 runs were performed by two operators with the *C. difficile* assay on one Revogene. Four (4) high positive samples and 4 negative samples were tested by alternating positive and negative samples in each run. For the between-run study, a run of 8 replicates of high positive samples followed by a run of 8 replicates of negative samples were performed by two operators, for a total of 10 runs on one Revogene.

Absence of carry-over and cross-contamination was demonstrated.

## REPRODUCIBILITY

Between-site reproducibility study was performed on three sites, by two operators per site, over five distinct days using one *C. difficile* assay kit lot.

Between-lot reproducibility study was performed on one site, by two operators, over 15 days using three *C. difficile* assay kit lots (5 days per kit lot). For each reproducibility study, a total of 120 replicates for negative samples and 90 replicates for each category of positive samples, all prepared in a toxigenic *C. difficile*-negative liquid stool matrix, were tested. Two toxigenic *C. difficile* strains were used for positive samples: ATCC® 43255™ (toxintype 0, Ribotype 087) and ATCC® BAA-1805™ (toxintype IIIb, NAP1, Ribotype 027).

Sample categories were described as follows:

1. Low Positive (LP): ATCC 43255 and ATCC BAA-1805 strains were spiked at 2 438 CFU/mL and 2 925 CFU/mL of SB, respectively
2. Moderate Positive (MP): ATCC 43255 and ATCC BAA-1805 strains were spiked at 3 750 and 4 500 CFU/mL of SB, respectively
3. Negative (NEG): samples without toxigenic *C. difficile* strain

Results for each sample category tested during the between-site and the between-lot reproducibility studies are shown in **Tables 11 and 12** respectively. For the between-site reproducibility, the overall percent agreement was 100% for NEG, LP and MP of toxigenic *C. difficile* ATCC® BAA-1805™ strain. The overall percent agreement for LP and MP of toxigenic *C. difficile* ATCC® 43255™ strain was 94,4% and 96,7% respectively (**Table 11**). For the between-lot reproducibility, the overall percent agreement was 100% for NEG and MP of toxigenic *C. difficile* ATCC® BAA-1805™ strain categories and 98,9% for LP. The overall percent agreement for LP and MP of toxigenic *C. difficile* ATCC® 43255™ strain was 91,1% and 96,7% respectively (**Table 12**).

Overall mean Cycle Threshold (Ct) values with variance components (SD and %CV) are shown in **Tables 11 and 12**.

**Table 11.** Between-site reproducibility study results using one *C. difficile* assay kit lot.

Category	Toxigenic <i>C. difficile</i> strain (ATCC® number)	Site 1		Site 2		Site 3		Overall Results/ Total	Overall Percent Agreement <sup>1</sup>	Overall 95% CI	Ct Values <sup>2</sup>		
		Results/ Total	Percent Agreement	Results/ Total	Percent Agreement	Results/ Total	Percent Agreement				Overall Mean	SD	%CV
LP	ATCC® 43255™	27/30	90,0%	28/30	93,3%	30/30	100%	85/90	94,4%	87,5%-98,2%	38,48	1,56	4,05
	ATCC® BAA-1805™	30/30	100%	30/30	100%	30/30	100%	90/90	100%	96,7%-100%	36,37	2,74	7,54
MP	ATCC® 43255™	29/30	96,7%	29/30	96,7%	29/30	96,7%	87/90	96,7%	90,6%-99,3%	37,68	1,63	4,32
	ATCC® BAA-1805™	30/30	100%	30/30	100%	30/30	100%	90/90	100%	96,7%-100%	36,57	1,65	4,52
NEG	N/A	40/40	100%	40/40	100%	40/40	100%	120/120	100%	97,5%-100%	33,09	1,02	3,09

<sup>1</sup>For the NEG category, percent agreement was calculated for negative results.

<sup>2</sup>For the LP and MP categories, Ct values reported are for the toxigenic *C. difficile* target. For the NEG category, Ct values reported are for the PrC.

**Table 12.** Between-lot reproducibility study results on one site using three *C. difficile* assay kit lots.

Category	Toxigenic <i>C. difficile</i> strain (ATCC® number)	Lot 1		Lot 2		Lot 3		Overall Results/ Total	Overall Percent Agreement <sup>1</sup>	Overall 95% CI	Ct Values <sup>2</sup>		
		Results/ Total	Percent Agreement	Results/ Total	Percent Agreement	Results/ Total	Percent Agreement				Overall Mean	SD	%CV
LP	ATCC® 43255™	27/30	90,0%	27/30	90,0%	28/30	93,3%	82/90	91,1%	83,2%-96,1%	38,57	1,72	4,46
	ATCC® BAA-1805™	30/30	100%	29/30	96,7%	30/30	100%	89/90	98,9%	94,0%-100%	37,02	1,98	5,35
MP	ATCC® 43255™	29/30	96,7%	30/30	100%	28/30	93,3%	87/90	96,7%	90,6%-99,3%	37,55	2,16	5,75
	ATCC® BAA-1805™	30/30	100%	30/30	100%	30/30	100%	90/90	100%	96,7%-100%	36,86	1,31	3,56
NEG	N/A	40/40	100%	40/40	100%	40/40	100%	120/120	100%	97,5%-100%	33,21	0,90	2,71

<sup>1</sup>For the NEG category, percent agreement was calculated for negative results.

<sup>2</sup>For the LP and MP categories, Ct values reported are for the toxigenic *C. difficile* target. For the NEG category, Ct values reported are for the PrC.

**PRECISION**

Precision study was performed on one site, by two operators, over 20 days using one *C. difficile* assay kit lot. For each sample category, 2 runs per day with 2 replicates per run were performed for a total of 80 replicates.

For each category, all samples were prepared in a toxigenic *C. difficile*-negative liquid stool matrix. Two toxigenic *C. difficile* strains were used for positive samples: ATCC® 43255™ (toxintype 0, Ribotype 087) and ATCC® BAA-1805™ (toxintype IIIb, NAP1, Ribotype 027).

Sample categories were described as follows:

1. Low Positive (LP): ATCC® 43255™ and ATCC® BAA-1805™ strains were spiked at 2 925 CFU/mL of SB
2. Moderate Positive (MP): ATCC® 43255™ and ATCC® BAA-1805™ strains were spiked at 4 500 CFU/mL of SB
3. Negative (NEG): samples without toxigenic *C. difficile* strain

Precision study results for NEG and LP of toxigenic *C. difficile* ATCC® 43255™ and *C. difficile* ATCC® BAA-1805™ strains demonstrated 100% agreement. Precision study results for MP of toxigenic *C. difficile* ATCC® 43255™ and *C. difficile* ATCC® BAA-1805™ strains demonstrated agreement of 98,8% (Table 13).

**Table 13.** Precision study percent agreement on one site using one *C. difficile* assay kit lot.

Category	Toxigenic <i>C. difficile</i> strain ATCC® number	Overall Percent Agreement <sup>1</sup>	Overall 95% CI
LP	ATCC® 43255™	100%	96,3%-100%
	ATCC® BAA-1805™	100%	96,3%-100%
MP	ATCC® 43255™	98,8%	93,2%-100%
	ATCC® BAA-1805™	98,8%	93,2%-100%
NEG	N/A	100%	96,3%-100%

<sup>1</sup>For the NEG category, percent agreement was calculated for negative results.

**E-LABELING**

Documentation related to this product can be accessed online at [www.meridianbioscience.com/pi](http://www.meridianbioscience.com/pi). Additionally, paper copies are available upon request by contacting your local distributor or via the phone number listed on the kit box.

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#### REFERENCES

1. Rousseau et al. 2012. *Clostridium difficile* carriage in healthy infants in the community: a potential reservoir for pathogenic strains. Clin Infect Dis. 55(9);1209-1215.
2. Collignon et al. 1993. Heterogeneity of *Clostridium difficile* isolates from infants. Eur J Pediatr. 152(4);319-322.
3. Kyne et al. 2002. Health care costs and mortality associated with nosocomial diarrhea due to *Clostridium difficile*. Clin Infect Dis. 34(3);346-353.
4. Loo et al. 2005. A predominantly clonal multi-institutional outbreak of *Clostridium difficile*-associated diarrhea with high morbidity and mortality. N Engl J Med. 353(23);2442-2449. (Erratum, N Engl J Med 354(20);2200, 2006.)
5. McDonald et al. 2005. An epidemic, toxin gene-variant strain of *Clostridium difficile*. N Engl J Med. 353(23);2433-2441.
6. O'Brien et al. 2007. The emerging infectious challenge of *Clostridium difficile*-associated disease in Massachusetts hospitals: clinical and economic consequences. Infect. Control Hosp. Epidemiol. 28(11);1219-1227.
7. Bartlett. 2006. Narrative review: the new epidemic of *Clostridium difficile*-associated enteric disease. Ann. Intern. Med. 145(10);758-764.
8. Drudy et al. 2007. Emergence and control of fluoroquinolone-resistant, toxin A-negative, toxin B-positive *Clostridium difficile*. Infect. Control Hosp. Epidemiol. 28(8);932-940.
9. Cohen et al. 2010. Clinical practice guidelines for *Clostridium difficile* infection in adults: 2010 Update by the Society for Healthcare Epidemiology of America (SHEA) and the Infectious Diseases Society of America (IDSA). Infect. Control Hosp. Epidemiol. 31(5);431-455.
10. Bartlett JG. 2002. Clinical practice. Antibiotic-associated diarrhea. N Engl J Med. 346(5);334-339.
11. Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health. Biosafety in microbiological and biomedical laboratories-5th ed. HHS Publication No. (CDC) 21- 1112. Revised December 2009
12. Clinical and Laboratory Standards Institute. Protection of laboratory workers from occupationally acquired infections; Approved Guideline-4th ed.
13. SN134822 Revogene® Operator's Manual.

The purchase of this product grants the purchaser rights under certain Roche patents to use it solely for providing human *in vitro* diagnostic services. No general patent or other license of any kind other than this specific right of use from purchase is granted hereby.

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REV. 01/20


Manufactured By

**Meridian Bioscience, Inc.**  
**Corporate Office**  
**3471 River Hills Drive**  
**Cincinnati, Ohio 45244 USA**  
**Telephone: 513.271.3700**  
**Orders/Customer Service:**  
**800.543.1980**  
**Technical Support Center:**  
**800.343.3858**  
**Information Fax: 513.272.5432**  
**Ordering Fax: 513.271.0124**

**INTERNATIONAL SYMBOLS USAGE**

You may see one or more of these symbols on the labeling/packaging of this product:

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For technical assistance call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1960.

TaqMan is a registered trademark of Roche Molecular Systems, Inc.  
 ATCC is a trademark of the American Type Culture Collection. BBL is a trademark of Becton, Dickinson and Company.  
 Revogene and associated logos are trademarks of Meridian Bioscience, Inc.  
 © 2020-01 Meridian Bioscience, Inc.  
 Made in Canada

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TaqMan est une marque commerciale déposée de Roche Molecular Systems, Inc.
Revogene et les logos connexes sont des marques déposées de Meridian Bioscience, Inc.