REF 410500

IVD For *in vitro* Diagnostic Use

Rx Only

INTENDED USE

The Revogene® Carba C assay, performed on the Revogene instrument, is a qualitative *in vitro* diagnostic test designed for the detection and differentiation of the *bla_{KPC}*, *bla_{NDM}*, *bla_{VIM}*, *bla_{OXA-48-like}*, and *bla_{IMP}* gene sequences associated with carbapenem-non-susceptible pure colonies of *Enterobacteriaceae*, *Acinetobacter baumannii*, or *Pseudomonas aeruginosa*, when grown on blood agar or MacConkey agar. The test utilizes automated real-time Polymerase Chain Reaction (PCR).

The Revogene Carba C assay should be used in conjunction with other laboratory tests including phenotypic antimicrobial susceptibility testing. A negative Revogene Carba C assay result does not preclude the presence of other resistance mechanisms.

The Revogene Carba C assay is intended as an aid for infection control in the detection of carbapenem-non-susceptible bacteria that colonize patients in healthcare settings. The identification of a *bla_{IMP}*, *bla_{NDM}* or *bla_{VIM}* metallo-β-lactamase gene (i.e., the genes that encode the IMP, NDM and VIM metallo-β-lactamases, respectively) may be used as an aid to clinicians in determining appropriate therapeutic strategies for patients with known or suspected carbapenem non-susceptible infections.

SUMMARY AND EXPLANATION

Resistance to carbapenems may be due to different mechanisms, and one of them is the production of carbapenemases. Carbapenemases are heterogeneous enzymes capable of hydrolyzing carbapenem antibiotics¹. Genes encoding carbapenemases are mainly located on mobile elements such as plasmids, facilitating their dissemination between different bacterial species. Plasmids carrying genes encoding carbapenemases of classes A (KPC), B (IMP, NDM, VIM) and D (OXA-48-like) have spread beyond healthcare institutions to cause critical community-acquired infections². Bacteria harboring those plasmid genes are often referred as carbapenemase-producing organisms (CPO) or carbapenem-resistant organisms (CRO).

Gram-negative *Enterobacteriaceae* bacilli are part of the normal gut flora and are reported as hosts of carbapenemase genes. *Enterobacteriaceae* that are non-susceptible to carbapenem antibiotics through the production of carbapenemase(s) are known as carbapenemase-producing *Enterobacteriaceae* (CPE). CPE can cause infections in almost any body part including bloodstream and urinary tract infections, ventilator-associated pneumonia and intra-abdominal abscesses³. The treatment of infections caused by these organisms is extremely difficult because of their multidrug resistance and hence, results in high mortality rates⁴.

Pseudomonas aeruginosa and *Acinetobacter baumannii* are two (2) major opportunistic pathogens of the non-fermentative Gram-negative bacilli group. As intrinsic antibiotic resistance to multiple antibiotic classes is common with these organisms, carbapenems are often used for treatment⁵. Non-fermentative Gram-negative bacilli are increasingly acquiring resistance to carbapenem antibiotics as well, leaving few therapeutic options. Non-fermentative Gram-negative bacteria that are non-susceptible to carbapenem(s) through the production of carbapenemase(s) are known as carbapenemase-producing non-fermenters (CP-NF). CP-NF can be recovered from all body sites but are most commonly associated with healthcare-acquired pneumonias⁶.

Early identification and isolation of colonized or infected patients are key infection control measures to limit the spread of these organisms in healthcare settings and to help in a better guidance of treatment.

PRINCIPLE OF THE PROCEDURE

The Revogene Carba C assay can provide results from one (1) up to eight (8) samples in approximately 70 minutes using characterized carbapenem-non-susceptible isolated colonies of *Enterobacteriaceae*, *Acinetobacter baumannii*, or *Pseudomonas aeruginosa*. The assay minimizes operator intervention from the loading of the sample into the single-use microfluidic cartridge (PIE) to final result.

The test is performed using the Revogene instrument which automates sample homogenization and dilution, cells lysis, DNA amplification and detection of the amplified PCR products. User intervention is only required for the preparation of a standardized bacterial suspension from characterized carbapenem-non-susceptible isolated colonies, inoculation of the sample into the Sample Buffer Tube (SBT), transfer of the sample from the SBT into the PIE, and insertion of the PIES into the Revogene carousel.

Each PIE is a fully integrated and closed device, in which a sample is dispensed and processed through different microfluidic chambers and channels which allow for the sample processing (i.e. sample homogenization, sample dilution, cells lysis and DNA extraction) and subsequent real-time PCR steps (Figure 1). The liquid from a single sample is transferred by centrifugation from one (1) chamber to the next in sequence and all reagents specific for the PCR are incorporated and dried within the PCR wells. A Process Control (PrC) is incorporated into each PIE to verify sample processing and amplification steps including the verification of potential inhibitory substances as well as microfluidic, instrument or reagent failure. The amplified products are detected in real-time using target-specific TaqMan® chemistry-based probes.


Figure 1. Top view of a PIE.
1: Sample Loading Chamber, 2: Overflow Chamber,
3: Homogenization Chamber containing PrC, 4: Dilution/Lysis Chamber,
5: Three (3) PCR Wells (A to C from left to right) and one (1) Waste Chamber (at the right end).

The Revogene can process from one (1) up to eight (8) samples simultaneously in the same run. The carousel must contain eight (8) PIES to maintain thermodynamic balance throughout the run. At run completion, the results are computed by the system from measured fluorescent signals and embedded calculation algorithms. Results are displayed on the touchscreen, and may be printed, transferred and stored by the user using the USB port or the connectivity option.

REAGENTS AND MATERIALS PROVIDED

The Revogene Carba C kit contains sufficient reagents and materials to process 24 samples. The kit contains the following materials:

1. 24 Disposable Transfer Tools (DTT): Plastic pipette with minimal and maximal volume marks for transferring the sample from the SBT to the PIE.
2. 24 Sample Buffer Tubes (SBT): Tube containing Tris-HCl pH 8.0/EDTA.Na₂ (TE 1X) buffered solution as a dilution and preservation buffer for sample.
3. 24 individual pouches containing one (1) Carba C microfluidic cartridge (PIE): Integrated device which comprises dried reagents allowing sample process and real-time PCR steps for *bla_{KPC}*, *bla_{NDM}*, *bla_{VIM}*, *bla_{OXA-48-like}*, and *bla_{IMP}* gene sequence amplification and detection. Each PIE contains PrC, PrC-specific primers and probe, *bla_{KPC}*, *bla_{NDM}*, *bla_{VIM}*, *bla_{OXA-48-like}*, and *bla_{IMP}* gene-specific primers and probes, dNTPs, buffer and DNA polymerase.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Revogene® (cat# 610210)
2. Disposable powderless gloves
3. Disposable sterile swabs
4. Saline solution (0.45% to 0.9% NaCl recommended)
5. Blood agar or MacConkey agar plates
6. 10 µg meropenem disc (BD BBL™ Sensi-Disc™ Antimicrobial Susceptibility Test Discs, cat# 231704 or equivalent)
7. Vortex mixer with a maximal speed of at least 3200 rpm (VWR cat# 58816-121 or equivalent)
8. Densitometer (VWR cat# 89402-910 or equivalent) or card with a white background and contrasting black lines (Wickerham Card for turbidity reading, 2x3 inches, black stripes recommended; cat# Z08 or equivalent)
9. Calibrated micropipette (P20 recommended, VWR cat# 89079-964 or equivalent)
10. DNase/RNase-free, aerosol resistant micropipette tips (extended length tips recommended; Sarstedt cat# 70.1189.215 or equivalent)
11. Sample rack (cat# 132539; optional)
12. MOCK PIE(s) (cat# 610208; optional)

WARNING AND PRECAUTIONS

1. The Revogene Carba C assay can only be used on the Revogene.
2. Do not use the kit if the label that seals the outer box is broken upon arrival.
3. Do not use PIEs if the protective pouches are open or broken upon arrival.
4. Do not interchange DTT, SBT, and PIE between kit lots.
5. Each single-use PIE and DTT are used to process one (1) sample. Do not reuse PIE and DTT.
6. Always handle samples as if they are infectious and in accordance with Good Laboratory Practices such as those described in Biosafety in Microbiological and Biomedical Laboratories⁷ and in Clinical and Laboratory Standards Institute (CLSI) Document M29-A4⁸.
7. Wear disposable powderless gloves while handling samples and thoroughly wash hands afterwards.
8. The PIE contains dried reagents. The protective pouch should not be opened until ready to perform the test.
9. Dispose of unused reagents and waste in accordance with country, federal, provincial, state and local regulations.
10. Do not open or break apart the PIE after use. The cap and the seals in the PIE prevent contamination with amplification products and/or infectious particles.
11. Do not use a PIE that has been dropped, shaken or inverted after the sample has been loaded as this may cause invalid results.
12. Do not use a kit that has passed its stated expiration date.
13. Do not refrigerate the loaded PIE.
14. Each run must be performed with eight (8) PIEs in the Revogene carousel to maintain thermodynamic and mechanical balance within the run. Place MOCK PIEs or assay PIEs loaded with Sample Buffer in empty positions if less than eight (8) samples are tested.
15. Local, state, and federal rules and regulations for notification of reportable diseases are continually updated and specify a number of organisms that are important for surveillance and outbreak investigation^{9, 10}. Laboratories are responsible for following their state and/or local rules pertaining to reportable pathogens and should consult their local and/or state public health laboratories for isolate and/or clinical sample submission guidelines.

HAZARD AND PRECAUTIONARY STATEMENTS

There are no known hazards associated with this product.

STORAGE AND STABILITY

1. Store the Revogene Carba C kit at 2-25°C. The expiration date is indicated on the kit box label.
2. Do not open a pouch until ready to perform testing. Use the PIE within one (1) hour after opening the pouch.
3. Inoculated SBT can be stored at 25°C for up to four (4) days or, at 2-8°C for up to seven (7) days.

INSTRUCTION FOR USE

SAMPLE PREPARATION

1. Organisms must be identified as *Enterobacteriaceae*, *Pseudomonas aeruginosa* or *Acinetobacter baumannii*, and carbapenem non-susceptibility status must be determined in accordance with the antibiotic disc package insert or the latest version of CLSI M100 Standards¹¹ prior to testing with the Revogene Carba C assay.
2. Inoculate the organism onto a blood agar or a MacConkey agar plate, streak for isolation, and place a 10 µg meropenem disc in the first quadrant as a means to ensure that the isolate retains its non-susceptibility to carbapenem.
3. Incubate the plate at 35°C for 18–24 hours in aerobic conditions.
4. Prepare a standardized 0.5 McFarland suspension in saline solution by collecting carbapenem-non-susceptible isolated colonies selected from the incubated agar plates close to the meropenem disc using a sterile swab.
5. Adjust the suspension to achieve a turbidity equivalent to a 0.5 McFarland standard. Use either a photometric device or, if performed visually, adequate light to compare the suspension and the 0.5 McFarland standard against a card with a white background and contrasting black lines. This results in a suspension containing approximately 1 to 2 x 10⁶ Colony-Forming Units (CFU)/mL.

SBT PREPARATION

1. For each sample to be tested, obtain one (1) SBT from the kit box.
2. Identify (or label) the SBT with the appropriate sample identification without obscuring or writing over the barcode. Place the SBT on the Sample rack, if used.
3. Vortex the 0.5 McFarland suspension at maximal speed for 15 seconds.
4. Using a micropipette, aspirate 15 µL of the standardized 0.5 McFarland suspension.
5. Remove the cap from the SBT and dispense the McFarland aliquot into the SBT, taking care not to produce aerosols. Pipet liquid up and down to ensure complete transfer of the aliquot. Make sure only one (1) SBT is open at any given moment.
6. Replace the cap on the SBT, close tightly, and replace the SBT on the Sample rack (if used).
7. Prepare any additional samples for testing by repeating steps 1 to 6.
8. When all samples are processed, proceed to step 9 (PIE PREPARATION).

PIE PREPARATION

NOTE: Process one (1) sample at a time.

9. Mix the SBT for a minimum of 15 seconds at maximal speed (≥ 3200 rpm) using a vortex mixer.
10. Unseal the PIE pouch and remove the PIE. Once unsealed, the PIE must be used within one (1) hour.
11. Obtain one DTT from the kit box and aspirate the inoculated Sample Buffer (SB) by squeezing the entire bulb. The liquid level into the DTT must be anywhere between the two marks (Figure 2). If the liquid level is not between the two marks, discharge the SB completely into the SBT and repeat this step.
12. Completely discharge the SB volume contained in the DTT into the Sample Loading Chamber of the PIE (Figure 1). Make sure not to touch the outer edges or the bottom of the Sample Loading Chamber with the DTT.
13. Close the cap of the PIE tightly. Place the PIE on the Sample rack, if used. Do not refrigerate the loaded PIE.
14. Prepare all additional samples for testing by repeating steps 9 to 13 then proceed to step 1 of the REVOCENE SYSTEM OPERATION section.


Figure 2.
Representation of an appropriate Sample Buffer (SB) level using the Disposable Transfer Tool (DTT).

REVOCENE SYSTEM OPERATION

NOTE 1: Each run must be performed with eight (8) PIEs in the Revogene. When less than eight (8) samples are processed, fill the empty places with MOCK PIEs*.

NOTE 2: Refer to the Revogene Operator's Manual¹² for further information regarding Revogene set-up and operation.

1. Power on the Revogene (if not already done). The software will launch automatically.
2. Log in by entering <Username> and <Password> and tap <Login>. The main menu will appear automatically.
3. Tap <Setup Run>.
4. Enter the sample identification number using either the barcode scanner or manual entry. Manual entry can be done by tapping the pencil icon of the <Scan or Enter Sample ID> line.
5. Enter the SBT and PIE barcodes using the barcode scanner by gently positioning the PIE almost vertically in front of the scanner. Alternatively, SBT and PIE barcodes may be entered manually (tap the pencil icon of their respective lines). Handle the PIE carefully without shaking or dropping it.
6. (Optional) Tap the pencil icon of the <Add Comments> line and type a comment.
7. Insert the PIE into the Revogene at any position of the carousel. The software will automatically associate sample and SBT to the correct PIE.
8. Confirm that the PIE is inserted by tapping <OK> on the <Insert PIE into instrument> line and repeat steps 4 to 8 for all samples. If less than eight (8) PIEs are being tested, load MOCK PIEs* in the carousel remaining positions. No scan is required when inserting MOCK PIEs into the Revogene.
9. When all PIEs are inserted into the carousel, tap <Next>.
10. If you have entered MOCK PIEs into the carousel, follow the instructions on the screen.
11. Scan the retention ring and place it on the carousel. Close the instrument lid.
12. Initiate the test run by tapping <Start>.
13. Store inoculated SBT in appropriate conditions (refer to STORAGE AND STABILITY section) for further repeat testing, if needed.

*If MOCK PIEs are not available, use new assay PIEs filled with SB (BLANK).

VIEWING AND EXPORTING RESULTS

- Once the run is completed, the lid opens automatically.
- If the Revogene has logged out, re-enter <Username> and <Password> and tap <Login>. The main menu will appear automatically.
- Tap Results icon to access test results. The Results window shows the reported result(s) for each sample (**Figure 3**).

Sample ID	Patient ID	Started	Assay	Result
1295227899		11/23/2018 2:53 PM	CARBA	
1295227900		11/23/2018 2:53 PM	CARBA	
1295227901		11/23/2018 2:53 PM	CARBA	
1295227902		11/23/2018 2:53 PM	CARBA	
1295227903		11/23/2018 2:53 PM	CARBA	
1295227904		11/23/2018 2:53 PM	CARBA	
1295227905		11/23/2018 2:53 PM	CARBA	
1295227906		11/23/2018 2:53 PM	CARBA	

Results:

IMP	Negative
KPC	Positive
NDM	Negative
OXA-48-like	Negative
VIM	Negative

Figure 3.
Results window showing the reported result(s) for each sample.

- Tap <Last Run> to see the latest test results.
- Tap <Reports> to see the specific target result for each sample. Refer to the RESULTS INTERPRETATION section in order to confirm if any additional action needs to be taken according to the target(s) result(s) obtained.
NOTE: In the case of a repeat test of a sample, all positive results from both the initial and the repeat test should be reported, even if there is a discrepancy between the results for the same target.
- From the <Last Run>, select samples for which results report(s) has (have) to be exported. All samples can be selected simultaneously by checking the box on the upper left corner of the screen.
- Tap <Exports> and save where appropriate (e.g. USB key).
- (Optional) Tap <Search> to find a specific sample and its result.
- Remove the retention ring and PIEs from the Revogene. The used PIEs should be disposed in the appropriate waste containers according to your institution's standard practices.

QUALITY CONTROL

Quality control procedures monitor the accuracy and precision of the analytical process. Each laboratory must establish the number, type and frequency of testing control materials per applicable regulations or accrediting agencies. The procedure described below may be employed, if appropriate, based on local policies and procedures.

PROCESS CONTROL

Each PIE contains a Process Control (PrC) that verifies the sample homogenization, sample dilution, cell lysis, DNA amplification inhibition and assay reagents failure.

EXTERNAL CONTROLS

NOTE: Separate DTT, SBT and PIE must be used for each External Control.

- Good laboratory practice recommends the use of control materials. User should follow the appropriate guidelines concerning the running of External Controls. It is recommended that one (1) Positive External Control and one (1) Negative External Control should be run at least on a daily basis until adequate process validation is achieved with the Revogene Carba C assay on the Revogene in each laboratory setting.
- External Control materials are not provided by Meridian Bioscience, Inc. External Controls are not used by the Revogene software for the purpose of sample test result interpretation. External Controls are treated as if they are specimens.

Five (5) different Positive External Controls and one (1) Negative External Control are recommended and listed below.

Positive External Controls

A standardized 0.5 McFarland suspension of isolated colonies from characterized carbapenem-non-susceptible commercially available strains, diluted in saline solution, is recommended as a Positive External Control. Recommended strains are:

- K. pneumoniae* bla_{KPC-2} (CCUG 59348),
- K. pneumoniae* bla_{NDM-1} (ATCC® BAA-2146™),
- K. pneumoniae* bla_{VIM-1} (NCTC 13440),
- E. coli* bla_{OXA-48} (ATCC® BAA-2523™),
- E. coli* bla_{IMP-1} (NCTC 13476).

Each of the five (5) recommended Positive External Controls may be tested alternately.

Negative External control

A standardized 0.5 McFarland suspension in saline of isolated colonies from a representative, carbapenem-non-susceptible strain of *Enterobacteriaceae*, *A. baumannii* or *P. aeruginosa* that does not carry any of the resistance genes targeted by the Revogene Carba C assay is recommended as a Negative External Control.

EXTERNAL CONTROLS PREPARATION

- Inoculate the organism onto a blood agar or a MacConkey agar plate, streak for isolation, and place a 10 µg meropenem disc in the first quadrant as a means to ensure that the isolate retains its non-susceptibility to carbapenem.
- Incubate the plate at 35 °C for 18–24 hours in aerobic conditions.
- Prepare a standardized 0.5 McFarland suspension in saline solution by collecting isolated colonies selected from the incubated agar plate close to the meropenem disc using a sterile swab.
- Adjust the suspension to achieve a turbidity equivalent to a 0.5 McFarland standard. Use either a photometric device or, if performed visually, adequate light to compare the suspension and the 0.5 McFarland standard against a card with a white background and contrasting black lines. This results in a suspension containing approximately 1 to 2 × 10⁶ Colony-Forming Units (CFU)/mL.

Process and test External Control suspensions according to SBT PREPARATION/PIE PREPARATION sections, then follow the REVOGENE SYSTEM OPERATION section.

REPEAT TESTING PROCEDURE

INDETERMINATE OR UNRESOLVED RESULT

When an Indeterminate (IND) or an Unresolved (UNR) target result is obtained for a sample, a repeat test from the corresponding inoculated SBT must be performed within the specified timeframe described in the STORAGE AND STABILITY section.

Use the original corresponding inoculated SBT to load a new PIE from the same kit lot and discard the unused SBT. Follow steps 9 to 13 of the PIE PREPARATION section for each sample and then, the REVOGENE SYSTEM OPERATION section. If needed, a new SBT can be inoculated from a new standardized 0.5 McFarland suspension.

Indeterminate or Unresolved Results for an External Control

When an Indeterminate (IND) or an Unresolved (UNR) target result is obtained for an External Control, the run should be considered invalid. A repeat test of Positive and Negative External Controls, as well as of the samples included in the same run, should be performed from the corresponding inoculated SBT. All inoculated SBT must be tested within the specified timeframe described in the STORAGE AND STABILITY section.

Use the corresponding inoculated SBT to load a new PIE from the same kit lot. Follow the PIE PREPARATION section from step 9, then follow the REVOGENE SYSTEM OPERATION section.

Unexpected Negative or Unexpected Positive Result for an External Control

When an unexpected Negative or an unexpected Positive target result is obtained for an External Control, the run should be considered invalid. A repeat test of Positive and Negative External Controls should be performed with a new set of External Controls, as described in the QUALITY CONTROL section. Handling and preparation techniques should be reviewed.

In addition, repeat testing should be performed for all samples included in the run, using the original corresponding inoculated SBT within the timeframe defined in the STORAGE AND STABILITY section with a new set of External Controls. Refer to the QUALITY CONTROL section for preparation of a new set of External controls.

Use the corresponding inoculated SBT to load a new PIE from the same kit lot. Follow the PIE PREPARATION section from step 9, then follow the REVOGENE SYSTEM OPERATION section.

RESULTS INTERPRETATION

The results are computed by the Revogene from measured fluorescent signals and embedded calculation algorithms and are available on the "Results" window. Possible reported results are listed below.

Sample type	Symbol(s) displayed on user screen	Reported result(s)	Interpretation	Additional step
Isolated Colony		Positive	KPC, NDM, VIM, OXA-48-like and IMP target DNA detected.	No additional step needed.
		Negative	KPC, NDM, VIM, OXA-48-like and IMP target DNA not detected.	No additional step needed.
		Positive /Negative	KPC and/or NDM and/or VIM and/or OXA-48-like and/or IMP target DNA detected AND Whether KPC and/or NDM and/or VIM and/or OXA-48-like and/or IMP target DNA not detected.	No additional step needed.
		Unresolved	Amplification/detection failure of the Process Control. This could be caused by inhibitory samples, microfluidic or reagent failure.	Repeat testing must be performed (refer to the REPEAT TESTING PROCEDURE section).
		Positive/ Unresolved	KPC and/or NDM and/or VIM and/or OXA-48-like and/or IMP target DNA detected AND Amplification/detection failure of the Process Control. This could be caused by inhibitory samples, microfluidic or reagent failure.	Repeat testing must be performed (refer to the REPEAT TESTING PROCEDURE section). All initial Positive results must be considered as Positive even if the repeat test reports a different result for those targets.
		Indeterminate	No reportable result due to possible Revogene detection error during the assay processing or the data analysis, or if the run is interrupted.	Repeat testing must be performed (refer to the REPEAT TESTING PROCEDURE section).

LIMITATIONS OF THE PROCEDURE

1. The Revogene Carba C assay detects *bla*_{KPC}, *bla*_{NDM}, *bla*_{VIM}, *bla*_{OXA-48-like}, and *bla*_{IMP} gene sequences, and is not for bacterial identification.
2. The performance of the Revogene Carba C assay with bacteria other than *Enterobacteriaceae*, *Acinetobacter baumannii* or *Pseudomonas aeruginosa* has not been evaluated. Organisms should be identified, and carbapenem non-susceptibility status should be determined, prior to testing with the Revogene Carba C assay.
3. *Pseudomonas aeruginosa* and *Acinetobacter baumannii* exhibit intrinsic resistance to ertapenem. Determination of the carbapenem susceptibility of these species and applicability of the Revogene Carba C assay should be based on testing with meropenem, imipenem and/or doripenem.
4. The Revogene Carba C assay is not a sub-typing tool and does not report variants of the *bla*_{KPC}, *bla*_{NDM}, *bla*_{VIM}, *bla*_{OXA-48-like}, and *bla*_{IMP} genes.
5. The Revogene Carba C assay must only be used with the Revogene instrument by trained personnel.
6. Only Revogene Carba C PIEs can be used within a single run. PIEs from other Revogene assays are not compatible with the Revogene Carba C assay and should not be used.
7. Results from the Revogene Carba C assay should be used as an adjunct to clinical observations and other information available to the physician.
8. Erroneous test results may occur from improper culture preparation, handling or storage, technical error, or sample mix-up. Careful compliance with the instructions of this insert, the Revogene Operator's Manual and to established guidelines is necessary to avoid erroneous results.
9. Contamination or unexpected positive results may occur if a PIE cap is incorrectly closed or if a droplet has been dropped on the edges of the Sample Loading Chamber.
10. Storage longer than the time period indicated in the STORAGE AND STABILITY section for the SBT samples may yield unexpected or invalid results.
11. Mutations or polymorphisms in current, new or unknown *bla*_{KPC}, *bla*_{NDM}, *bla*_{VIM}, *bla*_{OXA-48-like} and *bla*_{IMP} gene sequences may affect detection, resulting in false negative results.
12. *In silico* predictions of the inclusivity of the Revogene Carba C primers and probes were performed using sequence data available in GenBank in 2018. Analysis of new variant sequences of the targeted carbapenemase genes deposited in GenBank after 2018 has not been performed.
13. Performance of the Revogene Carba C assay with non-target carbapenemase genes, other than *bla*_{CTX-M}, *bla*_{AmpC}, *bla*_{SHV}, *bla*_{SME}, *bla*_{TEM}, and *bla*_{SPM} is unknown.

EXPECTED VALUES

A total of 512 carbapenem non-susceptible bacterial isolates were evaluated during the Carba C clinical study at three (3) sites in United States and Canada. Revogene Carba C assay results for each gene target in comparison to the Reference Method (a marketed Nucleic Acid Amplification Test [NAAT]) are presented in **Tables 1 to 4**.

CLINICAL PERFORMANCE CHARACTERISTICS

The performance of the Revogene Carba C assay was established in a clinical study conducted at three (3) sites in United States and Canada.

Bacterial isolates characterized as *Enterobacteriaceae*, *Acinetobacter baumannii* or *Pseudomonas aeruginosa* were grown on both blood agar and MacConkey agar plates. *Enterobacteriaceae* isolates were intermediate or resistant to doripenem, ertapenem, imipenem and/or meropenem according to CLSI M02¹³ Standards. *Acinetobacter baumannii* and *Pseudomonas aeruginosa* isolates were intermediate or resistant to doripenem, imipenem and/or meropenem, since they are intrinsically resistant to ertapenem. Colonies isolated on blood agar and MacConkey agar were diluted in saline to obtain a 0.5 McFarland standard equivalent suspension and tested with the Revogene Carba C assay. The Reference Method consisted of a high-performing FDA-cleared NAAT for the detection of the targeted carbapenemase genes from isolated colonies of known carbapenem non-susceptible organisms, the performance of which was established in comparison to PCR/bidirectional sequencing and which was used according to the manufacturer's instructions.

The sensitivity and specificity were calculated by comparing Revogene Carba C assay results with Reference Method. Discrepant analysis was performed on samples with discordant results between the Revogene Carba C assay and reference NAAT using an alternative PCR method for the amplification and detection of *bla*_{KPC}, *bla*_{NDM}, *bla*_{VIM}, *bla*_{OXA-48-like}, and *bla*_{IMP} genes followed by bi-directional sequencing of the amplified product.

RESULTS

A total of 532 bacterial isolates (475 clinical stock isolates and 57 prospectively collected fresh isolates) were tested during the clinical study. Twenty isolates were excluded from the analysis of performance because either they belonged to a non-target species, their identity could not be confirmed, they exhibited susceptibility to carbapenems, there was a laboratory error or valid PCR results were not obtained. Therefore, a total of 512 bacterial isolates were used to establish the performance of the Revogene Carba C assay in comparison to the Reference Method.

The sensitivity and specificity of the Revogene Carba C assay for the detection of the five (5) targeted carbapenemase genes from colonies grown on blood agar and MacConkey agar are shown in **Tables 1 and 2**, respectively.

A total of 24 isolates contained more than one carbapenemase gene as determined by the Revogene Carba C assay. The reported carbapenemase gene content of 21 of these isolates was confirmed by the Reference Method. However, the Reference Method only detected one carbapenemase gene in the other three (3) isolates.

The initial non-reportable rates (combining Unresolved and Indeterminate rates), at the media level, were 1.7% (9/516) from blood agar, and 1.4% (7/516) from MacConkey agar. At the target level, from blood agar, the rates were 1.7% (9/516) for *bla*_{KPC}, *bla*_{NDM}, *bla*_{IMP} and *bla*_{OXA-48-like} gene targets and 1.6% (8/516) for *bla*_{VIM} gene target, with an overall rate of 1.7% (9/516). From MacConkey agar, the non-reportable rates were 1.4% (7/516) across all targets. After repeat testing, the non-reportable rates were 0.8% (4/516) across all media and gene targets.

The Revogene Carba C assay performance for each organism group for blood agar and MacConkey agar are presented in **Tables 3 and 4**.

Table 1. Performance of Revogene Carba C Assay on Isolated Colonies Grown on Blood Agar

Target	N	TP	FP	TN	FN	Sensitivity (95% CI)	Specificity (95% CI)
NDM	512	186	2 ^a	322	2 ^f	98.9% (96.2 – 99.7)	99.4% (97.8 – 99.8)
KPC	512	113	3 ^b	395	1	99.1% (95.2 – 99.8)	99.2% (97.8 – 99.7)
OXA-48-like	512	65	3 ^c	444	0	100.0% (94.4 – 100.0)	99.3% (98.0 – 99.8)
IMP	512	27	19 ^d	466	0	100.0% (87.5 – 100.0)	96.1% (94.0 – 97.5)
VIM	512	52	1 ^e	459	0	100.0% (93.1 – 100.0)	99.8% (98.8 – 100.0)

N: Number; TP: True Positive; FP: False Positive; TN: True Negative; FN: False Negative; CI: Confidence Interval

Note: Discrepant testing consisted of five (5) alternative PCR followed by bi-directional sequencing and was performed for every discrepant target result. The results from discrepant testing of 21/31 samples agreed with those of the Revogene Carba C assay.

^a One (1) out of two (2) was NDM-1 positive.

^b Two (2) out of three (3) were KPC-3/KPC-38 positive.

^c One (1) out of three (3) was OXA-48 positive. Investigation suggested an OXA-48 cross-contamination at the step of sample preparation in one (1) out of three (3) isolates. Discrepant testing did not produce a sequence match with the OXA-48-like target.

^d Seventeen (17) out of 19 were found IMP positive including 1 (one) variant of IMP-4 (from Australia [2010]).

11 variants IMP-13/IMP-37 one (1) from Argentina [2006], one (1) from North America [2014], nine (9) from Europe [2005–2015], one (1) variant of IMP-27/IMP-64 (from Canada [2017]), one (1) variant of IMP-15 (from Argentina [2004]), one (1) variant of IMP-16 (from Brazil [2004]) and two (2) variants of IMP-62 (from Argentina [2006]). The discrepant analysis pointed out potential differences in IMP variant coverage between the Revogene Carba C assay and the Reference Method. Investigation suggested an IMP cross-contamination at the step of sample preparation in two (2) out of 19 isolates for which discrepant testing did not produce a sequence match with the IMP target.

^e Investigation suggested a VIM cross-contamination at the step of sample preparation. Discrepant testing did not produce a sequence match with the VIM target but produced a sequence match for the NDM target.

^f Discrepant testing did produce a sequence match with the NDM-1 target in one (1) out of two (2) isolates and produced a sequence match for the OXA-48-like target in one (1) out of two (2) isolates. The OXA-48 positive isolate was classified as FP.

Table 2. Performance of Revogene Carba C Assay on Isolated Colonies Grown on MacConkey Agar

Target	N	TP	FP	TN	FN	Sensitivity (95% CI)	Specificity (95% CI)
NDM	512	186	2 ^a	322	2 ^f	98.9% (96.2 – 99.7)	99.4% (97.8 – 99.8)
KPC	512	114	3 ^b	395	0	100.0% (96.7 – 100.0)	99.2% (97.8 – 99.7)
OXA-48-like	512	65	2 ^c	445	0	100.0% (94.4 – 100.0)	99.6% (98.4 – 99.9)
IMP	512	27	21 ^d	464	0	100.0% (87.5 – 100.0)	95.7% (93.5 – 97.2)
VIM	512	52	1 ^e	459	0	100.0% (93.1 – 100.0)	99.8% (98.8 – 100.0)

N: Number; TP: True Positive; FP: False Positive; TN: True Negative; FN: False Negative; CI: Confidence Interval

Note: Discrepant testing consisted of five (5) alternative PCR followed by bi-directional sequencing and was performed for every discrepant target result. The results from discrepant testing of 21/31 samples agreed with those of the Revogene Carba C assay.

^a One (1) out of two (2) was NDM-1 positive.^b Two (2) out of three (3) were KPC-3/KPC-38 positive.^c One (1) out of two (2) was OXA-48 positive. Investigation suggested an OXA-48 cross-contamination at the step of sample preparation in one (1) out of two (2) isolates. Discrepant testing did not produce a sequence match with the OXA-48-like target.^d Seventeen (17) out of 21 were found IMP positive including one (1) variant of IMP-4 (from Australia [2010]).

11 variants IMP-13/IMP-37 (one (1) from Argentina [2006], one (1) from North America [2014], nine (9) from Europe [2005–2015]), one (1) variant of IMP-27/IMP-64 (from Canada [2017]), one (1) variant of IMP-15 (from Argentina [2004]), one (1) variant of IMP-16 (from Brazil [2004]) and two (2) variants of IMP-62 (from Argentina [2006]). The discrepant analysis pointed out potential differences in IMP variant coverage between the Revogene Carba C assay and the Reference Method. Investigation suggested an IMP cross-contamination at the step of sample preparation in four (4) out of 21 isolates for which discrepant testing did not produce a sequence match with the IMP target.

^e Investigation suggested a VIM cross-contamination at the step of sample preparation. Discrepant testing did not produce a sequence match with the VIM target but produced a sequence match for the KPC target.^f Discrepant testing did produce a sequence match with the NDM-1 target in one (1) out of two (2) isolates and produced a sequence match for the OXA-48-like target in one (1) out of two (2) isolates. The OXA-48 positive isolate was classified as FP.

Table 3. Performance of Revogene Carba C Assay by Organism Category and by Target Gene, for Isolated Colonies Grown on Blood Agar

Medium	Organisms	Target	N	TP	FP	TN	FN	Sensitivity (95% CI)	Specificity (95% CI)
Blood agar	Enterobacteriaceae	NDM	306	85	1	219	1	98.8% (93.7 – 99.8)	99.5% (97.5 – 99.9)
		KPC	306	112	3	190	1	99.1% (95.2 – 99.8)	98.4% (95.5 – 99.5)
		OXA-48-like	306	64	3	239	0	100.0% (94.3 – 100.0)	98.8% (96.4 – 99.6)
		IMP	306	14	5	287	0	100.0% (78.5 – 100.0)	98.3% (96.1 – 99.3)
		VIM	306	12	0	294	0	100.0% (75.8 – 100.0)	100.0% (98.7 – 100.0)
	Pseudomonas aeruginosa	NDM	107	26	1	80	0	100.0% (87.1 – 100.0)	98.8% (93.3 – 99.8)
		KPC	107	0	0	107	0	N/A	100.0% (96.5 – 100.0)
		OXA-48-like	107	1	0	106	0	100.0% (20.7 – 100.0)	100.0% (96.5 – 100.0)
		IMP	107	5	14	88	0	100.0% (56.6 – 100.0)	86.3% (78.3 – 91.6)
		VIM	107	39	0	68	0	100.0% (91.0 – 100.0)	100.0% (94.7 – 100.0)
	Acinetobacter baumannii	NDM	99	75	0	23	1	98.7% (92.9 – 99.8)	100.0% (85.7 – 100.0)
		KPC	99	1	0	98	0	100.0% (20.7 – 100.0)	100.0% (96.2 – 100.0)
		OXA-48-like	99	0	0	99	0	N/A	100.0% (96.3 – 100.0)
		IMP	99	8	0	91	0	100.0% (67.6 – 100.0)	100.0% (96.0 – 100.0)
		VIM	99	1	1	97	0	100.0% (20.7 – 100.0)	99.0% (94.4 – 99.8)

N: Number; TP: True Positive; FP: False Positive; TN: True Negative; FN: False Negative; CI: Confidence Interval

More than one (1) targeted genes were detected for some isolates.

Table 4. Performance of Revogene Carba C Assay by Organism Category and by Target Gene, on Isolated Colonies Grown on MacConkey Agar

Medium	Organisms	Target	N	TP	FP	TN	FN	Sensitivity (95% CI)	Specificity (95% CI)
MacConkey agar	Enterobacteriaceae	NDM	306	85	1	219	1	98.8% (93.7 – 99.8)	99.5% (97.5 – 99.9)
		KPC	306	113	3	190	0	100.0% (96.7 – 100.0)	98.4% (95.5 – 99.5)
		OXA-48-like	306	64	2	240	0	100.0% (94.3 – 100.0)	99.2% (97.0 – 99.8)
		IMP	306	14	5	287	0	100.0% (78.5 – 100.0)	98.3% (96.1 – 99.3)
		VIM	306	12	1	293	0	100.0% (75.8 – 100.0)	99.7% (98.1 – 99.9)
	Pseudomonas aeruginosa	NDM	107	26	1	80	0	100.0% (87.1 – 100.0)	98.8% (93.3 – 99.8)
		KPC	107	0	0	107	0	N/A	100.0% (96.5 – 100.0)
		OXA-48-like	107	1	0	106	0	100.0% (20.7 – 100.0)	100.0% (96.5 – 100.0)
		IMP	107	5	14	88	0	100.0% (56.6 – 100.0)	86.3% (78.3 – 91.6)
		VIM	107	39	0	68	0	100.0% (91.0 – 100.0)	100.0% (94.7 – 100.0)
	Acinetobacter baumannii	NDM	99	75	0	23	1	98.7% (92.9 – 99.8)	100.0% (85.7 – 100.0)
		KPC	99	1	0	98	0	100.0% (20.7 – 100.0)	100.0% (96.2 – 100.0)
		OXA-48-like	99	0	0	99	0	N/A	100.0% (96.3 – 100.0)
		IMP	99	8	2	89	0	100.0% (67.6 – 100.0)	97.8% (92.3 – 99.4)
		VIM	99	1	0	98	0	100.0% (20.7 – 100.0)	100.0% (96.2 – 100.0)

N: Number; TP: True Positive; FP: False Positive; TN: True Negative; FN: False Negative; CI: Confidence Interval

More than one (1) targeted genes were detected for some isolates.

ANALYTICAL PERFORMANCE CHARACTERISTICS
INCLUSIVITY

Inclusivity of the Revogene Carba C assay was determined for 58 carbapenem non-susceptible isolates of *Enterobacteriaceae*, *Pseudomonas aeruginosa* and *Acinetobacter baumannii* from various geographic and temporal origins including:

- Two (2) strains with two (2) resistance genes;
- Eleven (11) strains including five (5) variants with the *bla_{KPC}* resistance gene;
- Eleven (11) strains including at least three (3) variants with the *bla_{KPC}* resistance gene;
- Fourteen (14) strains including five (5) variants with the *bla_{NDM}* resistance gene;
- Eleven (11) strains including three (3) variants with the *bla_{OXA-48-like}* resistance gene;
- Nine (9) strains including four (4) variants with the *bla_{VIM}* resistance gene.

Each strain was tested from a standardized 0.5 McFarland bacterial suspension corresponding to 1.5 to 3x10⁶ CFU/mL of SB. Three (3) replicates per strain were tested using three (3) different Revogene Carba C kit lots (1 replicate per kit lot). The 58 strains were detected by the Revogene Carba C assay and are described in Table 5.

Table 5. Carbapenem Non-susceptible Strains Tested for Inclusivity with the Revogene Carba C Assay

Species	Collection Number	Resistance Gene and Variant
<i>Pseudomonas aeruginosa</i>	ATCC® BAA-2793™	KPC-2 VIM-2
<i>Klebsiella pneumoniae</i>	CCRI-23061	OXA-232 NDM-1
<i>Klebsiella pneumoniae</i>	NCTC 13438	KPC-3
<i>Klebsiella pneumoniae</i>	ATCC® BAA-1705™	KPC-2
<i>Pseudomonas aeruginosa</i>	CCRI-21587	KPC-2
<i>Enterobacter cloacae</i>	CCRI-21578	KPC-4
<i>Klebsiella oxytoca</i>	CCRI-21581	KPC-3
<i>Klebsiella pneumoniae</i>	CCRI-19587	KPC-3
<i>Klebsiella pneumoniae</i>	CCRI-19570	KPC-2
<i>Klebsiella pneumoniae</i>	CCUG 59413	KPC-3
<i>Klebsiella pneumoniae</i>	CCUG 59348	KPC-2
<i>Klebsiella pneumoniae</i>	CCUG 56233	KPC-2
<i>Escherichia coli</i>	ATCC® BAA-2340™	KPC ^a
<i>Klebsiella pneumoniae</i>	NCTC 13443	NDM-1
<i>Klebsiella pneumoniae</i>	ATCC® BAA-2146™	NDM-1
<i>Escherichia coli</i>	CCRI-22255	NDM-1
<i>Klebsiella pneumoniae</i>	CCRI-21711	NDM-1
<i>Klebsiella pneumoniae</i>	CCRI-22199	NDM-1
<i>Providencia rettgeri</i>	CCRI-22257	NDM-1
<i>Providencia stuartii</i>	CCRI-22256	NDM-1
<i>Klebsiella pneumoniae</i>	CCRI-22254	NDM-4
<i>Escherichia coli</i>	CCRI-23064	NDM-5
<i>Escherichia coli</i>	CCRI-23464	NDM-5
<i>Escherichia coli</i>	CCRI-23065	NDM-6
<i>Escherichia coli</i>	CCRI-23066	NDM-7
<i>Enterobacter cloacae</i>	ATCC® BAA-2468™	NDM-1
<i>Klebsiella pneumoniae</i>	CCUG 60138	NDM-1
<i>Pseudomonas aeruginosa</i>	NCTC 13437	VIM-10
<i>Klebsiella pneumoniae</i>	NCTC 13439	VIM-1
<i>Klebsiella pneumoniae</i>	NCTC 13440	VIM-1
<i>Klebsiella pneumoniae</i>	CCRI-19585	VIM-1
<i>Klebsiella pneumoniae</i>	CCRI-22258	VIM-1
<i>Pseudomonas aeruginosa</i>	CCRI-21588	VIM-2
<i>Serratia marcescens</i>	CCRI-22261	VIM-2
<i>Pseudomonas aeruginosa</i>	CCRI-22720	VIM-2
<i>Klebsiella pneumoniae</i>	CCRI-22259	VIM-19
<i>Klebsiella pneumoniae</i>	NCTC 13442	OXA-48
<i>Klebsiella pneumoniae</i>	CCRI-22263	OXA-48
<i>Escherichia coli</i>	CCRI-22265	OXA-48
<i>Enterobacter cloacae</i>	CCRI-22266	OXA-48
<i>Klebsiella pneumoniae</i>	CCRI-22264	OXA-181
<i>Providencia rettgeri</i>	CCRI-22267	OXA-181
<i>Klebsiella pneumoniae</i>	CCRI-23060	OXA-204
<i>Citrobacter freundii</i>	CCRI-23374	OXA-204
<i>Escherichia coli</i>	ATCC® BAA-2523™	OXA-48
<i>Klebsiella pneumoniae</i>	ATCC® BAA-2524™	OXA-48
<i>Klebsiella pneumoniae</i>	CCUG 64452	OXA-48
<i>Escherichia coli</i>	NCTC 13476	IMP-1
<i>Acinetobacter baumannii</i>	CCRI-19488	IMP-1
<i>Klebsiella pneumoniae</i>	CCRI-19569	IMP-1
<i>Klebsiella pneumoniae</i>	CCRI-19582	IMP-1
<i>Pseudomonas aeruginosa</i>	CCRI-21589	IMP-1
<i>Klebsiella pneumoniae</i>	CCRI-19583	IMP-4
<i>Klebsiella pneumoniae</i>	CCRI-19588	IMP-4
<i>Citrobacter youngae</i>	CCRI-21591	IMP-4
<i>Klebsiella pneumoniae</i>	CCRI-19584	IMP-8
<i>Pseudomonas aeruginosa</i>	CCRI-21590	IMP-9
<i>Serratia marcescens</i>	CCRI-22262	IMP-11

^a Identification of the resistance gene variant is not available.

In addition, an *in silico* analysis was performed on November 7th, 2018 to assess inclusivity of the primers and probes of the Revogene Carba C assay targets. For each target, all the *bla_{KPC}*, *bla_{NDM}*, *bla_{VIM}*, *bla_{OXA-48-like}* and *bla_{IMP}* resistance gene sequences listed in the National Center for Biotechnology Information (NCBI) database were analyzed. One (1) representative sequence of each known variant was aligned. The number of variants analyzed and predicted to be detected is summarized for each target in Table 6.

Table 6. Summary of Targeted Resistance Gene Variants Detection Based on the *in silico* Prediction.

Target	Tested with Revogene® Carba C			<i>in silico</i> Prediction		
	Nb. of Samples	Variants Detected	Variants not Detected	Detectable ^a	Potentially Detectable ^b	Not Detectable
KPC	12	2, 3, 4	None	2 to 38	None	N/A
NDM	15	1, 4, 5, 6, 7	None	1 to 24	None	N/A
IMP	11	1, 4, 8, 9, 11	None	1, 2, 4 to 6, 8 to 10, 13 to 20, 23 to 26, 28 to 30, 32, 33, 37, 38, 40, 42, 45, 47 to 49, 53 to 56, 59, 60, 62, 66, 69 to 72, 74 to 79	3, 7 ^c , 11 ^d , 21, 22, 27, 34, 41, 43, 44, 51, 52, 58, 61, 64, 67, 68, 73	12, 31, 35, 63
OXA-48-like	12	48, 181, 204, 232	None	48, 162, 181, 199, 204, 232, 244, 245, 252 ^e , 370, 484, 505, 514 ^f , 515 ^e , 519, 546 ^e , 547 ^e , 566	None	54 ^e , 163 ⁱ , 247 ^j , 405 ⁱ , 416 ^e , 436, 438 ⁱ , 439 ⁱ , 517, 535 ^e , 538 ^e , 567
VIM	10	1, 2, 10, 19	None	1 to 6, 8 to 12, 14 to 20, 23 to 46, 48 to 50, 52 to 55, 57, 59, 60	51, 56, 58	7, 13, 47

^a Based on alignments with identity and query cover ≥ 95% and E-values < 0.01.^b Based on alignments presenting not more than two (2) nucleotides mismatches.^c A recombinant strain carrying an IMP-7 gene was tested with the Revogene Carba C assay, in complement to the *in silico* study, and its detection was confirmed.^d A clinical strain carrying an IMP-11 gene was tested with the Carba C assay in the analytical inclusivity study and its detection was confirmed.^e Variants only identified in the rare opportunistic human pathogens of the *Shewanella* genus.^f Variants presenting a deletion in their sequence resulting in no carbapenemase activity.**CROSS-REACTIVITY**

The cross-reactivity of the Revogene Carba C assay was assessed with high loads of carbapenem-non-susceptible strains carrying different β-lactamase resistance genes that are not targeted by the assay.

The study included 50 *Enterobacteriaceae*, *Acinetobacter baumannii* and *Pseudomonas aeruginosa* strains carrying different β-lactamase resistance genes (i.e., *bla*_{CTX-M}, *bla*_{AMP}, *bla*_{SHV}, *bla*_{SME}, *bla*_{TEM}, and *bla*_{SPM}) (**Table 7**). Strains were tested at a minimal load of 8.18x10⁶ CFU/mL of SB. Three (3) replicates per strain were tested using three (3) different Revogene Carba C kit lots (1 replicate per kit lot).Under the conditions of the study, 49 strains were found non-reactive with the Revogene Carba C assay. *Klebsiella pneumoniae* CCUG 59359 *bla*_{TEM-52} strain was found reactive at a final concentration of 1.13x10⁷ CFU/mL of SB but was found non-reactive when diluted 10-fold at a final concentration of 1.13x10⁶ CFU/mL of SB.The cross-reactivity with primers and probes of the Revogene Carba C assay was evaluated by an *in silico* analysis performed on sequences contained in the National Center for Biotechnology Information (NCBI) database on September 0th, 2018. No other resistance gene than the assay targets (i.e. *bla*_{KPC}, *bla*_{NDM}, *bla*_{VIM}, *bla*_{OXA-48-like} and *bla*_{IMP} resistance genes) was found to have a homology with the primers and probes of the Revogene Carba C assay, including *bla*_{MRI} and *bla*_{VEB}, which are two (2) other β-lactamase resistance genes that are not targeted by the assay.**Table 7. List of Strains Tested for Cross-Reactivity with the Revogene Carba C Assay.**

Species	Collection Number	Non-Targeted Resistance Gene(s) ¹
<i>Enterobacter cloacae</i>	CCRI-22353	ACT-15
<i>Enterobacter cloacae</i>	CCRI-23473	No β-lactamase Gene Identified
<i>Enterobacter cloacae</i>	CCRI-21540	ACT-7
<i>Enterobacter cloacae</i>	CCRI-22075	ACT-7
<i>Enterobacter cloacae</i>	CCRI-22097	ACT-16
<i>Enterobacter cloacae</i>	CCRI-23318	TEM-206, CTX-M-15, CMH-1
<i>Escherichia coli</i>	CCRI-21970	AmpC1, AmpC2, MrdA, AmpH, CMY-44
<i>Klebsiella aerogenes</i>	CCRI-19495	SHV-5, AmpC
<i>Serratia marcescens</i>	CCRI-21537	SRT-1
<i>Serratia marcescens</i>	CCRI-23334	SME-4, SRT-1
<i>Enterobacter cloacae</i>	CCRI-21536	ACT-5
<i>Enterobacter cloacae</i>	CCRI-21603	ACT-7
<i>Enterobacter cloacae</i>	CCRI-21692	ACT-14
<i>Acinetobacter baumannii</i>	CCRI-1016	TEM-90, Mbl, BlaA2, OXA-65 ³ , Zn-dependent hydrolase
<i>Klebsiella pneumoniae</i>	CCRI-1015	TEM-171, SCO-1, PER-2, OXA-9 ³ , SHV-39, AmpH
<i>Acinetobacter baumannii</i>	CCRI-1017	TEM-206, SCO-1, Mbl, BlaA2, OXA-67 ³ , Zn-dependent hydrolase
<i>Proteus mirabilis</i>	CCRI-831	TEM-206, CTX-M-2, OXA-2 ³
<i>Proteus mirabilis</i>	CCRI-825	TEM-33, CTX-M-2, OXA-2 ³
<i>Pseudomonas aeruginosa</i>	CCRI-873	OXA-50 ³
<i>Pseudomonas aeruginosa</i>	CCRI-1228	OXA-50 ³
<i>Salmonella enterica</i>	CCRI-8892	TEM-166, CTX-M-5, OXA-1 ³
<i>Salmonella enterica</i>	CCRI-8893	TEM-95, CTX-M-5, OXA-1 ³
<i>Escherichia coli</i>	CCRI-785	TEM-206, AmpC1, MrdA, AmpC2, AmpH
<i>Escherichia coli</i>	CCRI-779	TEM-206, AmpC1, MrdA, AmpC2, AmpH
<i>Escherichia coli</i>	CCRI-778	AmpC2, MrdA
<i>Enterobacter cloacae</i>	CCRI-3854	ACT-4 ²
<i>Enterobacter aerogenes</i>	CCRI-3879	AmpC
<i>Enterobacter aerogenes</i>	CCRI-3853	AmpC
<i>Enterobacter cloacae</i>	CCRI-3852	ACT-7
<i>Klebsiella quasipneumoniae</i>	CCRI-806	OKP-B-11
<i>Klebsiella pneumoniae</i>	CCRI-784	SHV-27, AmpH
<i>Escherichia coli</i>	CCRI-878	AmpC2, MrdA, AmpH
<i>Proteus mirabilis</i>	CCRI-826	TEM-215
<i>Klebsiella pneumoniae</i>	CCUG 54718	TEM-33, CTX-M-15, OXA-1 ³ , AmpH
<i>Klebsiella pneumoniae</i>	CCUG 59358	SHV-14, OXA-1 ³ , LAP-2
<i>Escherichia coli</i>	NCTC 13441	CTX-M-15, TEM-198, MrdA, OXA-1 ³ , AmpC2
<i>Klebsiella pneumoniae</i>	CCUG 58546	SHV-44, AmpH
<i>Klebsiella pneumoniae</i>	CCUG 59349	CTX-M-15, OXA-1 ³ , TEM-105, SHV-11, AmpH
<i>Klebsiella pneumoniae</i>	CCUG 59359	SHV-70, TEM-15, AmpH
<i>Klebsiella pneumoniae</i>	CCUG 59360	TEM-168, SHV-12, OXA-9 ³ , AmpH
<i>Escherichia coli</i>	CCUG 55970	CTX-M-9, AmpC2, TEM-206, MrdA, AmpC1, AmpH
<i>Escherichia coli</i>	CCUG 55971	TEM-143, CTX-M-15, AmpC2
<i>Escherichia coli</i>	CCUG 55972	AmpC1, CTX-M-2, AmpC2, AmpH
<i>Escherichia coli</i>	CCUG 58540	AmpC2, TEM-206, CTX-M-15, OXA-1 ³ , MrdA, AmpH
<i>Escherichia coli</i>	CCUG 58541	CTX-M-14, TEM-104, MrdA, AmpC2, AmpH
<i>Escherichia coli</i>	CCUG 58542	CTX-M-15, AmpC2, OXA-1 ³ , MrdA
<i>Proteus mirabilis</i>	CCRI-21789	No β-lactamase Gene Identified
<i>Klebsiella pneumoniae</i>	NCTC 13465	SHV-85, TEM-206, AmpH
<i>Escherichia coli</i>	CCRI-21710	AmpC2, MrdA, CTX-M-15, AmpH, OXA-1 ³
<i>Pseudomonas aeruginosa</i>	C72	SPM

¹ Analyzed by Whole Genome Sequencing, except for *P. aeruginosa* C72 strain.² Resistance gene identified with 99.9% homology.³ These *bla*_{OXA} variants are not part of the *bla*_{OXA-48-like} family, but are part of the Ambler class D.

INTERFERING SUBSTANCES

The potentially inhibitory effect of 12 combinations of agar plates and sterile saline solutions that may be used for isolation and 0.5 McFarland (McF) suspension preparation was assessed using five (5) positive bacterial strains, each harboring one (1) of the resistance genes detected by the Revogene Carba C assay, and one (1) negative carbapenem-non-susceptible strain that does not carry any of the resistance genes targeted by the assay. Each strain was first cultured on Blood Agar Plates (BAP) and MacConkey Agar Plates (MAP) in presence of a 10 µg meropenem disc, and then sub-cultured on fresh BAP or MAP. Isolated colonies obtained from the fresh BAP or MAP were sampled using a sterile dry swab and diluted in a sterile saline solution at a concentration equivalent to a standardized 0.5 McFarland suspension. A total of 12 combinations were tested (**Table 8**). For each combination, one (1) replicate of each strain was tested with three (3) Revogene Carba C kit lots. None of the 12 combinations showed reportable interference with the Revogene Carba C assay.

Table 8. List of Combinations of Agar Plates and Saline Solutions Tested with the Revogene Carba C Assay.

Agar Plate / Saline Solution (Manufacturer)
Colombia Blood Agar 5% (Hardy Diagnostics) / BBL™ Prepared Saline Solution (BD)
Colombia Blood Agar 5% (Hardy Diagnostics) / Saline 0.85% (Thermo Scientific™ Oxoid™)
Prepared Media BD BBL™ Columbia Agar with 5% Sheep Blood (BD) / BBL™ Prepared Saline Solution (BD)
Prepared Media BD BBL™ Columbia Agar with 5% Sheep Blood (BD) / Saline 0.85% (Thermo Scientific™ Oxoid™)
Columbia Agar with 5% Sheep Blood (bioMérieux) / BBL™ Prepared Saline Solution (BD)
Columbia Agar with 5% Sheep Blood (bioMérieux) / Saline 0.85% (Thermo Scientific™ Oxoid™)
MacConkey Agar (Hardy Diagnostics) / BBL™ Prepared Saline Solution (BD)
MacConkey Agar (Hardy Diagnostics) / Saline 0.85% (Thermo Scientific™ Oxoid™)
MacConkey Agar Medium (Thermo Scientific™ Remel™) / BBL™ Prepared Saline Solution (BD)
MacConkey Agar Medium (Thermo Scientific™ Remel™) / Saline 0.85% (Thermo Scientific™ Oxoid™)
MacConkey Agar (Thermo Scientific™ Oxoid™) / BBL™ Prepared Saline Solution (BD)
MacConkey Agar (Thermo Scientific™ Oxoid™) / Saline 0.85% (Thermo Scientific™ Oxoid™)

CARRY-OVER AND CROSS-CONTAMINATION

The carry-over (between-run) and cross-contamination (within-run) were assessed using high positive and negative samples. The concentration of each bacterial suspension was standardized to 4 McFarland ($\geq 1.14 \times 10^7$ CFU/mL of Sample Buffer), higher than the nominal concentration of 0.5 McFarland that is specified for use in the assay. Positive samples were prepared with the *Klebsiella pneumoniae* CCUG 59348 strain which harbors the *bla_{KPC}* gene. Negative samples were prepared with the carbapenem-non-susceptible *Enterobacter cloacae* CCR1-22760 strain that does not carry any of the carbapenemase resistance genes targeted by the assay.

For the carry-over study, a run of eight (8) replicates of high positive samples followed by a run of eight (8) replicates of negative samples was performed by two (2) operators with the Revogene Carba C assay, for a total of ten (10) runs on one (1) Revogene.

For the cross-contamination study, four (4) high positive samples and four (4) negative samples were tested by alternating sample preparation between positive and negative samples for each run. A total of ten (10) runs were performed with the Revogene Carba C assay by two (2) operators on one (1) Revogene.

No carry-over or cross-contamination of *bla_{KPC}* was observed.

REPRODUCIBILITY AND PRECISION

A study was conducted to assess the reproducibility and precision of the Revogene Carba C assay, including the between-laboratory reproducibility, the between-lot reproducibility and the within-laboratory precision study.

Between-laboratory reproducibility study was performed on three (3) sites, by two (2) operators per site, over five (5) distinct days using one (1) Revogene Carba C kit lot. A total of 120 replicates for each of the ten (10) positive bacterial carbapenem-non-susceptible strains, each one harboring one (1) resistance gene targeted by the Revogene Carba C assay (n=2 strains per resistance gene), were tested. A total of 240 replicates of the negative bacterial carbapenem-non-susceptible strain that does not carry any of the resistance genes targeted by the assay were tested.

Between-lot reproducibility study was performed on one (1) site, by two (2) operators, over 15 days using three (3) Revogene Carba C kit lots (5 days per kit lot). A total of 120 replicates for each of the ten (10) positive bacterial carbapenem-non-susceptible strains, each one harboring one (1) resistance gene targeted by the Revogene Carba C assay (n=2 strains per resistance gene), were tested. A total of 240 replicates of the negative bacterial carbapenem-non-susceptible strain that does not carry any of the resistance genes targeted by the assay were tested.

The within-laboratory precision was established with data obtained with kit lot #1 tested by two (2) operators exclusively on site #1 over a total of five (5) days. A total of 40 replicates for each of the ten (10) positive bacterial carbapenem-non-susceptible strains, each one harboring one (1) resistance gene targeted by the Revogene Carba C assay (n=2 strains per resistance gene), were considered for the analysis. A total of 80 replicates of the negative bacterial carbapenem-non-susceptible strain that does not carry any of the resistance genes targeted by the assay were also considered.

All strains were tested from standardized 0.5 McFarland (McF) suspensions.

For the between-laboratory reproducibility, the overall percent agreement was 100% (120/120) for each positive strain and was 98.8% (237/240) for the negative strain (**Table 9**). For the between-lot reproducibility, the overall percent agreement was 100% (120/120) for each positive strain and was 99.6% (239/240) for the negative strain (**Table 10**). Mean Cycle Threshold (Ct) values with variance components of Standard Deviation (SD) and Coefficient of Variation (CV) are also shown in **Tables 9 and 10**.

For the within-laboratory precision, the percent agreement was 100% (40/40) for each positive strain, and was 100% (80/80) for the negative strain (**Table 11**).

Table 9. Between-Laboratory Reproducibility Study Results Using One (1) Revogene Carba C Assay Kit Lot

Species	Resistance Gene and Variant	Site 1		Site 2		Site 3		Overall Results/Total	Overall Percent Agreement ¹	Overall 95% CI	Ct Values ²		
		Results/Total	Percent Agreement	Results/Total	Percent Agreement	Results/Total	Percent Agreement				Overall Mean	SD	CV%
<i>Escherichia coli</i>	IMP-1	40/40	100%	40/40	100%	40/40	100%	120/120	100%	96.9%-100.0%	29.7	1.6	5.5
<i>Klebsiella pneumoniae</i>	IMP-4	40/40	100%	40/40	100%	40/40	100%	120/120	100%	96.9%-100.0%	29.6	1.3	4.3
<i>Enterobacter cloacae</i>	KPC-4	40/40	100%	40/40	100%	40/40	100%	120/120	100%	96.9%-100.0%	30.5	2.1	7.0
<i>Klebsiella pneumoniae</i>	KPC-2	40/40	100%	40/40	100%	40/40	100%	120/120	100%	96.9%-100.0%	30.4	2.2	7.3
<i>Klebsiella pneumoniae</i>	NDM-1	40/40	100%	40/40	100%	40/40	100%	120/120	100%	96.9%-100.0%	30.3	1.0	3.1
<i>Escherichia coli</i>	NDM-5	40/40	100%	40/40	100%	40/40	100%	120/120	100%	96.9%-100.0%	30.2	1.5	4.8
<i>Klebsiella pneumoniae</i>	OXA-181	40/40	100%	40/40	100%	40/40	100%	120/120	100%	96.9%-100.0%	28.5	1.4	4.8
<i>Escherichia coli</i>	OXA-48	40/40	100%	40/40	100%	40/40	100%	120/120	100%	96.9%-100.0%	30.5	1.9	6.1
<i>Klebsiella pneumoniae</i>	VIM-1	40/40	100%	40/40	100%	40/40	100%	120/120	100%	96.9%-100.0%	27.6	1.0	3.5
<i>Klebsiella pneumoniae</i>	VIM-1	40/40	100%	40/40	100%	40/40	100%	120/120	100%	96.9%-100.0%	28.5	1.2	4.3
<i>Enterobacter cloacae</i>	Negative	80/80	100%	79/80	98.8%	78/80	97.5%	237/240	98.8%	96.4%-99.6%	34.9	1.9	5.5

¹ For the negative strain, percent agreement was calculated from the negative results.² For resistance genes and variants, Ct values reported are for the specified gene. For the negative strain, Ct values reported are for the PrC.

Table 10. Between-Lot Reproducibility Study Results on One (1) Site Using Three (3) Revogene Carba C Assay Kit Lots

Species	Resistance Gene and Variant	Lot 1		Lot 2		Lot 3		Overall Results/Total	Overall Percent Agreement ¹	Overall 95% CI	Ct Values ²		
		Results/Total	Percent Agreement	Results/Total	Percent Agreement	Results/Total	Percent Agreement				Overall Mean	SD	CV%
<i>Escherichia coli</i>	IMP-1	40/40	100%	40/40	100%	40/40	100%	120/120	100%	96.9%-100.0%	29.0	1.1	3.8
<i>Klebsiella pneumoniae</i>	IMP-4	40/40	100%	40/40	100%	40/40	100%	120/120	100%	96.9%-100.0%	29.2	1.4	4.7
<i>Enterobacter cloacae</i>	KPC-4	40/40	100%	40/40	100%	40/40	100%	120/120	100%	96.9%-100.0%	29.6	2.2	7.5
<i>Klebsiella pneumoniae</i>	KPC-2	40/40	100%	40/40	100%	40/40	100%	120/120	100%	96.9%-100.0%	30.8	2.3	7.4
<i>Klebsiella pneumoniae</i>	NDM-1	40/40	100%	40/40	100%	40/40	100%	120/120	100%	96.9%-100.0%	29.5	1.4	4.8
<i>Escherichia coli</i>	NDM-5	40/40	100%	40/40	100%	40/40	100%	120/120	100%	96.9%-100.0%	28.9	1.5	5.1
<i>Klebsiella pneumoniae</i>	OXA-181	40/40	100%	40/40	100%	40/40	100%	120/120	100%	96.9%-100.0%	28.2	1.5	5.2
<i>Escherichia coli</i>	OXA-48	40/40	100%	40/40	100%	40/40	100%	120/120	100%	96.9%-100.0%	28.6	1.6	5.7
<i>Klebsiella pneumoniae</i>	VIM-1	40/40	100%	40/40	100%	40/40	100%	120/120	100%	96.9%-100.0%	27.5	1.0	3.6
<i>Klebsiella pneumoniae</i>	VIM-1	40/40	100%	40/40	100%	40/40	100%	120/120	100%	96.9%-100.0%	28.1	1.4	4.8
<i>Enterobacter cloacae</i>	Negative	80/80	100%	79/80	98.8%	80/80	100%	239/240	99.6%	97.7%-99.9%	33.0	2.9	8.8

¹ For the negative strain, percent agreement was calculated from the negative results.² For resistance genes and variants, Ct values reported are for the specified gene. For the negative strain, Ct values reported are for the PrC.

Table 11. Within-Laboratory Precision Study Results on One (1) Site Using One (1) Revogene Carba C Assay Kit Lot

Species	Resistance Gene and Variant	Overall Percent Agreement ¹	Overall 95% CI
<i>Escherichia coli</i>	IMP-1	100% (40/40)	91.2%-100.0%
<i>Klebsiella pneumoniae</i>	IMP-4	100% (40/40)	91.2%-100.0%
<i>Enterobacter cloacae</i>	KPC-4	100% (40/40)	91.2%-100.0%
<i>Klebsiella pneumoniae</i>	KPC-2	100% (40/40)	91.2%-100.0%
<i>Klebsiella pneumoniae</i>	NDM-1	100% (40/40)	91.2%-100.0%
<i>Escherichia coli</i>	NDM-5	100% (40/40)	91.2%-100.0%
<i>Klebsiella pneumoniae</i>	OXA-181	100% (40/40)	91.2%-100.0%
<i>Escherichia coli</i>	OXA-48	100% (40/40)	91.2%-100.0%
<i>Klebsiella pneumoniae</i>	VIM-1	100% (40/40)	91.2%-100.0%
<i>Klebsiella pneumoniae</i>	VIM-1	100% (40/40)	91.2%-100.0%
<i>Enterobacter cloacae</i>	Negative	100% (80/80)	95.4%-100.0%

¹ For the negative strain, percent agreement was calculated from the negative results.**E-LABELING**Documentation related to this product can be accessed online at www.meridianbioscience.com/pi. Additionally, paper copies are available upon request by contacting your local distributor or via the phone number listed on the kit box.

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REFERENCES

- He, Q, Gallagher J. Pharmacodynamic and pharmacokinetic considerations in the treatment of critically ill patients infected with carbapenem-resistant *Enterobacteriaceae*. *Virulence*. 2017;8(4):440-452.
- Mollenkopf, Stull et al., Antimicrobial Agents and Chemotherapy, 2016.
- <https://www.cdc.gov/hai/organisms/cre/cre-clinicianfaq.html>
- Neuner, E, Gallagher J. Pharmacodynamic and pharmacokinetic considerations in the treatment of critically ill patients infected with carbapenem-resistant *Enterobacteriaceae*. *Virulence*. 2017;8(4):440-452.
- Gniadek, TJ, et al. Carbapenem-resistant non-glucose- fermenting Gram-negative bacilli: the missing piece to the puzzle. *J Clin microbiol*. 2016;54(7):1700-1710.
- Kaye, KS & Pogue JM. Infections caused by resistant Gram-negative bacteria: epidemiology and management. *Pharmacotherapy* 2015;35:949-962.
- Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health. Biosafety in microbiological and biomedical laboratories-5th ed. HHS Publication No. (CDC) 21-1112. Revised Dec 2009
- CLSI M29-A4, 2014, Protection of laboratory workers from occupationally acquired infections; Approved Guideline-4th ed.
- <https://www.cdc.gov/nndss/conditions/notifiable/2019/infectious-diseases/>
- CIFOR Analysis of State Legal Authorities, available at <http://www.cifor.us/>
- CLSI M100 Performance Standards for Antimicrobial Susceptibility Testing. 2018. 28th ed.
- SN134822 Revogene® Operator's Manual.
- CLSI M02 Performance Standards for Antimicrobial Disk Susceptibility Tests. 2018. 13th ed.

Meridian Bioscience Europe s.a./n.v.
2 Avenue du Japon - 1420 Braine l'Alleud
BELGIUM
Tel: +32 (0) 67 89 59 59
Fax: +32 (0) 67 89 59 58
Email: info.bn1@meridianbioscience.eu

Meridian Bioscience Europe France
34 rue de Ponthieu - 75008 Paris
FRANCE
Tel: +33 (0) 1 42 56 04 40
Fax: +33 (0) 9 70 06 62 10
Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe B.V.
Zekeringstraat 17 A
1014BM Amsterdam
NETHERLANDS
Tel: +31 (0) 411 62 11 66
Fax: +31 (0) 411 62 48 41
Email: Info.bn1@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labelling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de símbolos, Zeichenerklärung)

	Use-by date / Data di scadenza / Date de péremption / Fecha de caducidad / Verwendbar bis
LOT	Batch code / Codice di lotto / Code de lot / Código de lote / Chargennummer
IVD	In vitro diagnostic medical device / Dispositivo medico diagnostico <i>in vitro</i> / Instrument de test diagnostique <i>in vitro</i> / Producto sanitario para diagnóstico <i>In vitro</i> / <i>In-vitro-Diagnostikum</i>
	CE Mark / Marcatura CE / Symbole CE / Marcado CE / CE-Kennzeichen
REF	Catalog number / Numero di catalogo / Référence catalogue / Referencia / Bestellnummer
	Consult instructions for use / Consultare le istruzioni per l'uso / Consulter le mode d'emploi / Consultar las instrucciones de uso / Gebrauchsanweisung beachten
	Manufacturer / Produttore / Fabricant / Fabricante / Hersteller
	Contains sufficient for <n> tests / Contiene una quantità sufficiente per <n> test / Contient le matériel suffisant pour <n> tests / Contiene la cantidad suficiente para <n> ensayos / Inhalt ausreichend für <n> Tests
	Temperature limit / Limite di temperatura / Limite de température / Límite de temperatura / Temperaturgrenze
EC REP	Authorized representative in the European Community / Rappresentante autorizzato nella Comunità Europea / Représentant agréé dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter EU-Repräsentant
	Do not reuse / Non riutilizzare / Ne pas réutiliser / No reutilizar / Nicht wiederverwenden
	Keep dry / Conservare all'asciutto / Conserver au sec / Mantener seco / Vor Feuchtigkeit schützen
	Contains # pouches: 1 Disposable Transfer Tool (DTT), 1 Sample Buffer Tube (SBT), 1 PIE / Contiene # buste: 1 PIE, 1 Sample Buffer Tube (SBT), 1 Disposable Transfer Tool (DTT) / Contient # sachets: Une cartouche PIE, Un tube échantillon, Un outil de transfert jetable (OTJ) / Incluye # bolsitas: 1 PIE, 1 Sample Buffer Tube (SBT), 1 Disposable Transfer Tool (DTT) / Enthält # Beutel: 1 PIE, 1 Sample Buffer Tube (SBT), 1 Disposable Transfer Tool (DTT)
	Humidity Limitation / Limitazione dell'umidità / Limite d'humidité / Limitación de humedad / Feuchtigkeitsbegrenzung
Rx Only	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / Ne pas utiliser si l'emballage est endommagé / No utilizar si el envase está dañado / Nicht verwenden, wenn die Verpackung beschädigt ist
	Keep away from sunlight / Tenere lontano dalla luce del sole / Tenir à l'écart de la lumière du soleil / Mantener alejado de la luz solar / Vor Sonnenlicht schützen
24x	Contains 24 Disposable Transfer Loops (DTL) / Contiene 24 loop di trasferimento monouso (DTL) / Contient 24 boucles de transfert jetables (DTL) / Contiene 24 bucles de transferencia desechables (DTL) / Enthält 24 Einwegübertragungsschleifen (DTL)
DTT	Disposable Transfer Tool (DTT) / Disposable Transfer Tool 9DTT) / Un outil de transfert jetable (OTJ) / Disposable Transfer Tool (DTT) / Disposable Transfer Tool (DTT)
PIE	Revogene Test Device / Dispositivo Test Revogene / Dispositif de test Revogene / Dispositivo para la Prueba Revogene / Revogene-Analysegefäß
EUA	For Emergency Use Authorization only / Solo per l'autorizzazione per l'uso di emergenza / pour autorisation d'utilisation d'urgence uniquement / para autorización de uso de emergencia solamente / nur für Notfallverwendungsaufzierung
LOOP	Disposable Transfer Loop / Loop di trasferimento monouso / Boucles de transfert jetables / Bucle de transferencia desechables / Einwegübertragungsschleifen
SBT	Sample Buffer Tube / Sample Buffer Tube / Un tube échantillon / Smaple Buffer Tube / Sample Buffer Tube
	Revogene Test Device / Dispositivo Test Revogene / Dispositif de test Revogene / Dispositivo para la Prueba Revogene / Revogene-Analysegefäß

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order call Customer Service Department at 800-543-1980.

TaqMan is a registered trademark of Roche Molecular Systems, Inc.
ATCC is a trademark of the American Type Culture Collection.
BBL and Sensi-Disc are trademarks of Becton, Dickinson and Company.
Thermo Fisher Scientific is a trademark of Thermo Fisher Scientific Inc.
Oxoid is a trademark of Thermo Fisher Scientific Inc.
Remel is a trademark of Thermo Fisher Scientific Inc.
Revogene and associated logos are registered trademarks of Meridian Bioscience, Inc.
© 2020-12 Meridian Bioscience, Inc.
Made in Canada

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Revogene e i loghi associati sono marchi commerciali di Meridain Bioscience, Inc.
TaqMan est une marque commerciale déposée de Roche Molecular Systems, Inc.
Thermo Fisher Scientific est une marque commerciale de Thermo Fisher Scientific Inc.
Oxoid est une marque commerciale de Thermo Fisher Scientific Inc.
Remel est une marque commerciale de Thermo Fisher Scientific Inc.
Revogene et les logos connexes sont des marques déposées de Meridian Bioscience, Inc.