## Group A Streptococcus DNA Amplification Assay

DNA Amplification Assay for the Detection of  
*Streptococcus pyogenes* in throat swab samples

REF 480150

IVD

Rx Only

### INTENDED USE

The Alethia Group A *Streptococcus* (Group A Strep) assay, performed on the Alethia reader is a qualitative in vitro diagnostic test for the detection of *Streptococcus pyogenes* (Group A β-hemolytic *Streptococcus*) in throat swab specimens.

The Alethia Group A Strep assay utilizes loop-mediated isothermal DNA amplification (LAMP)<sup>1,2</sup> technology to detect *Streptococcus pyogenes* by targeting a segment of the *Streptococcus pyogenes* genome. Results from the Alethia Group A Strep assay can be used as an aid in the diagnosis of Group A Streptococcal pharyngitis. The assay is not intended to monitor treatment for Group A *Streptococcus* infections.

Alethia Group A Strep is intended for use in hospital, reference or state laboratory settings.

### SUMMARY AND EXPLANATION OF THE TEST

The Alethia Group A *Streptococcus* molecular assay is based on loop-mediated amplification (LAMP) technology. The assay targets a highly conserved 206 base pair (bp) sequence of the *Streptococcus pyogenes* genome residing in a region of the pyrogenic exotoxin B (*speB*) gene.

Loop-mediated amplification uses specially designed primers to provide for specific and continuous isothermal DNA amplification. A by-product of amplification is magnesium pyrophosphate, which forms a white precipitate leading to a turbid reaction solution. Reaction solution absorbance characteristics are monitored by the Meridian Alethia Incubator/Reader. Changes in reaction solution absorbance characteristics created by precipitation of magnesium pyrophosphate indicate the presence of target DNA. The absence of target DNA results in no significant change in sample absorbance.

The Alethia Group A Strep kit includes Alethia Sample Preparation Apparatuses and Alethia Group A *Streptococcus* Test Devices. The Alethia Sample Preparation Apparatus II, used for specimen dilution and preparation, is filled with a buffered solution containing formalin treated *E. coli* harboring *Staphylococcus aureus* DNA. The Alethia Group A Strep Test Device contains one lyophilized amplification reagent bead in each of two chambers: a TEST chamber with Group A *Streptococcus*-specific primers and a CONTROL chamber with *S. aureus*-specific primers. The *S. aureus* DNA in the Sample Preparation Apparatus and the *S. aureus*-specific primers in the CONTROL chamber function as the Internal Control for the assay. During specimen preparation, each patient specimen is added to the Sample Preparation Apparatus II and combined with the *S. aureus* DNA prior to amplification. Addition of *S. aureus* DNA to the patient sample allows for parallel processing of target DNA and Control DNA through amplification and detection. The Internal Control monitors amplification inhibition, assay reagent performance and sample processing effectiveness. The Control *S. aureus* target must be amplified and detected in the final reaction or the test is considered invalid and patient results are not reported.

The Alethia reader monitors changes in absorbance characteristics by measuring transmission of light through the Test and Control reaction solutions. Light transmission is checked at the assay Run Start (Signal<sub>initial</sub>, S<sub>i</sub>) and at the assay Run End (Signal<sub>final</sub>, S<sub>f</sub>). The Alethia reader calculates the change in light transmission between Run End and Run Start (S<sub>f</sub>:S<sub>i</sub>) and compares the ratio to a fixed cut-off value.

Fixed cut-off values for the TEST chamber are used to report sample results. TEST chamber S<sub>f</sub>:S<sub>i</sub> ratios less than 82% are reported as 'POSITIVE'; TEST chamber S<sub>f</sub>:S<sub>i</sub> ratios greater than or equal to 82% are reported as 'NEGATIVE'. Numerical values are not reported.

Fixed cut-off values for the CONTROL chamber are used to determine validity. CONTROL chamber S<sub>f</sub>:S<sub>i</sub> ratios less than 90% are considered valid and allow for reporting of TEST chamber results (POSITIVE, NEGATIVE). CONTROL chamber S<sub>f</sub>:S<sub>i</sub> ratios greater than or equal to 90% are considered invalid and prevent reporting of TEST chamber results. Invalid CONTROL chamber reactions are reported as 'INVALID'. Numerical values are not reported.

More stringent cut-off criteria are applied to the CONTROL chamber reaction to ensure amplification is not inhibited, reagents are performing as intended and that sample processing was performed appropriately.

### BIOLOGICAL PRINCIPLES

Group A *Streptococcus*, or *Streptococcus pyogenes*, is a bacterium commonly found in the human throat or on skin. *S. pyogenes* causes a wide variety of diseases in humans, the most common being acute pharyngitis or strep throat. Pharyngitis is diagnosed in approximately 11 million patients in the United States each year, with Group A *Streptococcus* accounting for 15% - 30% of cases in children and 5% - 20% of cases in adults.<sup>3</sup> Streptococcal pharyngitis is often accompanied by sore throat, fever, tonsillar exudates and enlarged cervical lymph nodes.

Group A *Streptococcus* infections may result in mild illness (e.g. pharyngitis, impetigo) or may lead to invasive, life-threatening illness (e.g. cellulitis, bacteremia, necrotizing fasciitis, Streptococcal toxic shock syndrome). Approximately 25% of patients with necrotizing fasciitis will die from the infection; more than 35% of patients with streptococcal toxic shock syndrome will die from the infection.<sup>4</sup> *Streptococcus pyogenes* may also present in healthy, asymptomatic patients. Throat culture surveys of healthy children taken during school outbreaks have shown a carrier prevalence of up to 20%.<sup>3</sup>

Signs and symptoms of Group A *Streptococcus* pharyngitis often overlap with other infections, making laboratory diagnosis an important tool in identification and treatment of infection. Specimens for laboratory testing should be obtained by vigorous swabbing of both tonsils and the posterior pharynx.<sup>5</sup> Traditional laboratory diagnosis is performed by rapid antigen testing or culture on sheep blood agar followed by differentiation with latex agglutination.<sup>5</sup> Specificities of rapid antigen tests is variable and confirmation of negative results with follow-up throat swab culture is recommended.<sup>5</sup> Throat swab cultures may be 90 to 95% sensitive when correct sampling and plating techniques are used.<sup>3</sup>

Diagnostic tests such chemiluminescent DNA probes may be as sensitive as standard throat cultures on sheep blood agar.<sup>5</sup>

### REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. Alethia Sample Preparation Apparatus II/Negative Control III (SMP PREP): Tris-buffered solution containing formalin-treated *E. coli* harboring *S. aureus* DNA and sodium azide (0.09%) as a preservative.
2. Alethia Group A Strep Test Device: Two-chambered device containing lyophilized amplification reagents (DNA polymerase, deoxynucleotide triphosphates) and either Group A *Streptococcus*-specific primers (TEST Chamber) or *Staphylococcus aureus*-specific primers (CONTROL Chamber).
3. Alethia Heat Treatment Tubes

### MATERIALS PROVIDED SEPARATELY

Alethia Group A *Streptococcus* (Group A Strep) External Control Kit, Catalog Number: 479910

### MATERIALS NOT PROVIDED

1. Disposable latex gloves, powder free
2. DNase/RNase-free, aerosol resistant pipette tips
3. Specimen collection and transport system
4. Swabs, breakable plastic shaft: Rayon, Polyester or Flocked Nylon
5. Non-nutritive Transport Medium: Liquid Amies, without charcoal; Liquid Stuart

### EQUIPMENT NOT PROVIDED

1. Dry-bath with 12 mm heat block capable of 95 C
2. Digital thermometer with Max/Min Temperature Memory (eg, Traceable® Lollipop™ Waterproof/Shockproof Thermometer)
3. Vortex mixer
4. Interval timer
5. Micropipette capable of dispensing 50 µL
6. Alethia reader Meridian Bioscience, Inc. Catalog Number: 610189

### PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Do not interchange Sample Preparation Apparatuses or Test Devices between lots. Heat Treatment Tubes are interchangeable provided they are within assigned expiration date when used.
3. Follow Biosafety Level 2 and Good Laboratory practices during testing.<sup>6</sup> Treat all specimens and used Test Devices as capable of transmitting infectious agents. Do not eat, drink or smoke in areas where specimens or kit reagents are handled.
4. Wear disposable gloves while handling specimens and thoroughly wash hands afterwards.
5. Quality Control Programs for Molecular Testing Laboratories, including proper use and care of equipment, should be employed.<sup>6</sup>
6. The Alethia Group A Strep Test Device contains lyophilized reagents. The protective pouch should not be opened until ready to perform the assay.
7. The Alethia Group A Strep Test Device includes a latch feature that is designed to prevent contamination of the test area with amplification product. Do NOT use Test Devices with broken latches.
8. Dispose of used Alethia Test Devices immediately after processing, leaving the device latch securely in place. Do NOT open the Test Device after processing. Opening the device after amplification may result in contamination of the test area with amplification product.

### HAZARD and PRECAUTIONARY STATEMENTS

There are no known hazards associated with this product.

### SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-27 C.

### REAGENT PREPARATION

Ensure kit reagents are at room temperature (21-27 C) before use. Incorrect results may be obtained if reagents are not brought to room temperature prior to use.

**SPECIMEN COLLECTION AND PREPARATION**

Sample type: Throat swabs

**Sample Collection:** Throat sample collection should be performed in accordance with institutional guidelines for collection of clinical specimens for culture of Group A *Streptococcus*. Samples should be collected by vigorously swabbing the tonsils and the posterior pharynx.

Place swab(s) in a non-nutritive transport medium (e.g. Liquid Amies, without charcoal or Liquid Stuart) and transport to the laboratory. Samples should be held at 2-27 C during transport.

Samples may be held at 21-27 C for up to 48 hours prior to testing. When testing will not be initiated within this time, the sample may be stored at 2-8 C for up to seven days.

**SPECIMEN PREPARATION**

**NOTE:** Ensure that the Alethia reader is powered on and required performance verifications have been completed prior to initiation of SPECIMEN PREPARATION. Refer to the Alethia incubator/reader Operator's Manual for further information regarding instrument set-up and operation.

1. Separate double swabs and remove one swab for Alethia testing.
2. Add the throat swab to the **Alethia SMP PREP**. Break off the handle of the swab. Replace and secure **SMP PREP** cap. Vortex the **SMP PREP** for a minimum of 10 seconds. Specimens in the **Alethia SMP PREP** may be held at 2-29 C for up to 2 hours prior to heat treatment.
3. Remove the tip cap from the **SMP PREP** and squeeze five to 10 drops of sample into a clean Alethia Heat Treatment Tube.
4. Repeat Sample Preparation Steps for all samples to be processed.
5. Heat each Sample/Control mixture in a dry-bath/heat block at  $95 \pm 5$  C for  $10 \pm 2$  minutes. Monitor heat-treatment step with digital thermometer and interval timer.
6. Remove each Heat Treatment Tube from the dry-bath/heat. Heat-treated samples may be held at 21-29 C for up to one hour prior to testing.

**TEST PROCEDURE**

**NOTE:** A maximum of 10 samples can be processed in a single Alethia run.

1. Vortex heat-treated samples for approximately 10 seconds.
2. Remove 1 Alethia Group A Strep Test Device from its protective pouch per sample. Carefully open the device by lifting the latch, holding the chambers such that the lyophilized reagent will not fall out upon opening. Place device on a flat surface or in a rack that can accommodate the device.
3. Using a micropipette, first transfer 50  $\mu$ L of the heat-treated sample to the TEST chamber (left/White Bead), and then transfer 50  $\mu$ L of the heat treated sample to the CONTROL chamber (right/Yellow Bead) of the Alethia Test Device. Take care not to introduce air to the reaction mixture. Do not mix reactions with pipette. Close the Alethia Test Device and fasten the latch securely.
4. Tap device on the bench top or mix to remove air bubbles. Carefully examine the test device for dissolution of the Control/Test Bead, for air bubbles left in the tube and liquid in the top of the device. If undissolved beads, air bubbles or liquid is noted, tap the device on the bench top and repeat the visual inspection.
5. Insert the Alethia Test Device into the Alethia reader and initiate amplification reaction and detection. Results will be displayed at the conclusion of the run.

**INTERPRETATION OF RESULTS**

Sample ID	Reported Result	Interpretation
Patient Specimen	POSITIVE	Sample contains <i>Streptococcus pyogenes</i> target DNA.
	NEGATIVE	No <i>Streptococcus pyogenes</i> DNA detected.
	INVALID	<b>No reportable result. Repeat the test using the original sample.</b> Inhibitory patient specimen, improper sample preparation, reagent failure, instrument failure or internal control failure.
Positive Control	POSITIVE	Valid positive control result. Reagents active at time of use, Alethia reader performing correctly.
	NEGATIVE	<b>Incorrect control result.</b> Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services at 1-800-343-3858 (US) or your local distributor.
	INVALID	<b>No reportable result. Repeat entire assay run using original samples.</b> Improper sample preparation, reagent failure, instrument failure or internal control failure.
Negative Control	POSITIVE	<b>Incorrect control result.</b> Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services at 1-800-343-3858 (US) or your local distributor.
	NEGATIVE	Valid negative control result. Reagents active at time of use, Alethia reader performing correctly.
	INVALID	<b>No reportable result. Repeat entire assay run using original samples.</b> Improper sample preparation, reagent failure, instrument failure or internal control failure.
EMPTY WELL	NONE	No Alethia Test Device in the Alethia reader Well. <b>OR</b> The Alethia Test Device present is compromised due to sample preparation failure, dirty device or improperly seated device. <b>Repeat the test using original sample.</b>

**QUALITY CONTROL**

**This test should be performed per applicable local, state, or federal regulations or accrediting agencies.**

1. Each device contains an internal control chamber that controls for amplification inhibition, assay reagents and sample processing effectiveness.
2. The heat treatment step is monitored with an external thermometer and interval timer. Use the max/min temperature memory of the thermometer to ensure that a temperature of  $95 \pm 5$  C is maintained. Use the interval timer to ensure that heat-treatment duration is  $10 \pm 2$  minutes.
3. Good laboratory practice recommends the use of control materials. Users should follow the appropriate federal, state and local guidelines concerning the running of external quality controls.
4. Alethia Group A Streptococcus External Control Reagents are supplied separately (Catalog 479910). It is recommended that the reactivity of each new lot and each new shipment of Alethia Group A Strep be verified on receipt or before use. External control tests should be performed thereafter in accordance with appropriate federal, state and local guidelines. The Alethia Group A Strep test kit should not be used in patient testing if the external controls do not produce the correct results.
5. A separate device must be used for each external control reagent.

**EXPECTED VALUES**

Overall incidence of *Streptococcus pyogenes* in patients tested during the 2012 clinical study was 14.5% (116/798).

**LIMITATIONS OF THE PROCEDURE**

1. The Alethia Group A Strep assay does not distinguish between viable and non-viable organisms.
2. Respiratory infections can be caused by Streptococci of serogroups other than A as well as other pathogens. This device does not differentiate between carriers and infected individuals.
3. The use of antibiotics or over-the-counter medications may suppress the *Streptococcus* growth in culture despite the presence of organisms detectable by nucleic acid tests.
4. As will all diagnostic tests, results should be interpreted together with other information available to the physician.
5. Zincum Acetum 2X, Zincum Gluconicum 1X, as found in Zicam® Homeopathic Cold Remedies, interfere with the Alethia Group A Strep Assay.

**SPECIFIC PERFORMANCE CHARACTERISTICS**

The Alethia Group A *Streptococcus* DNA Amplification Assay was evaluated in 2012 by independent clinical test sites representing geographically distinct regions throughout the United States. Compiled clinical study information is described below.

**Table 1: Alethia Group A Strep Performance Data Summary**

Group A <i>Streptococcus</i> Composite Culture	Alethia Group A Strep		
	Positive	Negative	Total
Positive	100	2	102
Negative	16	680*	696
Total	116	682	798
			95% CI
Sensitivity	100/102	98.0%	93.1 – 99.5%
Specificity	680/696	97.7%	96.3 – 98.6%
Invalid Rate	2/800	0.3%	0.1 – 0.9%

\* Two specimens produced initial invalid results that were negative upon re-test. A total of 798 qualified specimens were evaluated with the test device and composite bacterial culture methods.

Composite Culture Method testing included the Clinical Site Culture Method and a Reference Culture Method performed by Meridian Bioscience. Site culture testing was performed by direct plating of the throat swab specimen while reference culture testing was performed by plating the swab transport media. Specimens producing positive Group A *Streptococcus* results from either the Site Culture Method or the Reference Culture Method were considered positive.

Table 2 and Table 3 summarize performance information by Clinical Site. Site data was analyzed based on Composite Culture Methods and the Site Culture Method. Seven of the 102 Composite Culture-positive specimens were positive by Site Culture methods only; 26 of the 102 Composite Culture-positive specimens were positive only by the Reference Culture method; the two Alethia false-negative results were positive only by the Reference Culture performed by Meridian Bioscience. Statistical analysis of Site performance data was performed with no significant difference among Sites identified.

**Table 2: Alethia Group A Strep Assay Performance by Site; Composite Culture Method**

	Positive Specimens			Negative Specimens		
	Alethia / Composite Culture	% Sensitivity	95% CI	Alethia / Composite Culture	% Specificity	95% CI
Total	100/102	98.0%	93.1 - 99.5%	680/696	97.7%	96.3 - 98.6%
Site 1	47/47	100.0%	92.4 - 100.0%	287/291	98.6%	96.5 - 99.5%
Site 2	28/30	93.3%	78.7 - 98.2%	204/212	96.2%	92.7 - 98.1%
Site 3	25/25	100.0%	86.7 - 100.0%	189/193	97.9%	94.8 - 99.2%

**Table 3: Alethia Group A Strep Assay Performance by Site; Site Culture Method**

	Positive Specimens			Negative Specimens		
	Alethia / Site Culture	% Sensitivity	95% CI	Alethia / Site Culture	% Specificity	95% CI
Total	74/74	100.0%	95.1 - 100.0%	682/724	94.2%	92.3 - 95.7%
Site 1	40/40	100.0%	91.2 - 100.0%	287/298	96.3%	93.5 - 97.9%
Site 2	18/18	100.0%	82.4 - 100.0%	206/224	92.0%	87.7 - 94.9%
Site 3	16/16	100.0%	80.6 - 100.0%	189/202	93.6%	89.3 - 96.2%

Seventy two (9.0%) patients tested were two years of age or younger; 385 (48.2%) patients were between two and 12 years of age; while 259 (32.5%) patients were greater than 12 and less than 21 years of age. The remaining 82 (10.3%) study patients were greater than 21 years old. Age information was known for all patients included in the performance analysis. No performance differences were noted based on chronological age.

The study population included 410 (51.4%) female patients and 386 (48.4%) male patients. Gender was unknown for two (0.2%) of the study participants. No performance differences were noted based on gender.

**ANALYTICAL SENSITIVITY**

The analytical sensitivity or Limit of Detect for the Alethia Group A Strep assay was determined for common *S. pyogenes* strains.

Limit of Detect was determined using a minimum of 20 replicates for each measurand and a stated probability (eg. 95% where 19/20 replicates are positive) of obtaining positive responses. Analytical sensitivity testing is summarized below:

Streptococcus pyogenes Strain Description	CFU/Test
ATCC 12344	400
ATCC 19615	430

**ASSAY REACTIVITY**

The following *S. pyogenes* strains were tested and produced positive reactions at or below the stated assay limit of detect of 400 CFU/Test: ATCC12384, NCIMB 13285, CCUG 33061, CCUG 33409, CCUG 39158, ATCC 49399, CCUG 53553.

**REPRODUCIBILITY**

Blind-coded panels of 10 samples were supplied to three independent laboratories for reproducibility studies. Samples were randomly sorted within each panel to mask sample identities. The panels included contrived samples manufactured as low positive samples (ie. near the assay limit of detection, n=3) and high negative samples, (n=3). The panels also included contrived positive (n=3) samples and natural negative samples (n=1). Testing was performed by different operators at each site on the same day (intra-assay variability) for five days (inter-assay variability). Three lots of Alethia Group A Strep and five Alethia instruments were used in this study. Positive and Negative Controls were tested each day of testing. The results are given in the table below:

Sample Type	Clinical Site 1		Clinical Site 2		Clinical Site 3		Total	
	Percent agreement		Percent agreement		Percent agreement		Percent agreement	
Negative	10/10	100%	10/10	100%	10/10	100%	30/30	100%
High Negative	29/30	96.7%	30/30	100%	28/30	93.3%	87/90	96.7%
Low Positive	30/30	100%	30/30	100%	30/30	100%	90/90	100%
Positive	30/30	100%	30/30	100%	30/30	100%	90/90	100%
Negative Control	10/10	100%	10/10	100%	10/10	100%	30/30	100%
Positive Control	10/10	100%	10/10	100%	10/10	100%	30/30	100%

**CROSSREACTIVITY**

Crossreactivity studies were performed with positive and negative specimens inoculated with bacterial or fungal organisms to a final concentration of  $1.2 \times 10^8$  CFU/mL. Human DNA was tested at 0.02 mg/mL with no cross-reactivity observed. None of the following organisms reacted with Alethia Group A Strep:

*Acinetobacter baumannii*, *Acinetobacter lwoffii*, *Aeromonas hydrophila*, *Arcanobacterium haemolyticum*, *Bordetella bronchiseptica*, *Bordetella holmesii*, *Bordetella parapertussis*, *Bordetella pertussis*, *Burkholderia cepacia*, *Campylobacter jejuni*, *Candida albicans*, *Citrobacter freundii*, *Corynebacterium diphtheriae*, *Corynebacterium pseudodiphtheriticum*, *Enterococcus faecalis*, *Escherichia coli*, *Haemophilus influenzae*, *Klebsiella oxytoca*, *Klebsiella pneumoniae*, *Lactococcus lactis*, *Legionella jordanis*, *Legionella micdadei*, *Legionella pneumophila*, *Listeria monocytogenes*, *Moraxella catarrhalis*, *Morganella morganii*, *Neisseria gonorrhoeae*, *Neisseria meningitidis*, *Proteus mirabilis*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, *Pseudomonas fluorescens*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Stenotrophomonas maltophilia*, *Streptococcus agalactiae*, *Streptococcus anginosus*, *Streptococcus bovis*, *Streptococcus canis*, *Streptococcus dysgalactiae* (subspecies *equisimilis*), *Streptococcus equinus*, *Streptococcus intermedius*, *Streptococcus mitis*, *Streptococcus mutans*, *Streptococcus pneumoniae*, *Streptococcus salivarius*, *Streptococcus suis*, *Streptococcus uberis*, *Streptococcus sp. Vindans type*.

Two organisms produced unexpected results during initial testing that were not confirmed through further testing. *Arcanobacterium haemolyticum* produced false-negative results for one of 13 replicates tested. *Haemophilus influenzae* produced false-positive results for one of 13 replicates tested.

*Mycoplasma pneumoniae* was tested at  $1.5 \times 10^6$  CFU/mL with no reaction with the Alethia Group A Strep assay.

#### TESTS FOR INTERFERING SUBSTANCES

Interfering substance testing was performed by adding potentially interfering substances diluted in sterile saline or contrived positive samples (ATCC 12344, ATCC 19615) to the **Alethia SMP PREP**. Rayon swabs inoculated with negative throat flora and Liquid Amies Transport Media were added to the **Alethia SMP PREP** containing potentially interfering substances and tested.

The following biological substances, at the saturated solvent/diluent concentrations indicated, do not interfere with test results: Mucus (5.0mg/mL), Human saliva (10% v/v), and Whole Blood (2.5% v/v). Whole Blood at concentrations greater than 2.5% v/v may interfere with the Alethia Group A Strep assay

The following chemical substances, at the saturated solvent/diluents concentrations indicated, do not interfere with test results:  
Acetaminophen (19.5 mg/mL), Aspirin (12.3 mg/mL), Cepacol® Mouthwash, [Cetylpyridinium Chloride (0.005% v/v)], Cepacol® Sore Throat Lozenges [Benzocaine (0.09 mg/mL), Menthol (0.02 mg/mL)], Chloraseptic® Oral Anesthetic/Analgesic [Phenol (0.07% v/v)], Contac® Cold + Flu Tablets [Acetaminophen (16.2 mg/mL), Chlorpheniramine maleate (0.06 mg/mL), Phenylephrine HCl (0.16 mg/mL)], Crest® Complete Toothpaste [Sodium fluoride (0.1 mg/mL)], Diphenhydramine HCl (2.7 mg/mL), HALLS® Cough Drops [Menthol (0.08 mg/mL)], Ibuprofen (15.6 mg/mL), Listerine® Antiseptic Mouthwash [Eucalyptol (0.0092% v/v), Menthol (0.0042% v/v), Methyl salicylate (0.0060% v/v), Thymol (0.0064% v/v)], Robitussin® Cough/Chest Congestion Cough Syrup [Dextromethorphan HBr (0.2 mg/mL), Guaiifenesin (2.0 mg/mL)].

Zicam® Oral Mist [Zincum Aceticum 2X, Zincum Gluconicum 1X (0.625% v/v)] produced invalid results in all replicates tested.

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## REFERENCES

1. Nagamine K, Hase T, Notoni T. Accelerated reaction by loop-mediated isothermal amplifications using loop primers. Mol Cell Probes. 2002; 16:223-29.
2. Mori Y, Kitao M, Tomita N, Notoni T. Real-time turbidimetry of LAMP reaction for quantitating template DNA. J Biochem Biophys. 2004; 59:145-47.
3. Choby BA. Diagnosis and treatment of Streptococcal Pharyngitis. Am Fam Phys. 2009; 79(5):383-90.
4. US Department of Health and Human Services, Centers for Disease Control and Prevention. Group A Streptococcal Disease. [www.cdc.gov](http://www.cdc.gov).
5. Pickering LK, ed. American Academy of Pediatrics. Red Book: 2012 Report of the Committee on Infectious Diseases. 29th ed. Elk Grove Village, IL: American Academy of Pediatrics; 2012. US Department of Health and Human Services PHS/CDC/NIH. Biosafety in microbiology and biomedical laboratories. Washington DC: US Government Printing Office, 2007.
6. CLSI: MM3-A2 Molecular diagnostic methods for infectious disease; approved guideline, 2nd ed. Wayne PA: Clinical Laboratory Standards Institute. 2006.

SN11022

REV. 07/21

Manufactured By

**Meridian Bioscience, Inc.**  
Corporate Office  
3471 River Hills Drive  
Cincinnati, Ohio 45244 USA  
Telephone: 513.271.3700  
Orders/Customer Service:  
800.543.1980  
Technical Support Center:  
800.343.3858  
Information Fax: 513.272.5432  
Ordering Fax: 513.271.0124

Authorized Representative

**Meridian Bioscience Europe S. r. l**  
Via dell'Industria, 7  
20035 Villa Cortese, Milano  
ITALY  
Tel: +39 0331 43 36 36  
Fax: +39 0331 43 36 16  
Email: [info@meridianbioscience.eu](mailto:info@meridianbioscience.eu)  
WEB: [www.meridianbioscience.com/eu](http://www.meridianbioscience.com/eu)

**Australian Sponsor**  
**Emergo Australia**  
**Level 20, Tower II**  
**Darling Park**  
**201 Sussex Street**  
**Sydney, NSW 2000**  
**Australia**

Meridian Bioscience Europe s.a./n.v.  
2 Avenue du Japon - 1420 Braine l'Alleud  
BELGIUM  
Tel: +32 (0) 67 89 59 59  
Fax: +32 (0) 67 89 59 58  
Email: [info.bnl@meridianbioscience.eu](mailto:info.bnl@meridianbioscience.eu)

Meridian Bioscience Europe France  
34 rue de Ponthieu - 75008 Paris  
FRANCE  
Tel: +33 (0) 1 42 56 04 40  
Fax: +33 (0) 9 70 06 62 10  
Email: [info.fr@meridianbioscience.eu](mailto:info.fr@meridianbioscience.eu)

Meridian Bioscience Europe B.V.  
Zekeringstraat 17 A  
1014BM Amsterdam  
NETHERLANDS  
Tel: +31 (0) 411 62 11 66  
Fax: +31 (0) 411 62 48 41  
Email: [Info.bnl@meridianbioscience.eu](mailto:Info.bnl@meridianbioscience.eu)

## INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product.

**Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simblos, Zeichenerklärung)**

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For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.