EIA for the detection of the Toxins Produced by Enterohemorrhagic *E. coli* in Stool Specimens or Culture Systems

REF 608096

IVD In Vitro Diagnostic Medical Device

INTENDED USE

The Premier EHEC test is a rapid in vitro microwell EIA for the detection of Shiga Toxins I and II (Verotoxins) in stool specimens, broth cultures, from individual colonies or colony sweeps of agar plates. The Premier EHEC test is intended for use as an aid in the diagnosis of enterohemorrhagic *E. coli* (EHEC) infections.

SUMMARY AND EXPLANATION OF THE TEST

Enterohemorrhagic *E. coli* (EHEC) have been isolated from patients who have hemorrhagic colitis and hemolytic-uremic syndrome (HUS).^{5, 8, 14} To date, O157:H7 is the most frequently encountered EHEC strain in stools from these patients.^{9, 13} This is probably due to the fact that conventional diagnostic strategies such as O157 latex agglutination assays rely on the unique sorbitol negative fermentation property of this strain.^{2, 9, 13, 14} The major weakness in this approach is that at least 50 other EHEC serotypes have been reported to be associated with the development of hemorrhagic colitis and HUS.^{7-11, 16}

One virulence trait of all EHEC strains is the ability to produce cytotoxin(s) called Shiga toxin (ST) or verotoxin (VT).¹⁴ ST-I and ST-II are the two most common toxins and individual EHEC strains have the ability to produce both or either, in varying quantities. Therefore, ST production and not individual (O157:H7) serotype identification is a better diagnostic strategy for the determination of EHEC associated disease.^{7, 8, 16} Cytotoxin can be identified by a specific cytotoxin assay described by Karmali.⁸ The cytotoxin assay, however, is labor intensive, requires tissue culture facilities, has not been standardized and may take up to 72 hours to confirm the presence of cytotoxin.⁸

Exploiting the EHEC attribute to produce these toxins, Premier EHEC was developed for the direct detection of ST producing strains from stool specimens or culture systems.

BIOLOGICAL PRINCIPLES

The Premier EHEC test utilizes monoclonal anti-Shiga toxin capture antibodies adsorbed to microwells.⁴ Diluted samples are added to the wells and incubated at room temperature. A wash is performed to remove unbound material. Polyclonal anti-Shiga-like toxin antibodies are added for detection and incubated at room temperature. Another wash is used to remove unbound antibody. Enzyme conjugated anti-IgG polyclonal antibody is added and incubated at room temperature. If toxin is present, a reactive antibody-enzyme complex is formed. After washing to remove unbound conjugate, substrate is added and incubated for 10 minutes at room temperature. Color develops in the presence of bound enzyme. Stop Solution is added and the results are interpreted visually or spectrophotometrically.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this kit is listed on the outer box.

1. **Premier EHEC Antibody Coated Microwells** – Breakaway plastic wells coated with monoclonal antibodies specific for *E. coli* Shiga toxin I and II.
2. **Premier EHEC Positive Control** – Inactivated Shiga toxin in a buffered protein solution with preservative.
3. **Premier EHEC Negative Control** – Buffered solution with a preservative.
4. **Premier EHEC Sample Diluent** – Buffered Protein solution with a preservative.
5. **Premier 20X Wash Buffer II** – Concentrated wash buffer with a preservative.
6. **Premier EHEC Detection Antibody** – Rabbit antibodies specific for Shiga toxins in buffered protein solution containing preservative.
7. **Premier EHEC Enzyme Conjugate** – Goat anti-rabbit antibody conjugated to horseradish peroxidase in buffered protein solution containing preservative.
8. **Premier Substrate II** – Buffered solution containing urea peroxidase.
9. **Premier Stop Solution II** – 2N Sulfuric Acid. **CAUTION:** Avoid contact with skin. Flush with water if contact occurs.
10. **Transfer pipettes** - Each pipette is marked to indicate 50 µL, 100 µL, 200 µL and 300 µL volumes.
11. **Microwell strip holder**
12. **Microwell strip sealers**

MATERIALS NOT PROVIDED

1. Wooden applicator sticks
 2. Test tubes (12 x 75 mm) for dilution of sample
 3. Distilled or deionized water
 4. Squirt bottle
 5. Graduated cylinder for making 1X Wash Buffer
 6. EIA plate reader capable of reading absorbance of 450 nm or 450/630 nm*
 7. STAT Fax™ 2200 plate shaker or equivalent (Optional – Broth Method Only)*
- * NOTE: It is the operator's responsibility to validate readers and plate shakers prior to their use with this product.

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Patient specimens may contain infectious agents and should be handled at Biosafety Level 2 as recommended in the CDC/NIH manual, "Biosafety in Microbiology and Biomedical Laboratories".
3. All reagents should be mixed gently before use.
4. Do not interchange Microwells, Detection Antibody, Enzyme Conjugate or Positive Control Reagent between different kits. 20X Wash Buffer II, Sample Diluent, Substrate II and Stop Solution II can be interchanged provided they are within their expiration dates at the time of testing.
5. Allow reagents to warm to room temperature before use.
6. Hold reagent vials vertically at suitable distance above the well to insure proper drop size and delivery.
7. Do not use kit components beyond labeled expiration date.
8. Replace colored caps on correct vials.
9. Dispose of used wash buffer and all test materials in appropriate container. Treat as potential biohazard.
10. The Positive Control reagent contains inactivated Shiga toxin. It should be handled, however, as a potentially hazardous material.
11. Avoid skin contact with Stop Solution II (2N Sulfuric Acid). Flush with water immediately if contact occurs.
12. Do not reuse microwells.
13. Unused microwells must be placed back inside resealable pouch. It is important to protect strips from moisture.
14. Transfer pipettes provided must be used for specimen preparation and transfer. Use one per specimen.
15. Avoid splashing when dispensing diluted samples into microwell by placing transfer pipette tip about halfway down the well and dispensing slowly down the side of well.
16. Microwell washing is to be performed precisely as directed in assay procedure. **Inadequate washing will cause elevated background and false positive results.**
17. All reagents except the 20X Wash Buffer II are provided already diluted to the proper concentration.
18. Any deviation below or above set incubation times may affect sensitivity and specificity and should be avoided.
19. Do not use microwells with pouches that have been damaged (ie, show holes or punctures).

HAZARD and PRECAUTIONARY STATEMENTS


Refer to the SDS, available at www.meridianbioscience.com for Hazards and Precautionary Statements.

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-8 C and return the kit promptly to the refrigerator after each use.

PROCEDURAL NOTES

The Premier EHEC transfer pipette is diagrammed below:

REAGENT PREPARATION

1. Bring the entire kit, including microwell pouch, to room temperature before use.
 2. Prepare 1X Wash Buffer as needed.
- For example: 4.0 mL of 20X Wash Buffer II + 76.0 mL of distilled or deionized water is sufficient to wash one strip. Place in a clean squirt bottle. The 1X Wash Buffer can be stored at room temperature for up to three months.

SPECIMEN COLLECTION AND PREPARATION

1. CDC Recommendations for Stool Specimen Storage, Isolation and Identification of EHEC Organisms:
 - Raw stool specimens should be examined as soon as they are received in the laboratory. If not processed immediately, they should be placed at 2-8 C or frozen at ≤ -70 C.
 - Refrigerated raw stool specimens should be examined within 1-2 hours. If stools cannot be examined within this time, they should be placed in a Cary-Blair-based transport medium. All rectal swabs should be placed immediately in a Cary-Blair-based transport medium.
 - If specimens in transport medium will be examined in 2-3 days, they should be refrigerated. If specimens in transport medium are not examined in 2-3 days, they should be frozen immediately at ≤ -70 C. Specimens in transport medium should not be refrigerated for days, then frozen or left for any time at room temperature.
2. Stool and Broth Storage Prior to Premier EHEC Testing: Stool specimens and broths may be held up to seven days at 2-8 C before testing in the EIA. If testing is not performed within this time period, the stool and/or broth should be frozen at (≤ -70 C). Repeated freeze-thaws should be avoided.
 - A. Direct Stool:**
 1. Add 200 µL of Sample Diluent to a clean 12 x 75 mm test tube. (NOTE: Third marking from tip on transfer pipette represents 200 µL.)
 2. Mix stool as thoroughly as possible prior to pipetting.
 - a. Liquid, semi-solid or stools in transport medium: Using a transfer pipette, draw stool to the 50 µL calibration point (first mark from the tip of the pipette). Dispense the stool into the Sample Diluent. Using the same pipette, gently withdraw and expel the stool suspension several times, then vortex 15 seconds. Leave transfer pipette in tube for later use.
 - b. Non-pipetable stools: Using a wooden applicator stick, transfer a small "BB" sized portion (3-4 mm diameter) of thoroughly mixed stool into Sample Diluent. Emulsify stool using the wooden applicator stick, then vortex 15 seconds. Place a transfer pipette in the tube.
 - B. Plate Method:**
 1. Add 20 µL of stool or specimen to MacConkey or Sorbitol/MacConkey plate and spread with an inoculating loop.
 2. Incubate 16-24 hours at 37 C.
 3. Individual colonies or colony sweeps can be removed with a loop and placed in 200 µL of Sample Diluent for EIA testing.
 - C. Broth Method:**
 1. Add 10-50 µL of stool to 5 mL of MacConkey's Broth, or GN Broth. Vortex 10-15 seconds.
 2. Incubate 16-24 hours at 37 C.
 3. Add 50 µL of growth to 200 µL Sample Diluent for EIA testing.

TEST PROCEDURE

NOTE: With large numbers of specimens, repetitive or multichannel pipettes may be used for dispensing the reagents.

1. After the pouch has reached temperature, break off the required number of microwells (1 well for each specimen plus 1 positive and 1 negative control well per batch). Place the microwells in the microwell strip holder and record the location of all wells. Unused microwells must be resealed in the pouch immediately.
2. Add 100 µL of diluted specimen (second mark from the tip of the pipette) to the appropriate well (place pipette tip halfway into well and let sample slowly run down side of well).
3. Add 2 free-falling drops of Positive and Negative Control to the appropriate wells. Mix wells by firmly shaking/swirling the plate for 30 seconds.
4. Cut plate sealer to size and press firmly onto top of microwells to seal. Incubate the plate for 1 hour at room temperature (22-27 C). Broth Method Only: Alternatively, laboratories equipped with a heated plate shaker (STAT Fax™-2200) can incubate and rotate the plate for 30 minutes at 25 C at 1000 rpm (setting 5).
5. **Carefully**, remove the plate sealer and wash wells:
 - a. Dump plate contents firmly into a biohazard receptacle.
 - b. Bang the inverted plate on a clean stack of paper towels.
 - c. Fill all wells with 1X Wash Buffer. Use of a squirt bottle is recommended.
 - d. Repeat washing cycle (dump, bang on fresh towels, fill) 4 additional times. After the last fill, dump and bang plates on fresh towels hard enough to remove as much excess Wash Buffer as possible but do not allow wells to completely dry at any time.
6. Add 2 free-falling drops of Detection Antibody to each well. Firmly shake/swirl the plate for 30 seconds.
7. Press plate sealer firmly onto top of microwells to seal. Incubate the plate for 30 minutes at room temperature (22-27 C). Broth Method Only: Alternatively, incubate and rotate the plate for 15 minutes on the STAT Fax™ at 25 C and at 1000 rpm (setting 5).
8. Repeat wash procedure (Step #5)
9. Add 2 free-falling drops of Enzyme Conjugate to each well. Firmly shake/swirl the plate for 30 seconds.
10. Press plate sealer firmly onto top of microwells to seal. Incubate the plate for 30 minutes at room temperature (22-27 C). Broth Method Only: Alternatively, incubate and rotate the plate for 15 minutes on the STAT Fax™ at 25 C and at 1000 rpm (setting 5).
11. Repeat wash procedure (Step #5)
12. Clean the underside of all wells with a lint free tissue.
13. Add 2 free-falling drops of Substrate Solution II to each well. Firmly shake/swirl the plate for 30 seconds. Incubate for 10 minutes at room temperature.
14. Add 2 free-falling drops of Stop Solution II to each well. Firmly shake/swirl the plate for 30 seconds. **NOTE:** Initial color of positive reaction is blue, which changes to yellow upon addition of Stop Solution II.
15. Observe Reactions:
 - a. Visual Determination – Read within 15 minutes after adding Stop Solution II.
 - b. Spectrophotometric Determination – Zero EIA reader on air. Wipe underside of wells with a lint free tissue. Read absorbance at 450 nm or 450/630 nm within 30 minutes of adding Stop Solution II.

INTERPRETATION OF RESULTS

1. **Visual Reading**
 - Negative – colorless
 - Positive = definite yellow color
 - A faint yellow color must be evaluated spectrophotometrically
- NOTE:** In view of the epidemiological importance of obtaining bacterial isolates it is recommended that all toxin positive samples be submitted for isolation of Shiga toxin positive organisms. We suggest that individual microbiology laboratories coordinate bacterial isolation with their local state health laboratories.
2. **Spectrophotometric Single Wavelength (450 nm)**
 - Negative = OD₄₅₀ < 0.180
 - Positive = OD₄₅₀ ≥ 0.180
3. **Spectrophotometric Dual Wavelength (450/630 nm)**
 - Negative = OD_{450/630} < 0.150
 - Positive = OD_{450/630} ≥ 0.150

NOTE: A positive result indicates the presence of Shiga toxins. A negative result indicates the absence of Shiga toxins, or that the level of toxin is below that which can be detected by the test.

NOTE: See Quality Control section (Item #5) regarding low positive results. Extremely strong positive reactions may yield a purple precipitate within a few minutes of stopping the reaction.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

- The Positive and Negative Controls must be used with each assay run to provide quality assurance of the reagents.
- The Positive Control should read > 0.500 at either 450 nm or 450/630 nm. The Positive Control should be a definite yellow color when read visually.
- The Negative Control should read < 0.180 at 450 nm and <0.150 at 450/630 nm but greater than 0.00. If the Negative Control is < 0.00, re-blank the plate reader to air and re-read the plate. The Negative Control should be colorless when read visually.
- Any positive well without visible color should be repositioned, wiped on the underside of the well and re-read.
- If the frequency of low positive results (OD between 0.150 and 0.200 with dual wavelength or OD between 0.180 and 0.230 with single wavelength) is greater than 5% of specimens tested, this indicates insufficient washing. More vigorous washing or increasing the washes to seven washes in step 5 of the Procedure is recommended.
- If the expected control reactions are not observed and the reagents are still within their expiration date, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
- At the time of each use, kit components should be visually examined for obvious signs of microbial contamination, freezing or leakage.
- It is suggested that results of each quality control check be recorded in an appropriate log book to maintain high quality testing records.

EXPECTED VALUES

The Premier EHEC test detects the presence of Shiga toxin. Expected values for a given population should be determined for each laboratory. The rate of positivity may vary depending on patient age, geographic location, method of specimen collection, handling and transportation, test employed and general health of the patient population under study.^{4, 5, 9, 16}

Only *Campylobacter* and *Salmonella* ranked above *E. coli* O157:H7 as the most frequently encountered pathogen isolated from stool in the states of Washington and Minnesota during the years 1984-1987. Similar data was obtained during that same time period in Alberta, British Columbia and Ontario, Canada.^{2, 12, 13, 15, 17-19} Direct fecal cytotoxin test have shown that an additional 10% of all specimens may contain undiagnosed EHEC organisms. Our clinical studies show that this number may be conservative.^{3, 7, 10, 11, 16}

Isolation rates of O157:H7 have been shown to be the greatest during the summer months of June, July and August. The rates for non-O157 SLT producing strains still need to be determined.^{2, 12}

LIMITATIONS OF THE PROCEDURE

- The Premier EHEC test detects the presence of Shiga toxin from stool, broths or culture isolates. The level of toxin has not been shown to be correlated with either the presence or severity of disease. As with all *in vitro* diagnostic procedures, test results should be interpreted by a physician in conjunction with other clinical information.
- Shiga toxin I and the toxin produced by *Shigella dysenteriae* type 1 strains (shiga toxin) are nearly identical. Therefore, Premier EHEC will generate a positive signal when Shiga toxin is present in the specimen of other infectious organisms.
- A positive result does not preclude the presence of other infectious organisms.
- Toxin expression may be lost upon serial passage. Colony sweeps may increase the likelihood of detecting ST producing organisms.

SPECIFIC PERFORMANCE CHARACTERISTICS

The Premier EHEC test was evaluated at three major medical centers and one state health laboratory in the United States. Each study site tested diarrheal stools submitted for routine culture of enteric pathogens. Each stool was tested directly by and the broth method in the Premier EHEC assay. A specific ST cytotoxin neutralization assay was performed at Meridian Bioscience on these specimens. The results obtained for the direct stool assay and the broth assay are summarized in the following table.³

	Premier EHEC			
	Direct Stool		Broth Assay	
Cytotoxin	Positive	Negative	Positive	Negative
Positive	15	4	22	0
Negative	17	391	5	236
Indeterminate	2	17	0	3

Relative Sensitivity =	78.9% (56.6-91.5%)*	100% (85.1-100%)
Relative Specificity =	95.8% (93.4-97.4%)	97.9% (95.2-99.1%)
Relative Agreement =	95.1% (92.6-96.8%)	98.1% (95.6-99.2%)

ST producing EHEC were isolated from two of the 17 Premier EHEC positive/cytotoxin negative direct stool specimens: one O157:H7 and one O6:H untypeable. In addition, PCR** confirmed the presence of the ST gene in two Premier EHEC positive/cytotoxin negative and one Premier EHEC positive/cytotoxin indeterminate stool specimens. The 15 Premier EHEC positive/cytotoxin positive direct stool specimens were comprised of: nine – O157:H7, two – O6:H untypeable, one – O untypeable: H12, 41w or 51, and two – O111:NM EHEC strains. An isolate could not be obtained from the remaining stool specimen.

ST producing EHEC were isolated from two of five Premier EHEC positive/cytotoxin negative broths: one O157:H7 and one O111:NM. The 22 Premier EHEC positive/cytotoxin positive samples were comprised of 13 O157:H7, three O6:H untypeable, one O untypeable H:12, 41w or 51, and one O111:NM ST producing strains of EHEC. No isolates could be found in the remaining four broth specimens.

*95% confidence intervals calculated by the normal method as indicated in parenthesis.
**PCR is a research tool and not intended for *in vitro* diagnostic use.

Various *E. coli* ST producing organisms were tested in the EIA by both the Plate and Broth Methods (Specimen Preparation, Part B, and C.). Each strain is a clinical isolate and each was tested by a cytotoxin assay and by a polymerase chain reaction (PCR) to confirm the presence of the ST gene. All organisms generated positive results when tested in this manner. The number in parenthesis represents the number of different strains tested in that group type. Also, the type of toxin produced by each strain is indicated.³

Strain Type	ST Type	Strain Type	ST Type
O157:H7 (3)	I	O91:H21 (1)	II
O157:H7 (7)	II	O146:H21 (1)	I
O157:H7 (9)	I & II	O137:H41 (1)	I & II
O111:NM (3)	I & II	O111:H8 (1)	I & II
OX3:H21 (1)	II	O50:H7 (1)	I
O4:NM (1)	I & II	O145:NM (2)	II
O165:H25 (1)	II	O103:H2 (1)	I
O165:H25 (1)	I & II	O125:NM (2)	I
O45:H2 (1)	I	O26:H11 (1)	I
O126:H27 (1)	I	O5:NM (3)	I
O121:H19 (1)	I & II	O171:NM (1)	II
O121:H19 (1)	II	O83:H1 (1)	I & II

REPRODUCIBILITY

Intra-assay Variability – Ten replicates of two known positive broth specimens and two known stool positive specimens were tested in one assay.

Sample	Mean A _{450/630}	SD	% CV
Broth #1	0.362	0.013	3.5
Broth #2	1.920	0.094	4.9
Stool #1	0.317	0.010	3.1
Stool #2	0.774	0.032	4.2

Inter-assay Variability – Ten replicates of two known positive broth specimens and one negative broth specimen were run on three separate days. This protocol was repeated using one negative and two positive stool specimens.

Sample	Mean A _{450/630}	SD	% CV
Broth #1	0.369	0.032	8.7
Broth #2	1.546	0.293	18.9
Broth #3	0.074	0.007	9.2
Stool #1	0.335	0.027	8.1
Stool #2	0.818	0.079	9.6
Stool #3	0.058	0.007	11.4

CROSSREACTIVITY

The Premier EHEC EIA was tested for specificity using the clinical isolates (CI) or ATCC strains listed. Each strain was tested directly in the EIA by the Plate method using the colony isolate procedure (Specimen Preparation Section, Part B). Every strain was negative when tested in this manner. In an additional study, the bacterial strains were spiked into either a single positive EHEC stool or single negative EHEC stool at a concentration of approximately 2.4 x 10⁹ cfu/mL. The following clinical isolates (CI) or ATCC strains were all found to be negative after direct testing. In addition, the following strains did not alter the expected results when spiked into an EHEC positive or negative stool.³

Description		Description	
<i>Campylobacter coli</i>	CI	<i>Prevotella bivia</i>	ATCC 29303
<i>Campylobacter fetus</i>	CI	<i>Proteus vulgaris</i>	ATCC 6380
<i>Campylobacter jejuni</i>	CI	<i>Providencia alcalifaciens</i>	ATCC 9886
<i>Campylobacter lari</i>	CI	<i>Providencia stuartii</i>	ATCC 33672
<i>Candida albicans</i>	ATCC 10231	<i>Pseudomonas aeruginosa</i>	ATCC 9027
<i>Citrobacter freundii</i>	ATCC 8090	<i>Pseudomonas fluorescens</i>	ATCC 6570
<i>Clostridium perfringens</i>	ATCC 3624	<i>Pseudomonas fluorescens</i>	ATCC 6571
<i>Clostridium perfringens</i>	ATCC 14810	<i>Salmonella Group B</i>	CI
<i>Enterococcus faecalis</i>	CI	<i>Salmonella hilversum (Grp N)</i>	ATCC 15784
<i>Enterobacter cloacae</i>	ATCC 13048	<i>Salmonella minnesota</i>	ATCC 9700
<i>Escherichia coli</i>	ATCC 8739	<i>Serratia liquefaciens</i>	CI
<i>Escherichia coli</i>	ATCC 9637	<i>Shigella dysenteriae</i>	ATCC 11456
<i>Escherichia fergusonii</i>	ATCC 35469	<i>Shigella flexneri</i>	CI
<i>Escherichia faecalis</i>	ATCC 33650	<i>Shigella sonnei</i>	CI
<i>Gardnerella vaginalis</i>	CI	<i>Staphylococcus aureus</i>	ATCC 6538
<i>Helicobacter pylori</i>	CI	<i>Staphylococcus aureus (Cowan I)</i>	ATCC 12598
<i>Klebsiella pneumoniae</i>	ATCC 13883	<i>Staphylococcus epidermidis</i>	ATCC 12228
<i>Lactobacillus acidophilus</i>	ATCC 4356	<i>Streptococcus agalactiae</i>	CI
<i>Mycobacterium smegmatis</i>	CI	<i>Streptococcus faecalis</i>	CI
<i>Neisseria gonorrhoeae</i>	CI	<i>Yersinia enterocolitica</i>	CI
<i>Nocardia asteroides</i>	ATCC 3308	<i>Yersinia enterocolitica O:9</i>	ATCC 55075

ANALYTICAL SENSITIVITY

The limit of detection for ST-I and ST-II are approximately 7 and 15 pg/well respectively.

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REFERENCES

- Brian, MJ., Frosolono, M., Murray, BE., Miranda, A., et al. Polymerase chain reaction for diagnosis of enterohemorrhagic *Escherichia coli* infection and hemolytic-uremic syndrome. *J. Clin. Micro.* 1992;30:1801-1806.
- Centers for Disease Control and Prevention, National Center for Infectious Diseases, *E. Coli* O157:H7, What the clinical microbiologist should know. 1994.
- Data on file, Meridian Bioscience, Inc.
- Donohue-Rolfe, A., Acheson, DWK., Kane, AV. and GT. Keusch. Purification of shiga-toxin and shiga-like toxins I and II by receptor analog affinity chromatography with immobilized P1 glycoprotein and production of cross-reactive monoclonal antibodies. *Infect. Immun.* 1989;57:3888-3893.
- Griffin, PM., Ostroff, SM., Tauxe, RV., Green, KD., Wells, JG., Lewis, JH. and PA. Blake. Illnesses associated with *Escherichia coli* O157:H7 infections. A broad clinical spectrum. *Ann. Intern. Med.* 1988;109:705-712.
- Hull, AE., Acheson, DWK., Echeverria, P., Donohue-Rolfe, A., and GT. Keusch. Mitomycin immunoblot colony assay for detection of shiga-like toxin producing *Escherichia coli* in fecal samples: comparison with DNA probes. *J. Clin. Micro.* 1993;31:1167-1172.
- Karch, H., Böhm, J., Schmidt, J., Gunzer, F., Aleksic, S., and J. Heesemann. Clonal structure and pathogenicity of Shiga-like toxin producing, sorbitol-fermenting *Escherichia coli* O157:H7. *J. Clin. Micro.* 1993;31:1200-1205.
- Karmali, MB. Laboratory diagnosis of verotoxin-producing *Escherichia coli* infections. *Clin. Micro. Newsletter.* 1987;9:65-70
- Kay, BA., Griffin, PM., Stockbine, VA. and JG. Wells. Too fast food: Bloody diarrhea and death from *Escherichia coli* O157:H7. *Clin. Micro. Newsletter.* 1994;16:17-19.
- Mariani-Kurkdjian, P., Deamur, E., Mialn, A., Picard, B., et al. Identification of a clone of *Escherichia coli* O103:H2 as a potential agent of hemolytic-uremic syndrome in France. *J. Clin. Micro.* 1993;31:296-301.
- Maniar, AC., Williams, T., Anand, CM. and GW. Hammond. Detection of verotoxin in stool specimens. *J. Clin. Micro.* 1990;28:134-135.
- Martin, DL., MacDonald, KL., White, KE., Sober, JT. and MT. Osterholm. The epidemiology and clinical aspects of the hemolytic uremic syndrome in Minnesota. *New Engl. J. Med.* 1990;323:1161-1167.
- Neill, MA. *Escherichia coli* O157:H7. A pathogen of no small renown. *Infect. Dis. Newsletter.* 1991;10:19-24.
- O'Brien, AD. and RK. Holmes. Shiga and shiga-like toxins. *Micro Reviews.* 1987;51:206-219.
- Ostroff, SM., Kobayashi, JM. and JH. Lewis. Infections with *Escherichia coli* O157:H7 in Washington State. *JAMA.* 1989;262:355-359.
- Pai, CJ., Ahmed, N., Lior, J., Johnson, WM., Sims, HV., and DE. Woods. Epidemiology of sporadic diarrhea due to verocytotoxin-producing *Escherichia coli*. A two year prospective study. *J. Infect. Dis.* 1988;157:1054-1057.

- Ritchie, M., Partington, S., Jessop, J., and MT. Kelly. Comparison of a direct fecal shiga-like toxin assay and sorbitol-MacConkey Agar culture for laboratory diagnosis of enterohemorrhagic *Escherichia coli* infection. *J. Clin. Micro.* 1992;30:461-464.
- Rowe, PG., Ornbine, E., Lior, J., Wells, GA., McLaine, PN., et al. A prospective study on exposure to verotoxin-producing *Escherichia coli* among Canadian children with haemolytic uremic syndrome. *Epidemiol. Infect.* 1993;110:1-7.
- Wick, PD., Hamacher, ME., Archer, JR. and JS. Wintheiser. Relative incidence of *Escherichia coli* O157:H7 from stool specimens for routine culture in Wisconsin. Abstract. Americana Society for Microbiology, 1991 Annual Meeting.

SN11151

REV. 01/15

Manufactured By

Meridian Bioscience, Inc.
USA/Corporate Office
3471 River Hills Drive
Cincinnati, Ohio 45244
Telephone: 513.271.3700
Orders/Customer Service:
800.543.1980
Technical Support Center:
800.343.3858
Information Fax: 513.272.5432
Ordering Fax: 513.271.0124

Authorized Representative

Meridian Bioscience Europe S. r. l.
Via dell'Industria, 7
20020 Villa Cortese, Milano
ITALY
Tel: +39 0331 43 36 36
Fax: +39 0331 43 36 16
Email: info@meridianbioscience.eu
WEB: www.meridianbioscience.eu

Meridian Bioscience Europe s.a./n.v.
2 Avenue du Japon - 1420 Braine l'Alleud
BELGIUM
Tel: +32 (0) 67 89 59 59
Fax: +32 (0) 67 89 59 58
Email: info.bn@meridianbioscience.eu

Meridian Bioscience Europe France
34 rue de Ponthieu - 75008 Paris
FRANCE
Tel: +33 (0) 1 42 56 04 40
Fax: +33 (0) 9 70 06 62 10
Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.
Postbus 301 - 5460 AH Veghel
NETHERLANDS
Tel: +31 (0) 411 62 11 66
Fax: +31 (0) 411 62 48 41
Email: info.bn@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simbolos, Erläuterung der graphischen symbole)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	EC REP	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG.	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Apareto para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdünnungspuffer befindet
	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Solo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" tests / Contenido suficiente para "<n> ensayos / Inhalt ausreichend für <n> Prüfungen	SOLN STOP	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stoplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung	CONJ ENZ	Enzyme Conjugate / Coniugato enzimatico / Conjugué enzymatique / Conjugado enzimático / Enzymkonjugat
	Serial number / Numero di serie / Número de série / Número de serie / Seriennummer	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät	REAG	Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Date di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	BUF WASH	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Washpuffer
	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenza / Mise En Garde / Advertencia / Warnhinweise
	Conjugate / Coniugato / Conjugué / Conjugado / Konjugat	DIL SPE	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluent échantillon / Diluyente de muestra / Probenverdünnungspuffer
	Substrate / Substrato / Substrat / Substrato / Substrat	BUF WASH 20X	Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
		DET REAG	Detection Reagent / Reagente Diretto / Réactif de Détection / Reactivo de Detección / Nachweis Reagenz

For technical assistance, call Technical Support Services at (800) 343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at (800) 543-1980.