# PREMIER® ADENOCLONE®

## EIA for the Direct Detection of Adenovirus Antigen in Human Fecal Specimens and for the Confirmation of Adenovirus in Cell Culture

REF 696007

IVD In vitro diagnostic medical device

### INTENDED USE

**Premier Adenoclone** is an Enzyme immunoassay (EIA) for the qualitative detection of human adenoviruses directly from stool specimens and confirms the presence of adenovirus in cell culture isolates from respiratory, ophthalmic or enteric specimens.

### SUMMARY AND EXPLANATION OF THE TEST

Human adenoviruses consist of 41 known types, classified primarily by serology and DNA electrophoretotyping. The viruses are nonenveloped, DNA-containing icosahedrons composed of 252 capsomers. The virion contains a number of proteins associated with genus and type specificity. Those associated with genus (group) specificity are the alpha portion of the hexon, the penton base, the delta portion of the fiber and polypeptide IIIa. The monoclonal antibody in **Premier Adenoclone** is against the group-reactive hexon antigen shared by human, simian, canine, porcine, marine and bovine adenoviruses. This monoclonal antibody has been shown to react with all 41 human adenovirus types by EIA<sup>1</sup>.

Adenovirus may cause a variety of clinical illnesses in humans, including respiratory (acute febrile pharyngitis, pharyngoconjunctival fever, acute respiratory disease, pneumonia), ocular (epidemic keratoconjunctivitis, water-borne "swimming pool" conjunctivitis), and gastroenteritis<sup>2,4</sup>.

Identification of adenovirus is generally accomplished by cultivation in tissue culture followed by appearance of characteristic cytopathic effect (CPE). Cell lines susceptible to adenovirus are continuous cell lines of epithelial origin such as HeLa, HEK, HEP-2, KB and Graham 293<sup>2,3</sup>. Generally CPE occurs within 2 to 7 days with adenovirus types 1-7. However, other strains, especially subgroup D, can require up to 28 days before recognizable CPE develops<sup>16</sup>. Confirmation of adenovirus isolated in cell culture after CPE has included immunofluorescent assay, radioimmunoassay, DNA hybridization, enzyme immunoassay and DNA electrophoretotyping<sup>17,21</sup>. **Premier Adenoclone** offers a simple, rapid and accurate means of confirming cell culture isolates for the presence of adenovirus.

The second most commonly found virus in stools or rectal swabs from patients with diarrhea is adenovirus<sup>22</sup>. These fastidious adenovirus types 40 and 41 were first discovered by electron microscopic (EM) examination of stool extracts and are most often associated with gastroenteritis<sup>22</sup>. These adenovirus types can be isolated in Graham 293 cells<sup>5</sup>. **Premier Adenoclone** can be used for direct testing of stool specimens, thus reducing the time required by traditional cell culture isolation.

### BIOLOGICAL PRINCIPLES

**Premier Adenoclone** utilizes monoclonal antibodies in a solid phase sandwich type assay. Breakaway microwells are coated with a monoclonal anti-adenovirus antibody. An aliquot of fecal suspension or undiluted cell culture fluid is added to the microwell and incubated simultaneously with an anti-adenovirus monoclonal antibody conjugated to horseradish peroxidase, resulting in the adenovirus antigen being sandwiched between the solid phase and enzyme antibodies. After 60 minutes incubation at room temperature, the sample well is washed with purified water to remove unbound specimen and the excess enzyme labeled antibodies. Enzyme substrate (urea peroxide) and chromogen (tetramethyl-benzidine) are added to the wells and incubated for 10 minutes at room temperature. Any bound enzyme-conjugate in the wells converts the colorless substrate to a blue color.

### REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. **Microwells** - Coated with Group Specific Anti-Adenovirus Murine Monoclonal Antibody.
2. **Sample Diluent / Negative Control** - Phosphate buffered saline with 0.02% thimerosal as preservative
3. **Enzyme Conjugate** - Horseradish Peroxidase conjugated to anti-adenovirus murine monoclonal antibody in a buffered protein solution with gentamicin and 0.02% thimerosal as preservatives
4. **Positive Control** - Inactivated Adenovirus (AD41) in buffer with 0.02% thimerosal as preservative.
5. **Part A Substrate** - Substrate buffer contains urea peroxide
6. **Part B Substrate** - Substrate buffer contains tetramethylbenzidine
7. **Stop Solution** - Contains 1N sulfuric acid. **Caution:** Avoid contact with skin. If skin contact occurs, flush with water.
8. **Microwell Holder**
9. **Sample Transfer Pipettes**

### MATERIALS NOT PROVIDED

1. Test tubes, 12 x 75 mm
2. Deionized or distilled water
3. Absorbent paper
4. Precision micropipettes and micropipette tips to deliver 100 µL and 1000 µL (optional)
5. Waste container with a 1:10 dilution of household bleach. For autoclaving, use an iodophor disinfectant
6. Microwell EIA plate reader capable of reading absorbance at 450 nm (optional)
7. Device for dispensing wash solution such as multichannel pipette, syringe with manifold wash bottle,
8. Timer (minimum 1 hour)
9. Sterile swabs and viral transport medium (commonly a balanced salt solution containing a viral stabilizing protein and antibiotics). The medium must not be inhibitory to either adenovirus or to cell lines<sup>16</sup>
10. Adenovirus susceptible cell cultures and a light microscope to examine cultures. Susceptible cell cultures for the growth of adenovirus are continuous cell lines of epithelial origin such as HeLa, HEK, HEP-2, KB and Graham 293<sup>2,3</sup>. Only Graham 293 cells should be used for isolation of adenovirus types 40 and 41<sup>5</sup>
11. Adenovirus positive control cultures, i.e. ATCC#VR-2, Adenovirus type 2

### PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Do not mouth pipette samples or reagents. Avoid contact with broken skin or mucous membranes.
3. Do not smoke, eat or drink in areas where specimens or kit reagents are handled.
4. Wear disposable gloves while handling samples and wash hands after assay is complete. Patient specimens, assay controls and all materials coming into contact with them should be handled at Biosafety Level 2 as recommended in the CDC/NIH manual Biosafety in Microbiology and Biomedical Laboratories.
5. **CAUTION:** Because of the affinity of adenovirus for the ocular area, hand to eye contact should be carefully avoided during processing of specimens or performing diagnostic tests<sup>16</sup>.
6. Avoid skin contact with stop solution (1N sulfuric acid). It may cause irritation and burns. Flush with water if contact occurs.
7. Dispose of all materials used to perform the test by autoclaving at 121 C for at least one hour. Liquid waste may be disposed by mixing with a 1:10 dilution of household bleach for a minimum of 30 minutes. **CAUTION:** Liquid waste containing stop solution must be neutralized before addition of household bleach.
8. Avoid splashing or generation of aerosols.
9. Do not use **Premier Adenoclone** reagents beyond the kit expiration date. Each working reagent has been optimized for maximum performance. Dilution or adulteration of these reagents may result in a loss of sensitivity. Incubation time and temperatures other than those specified may give erroneous results.
10. Do not interchange or mix lots of **Premier Adenoclone** reagents.
11. Avoid microbial contamination of reagents or incorrect results may occur. Contamination of samples could cause false results.

12. Use aseptic technique and sterile equipment and materials for all tissue culture procedures.
13. Use separate pipettes or pipette tips for each sample, control and reagent. Avoid contamination with metal ions.
14. **DO NOT REUSE MICROWELLS.**

### HAZARD AND PRECAUTIONARY STATEMENTS

Refer to the SDS, available at [www.meridianbioscience.com](http://www.meridianbioscience.com) for Hazards and Precautionary Statements.

### SHELF LIFE AND STORAGE

Store kit reagents at 2-8 C. Bring kit reagents to room temperature (20-30 C) before use and promptly return to 2-8 C after use. Keep microwells in pouch until the pouch reaches room temperature. Return all unused microwells to their original foil pouch. The microwells are stable for the shelf life of the kit once opened.

### INDICATION OF INSTABILITY OR DETERIORATION

The following conditions may indicate reagent deterioration:

1. Any evidence of microbial contamination or heavy precipitation
  2. Any blue color in the substrate solutions before addition to microwells.
  3. A negative control value greater than 0.150 absorbance units at 450 nm may indicate deterioration of reagents.
  4. A positive control value of less than 0.3 absorbance units may indicate deterioration of reagents.
- If any of the above conditions are observed, contact our Technical Support Center at 1-800-343-3858.

### REAGENT PREPARATION

1. Bring all reagents to room temperature (20-30 C) before use.
2. Return all reagents to 2-8 C immediately after use.
3. Do not allow microwells to dry between steps.
4. Prepare decontamination vessel for discarding reagents and materials.
5. Reproducibility in any EIA assay is largely dependent upon the consistency with which the microwells are washed. Carefully follow the recommended washing sequence as outlined in the EIA test procedure.

### SPECIMEN COLLECTION

Proper specimen collection is critical for the isolation of adenovirus in cell culture and should only be attempted by experienced personnel. The optimum time for specimen collection is as soon as possible after onset of symptoms. Adenovirus is most commonly isolated from respiratory, ophthalmic and rectal specimens, but has been recovered from virtually every organ system in man<sup>16</sup>. The site of infection should be thoroughly swabbed to collect epithelial cells. Place the swab directly in approximately 3 mL of viral transport medium and swirl vigorously<sup>16</sup>. The specimen should be sent promptly to the laboratory on ice according to 42 CFR (Code of Federal Regulation) 72.

Alternatively, nasopharyngeal aspirates, secretions or washes can be collected<sup>16</sup>. It is generally more effective and convenient to collect 2-5 mL naso-pharyngeal secretions by aspiration, using a small tube connected to a specimen collection bottle with the tube inserted into the patient's nose. This system is commercially available as the mucous trap. In cases in which nasal secretions are insufficient, nasal washes with 1-2 mL of sterile saline through the nostril may be used. A minimum of 0.2 mL of secretion should be collected.

Rectal swabs or stool specimens collected for culture should be placed in 1 mL of viral transport media. Caution should be taken when using rectal swabs to assure that sufficient fecal sample (30 to 50 mg of raw stool) is obtained.

All specimens should be held at 2-8 C on ice and not frozen unless a delay of more than 72 hours before cell culture inoculation is anticipated. In such cases, rapidly freeze the specimen in a bath of dry ice and acetone or alcohol and store at -20 C or colder for up to one month. Avoid slow freezing (e.g. placing into a -20 C freezer) as this may destroy infectivity of adenovirus in clinical specimens. Do not store specimens in a self-defrosting freezer. Avoid repeated freezing and thawing of specimens. For direct testing of stool specimens: Stool specimens or rectal swabs should be collected in containers that do not contain media, preservatives, animal serum, metal ions, oxidizing agents or detergents as any of these additives may interfere with the **Premier Adenoclone** assay. Caution should be taken when using rectal swabs to assure that sufficient fecal sample (30 to 50 mg) is obtained. Diluted samples may be stored at 2-8 C for 3 days without interference with the assay performance. For long term storage of undiluted stool specimens, store at -20 C or colder. Do not store in self-defrosting freezers. Avoid repeated freezing and thawing of specimens.

### PREPARATION OF SPECIMENS FOR ISOLATION

Susceptible cell cultures for the growth of adenovirus are continuous cell lines of epithelial origin such as HeLa, HEK, HEP-2, KB and Graham 293<sup>2,3</sup>. Only Graham 293 cells should be used for the isolation of adenovirus types 40 and 41<sup>5</sup>. Different cell lines as well as different passages of the same cell lines may vary greatly in viral sensitivity. Each lab should establish their own QC procedures for ensuring the acceptability of tissue culture cells for routine tissue culture isolation.

Swabs arriving in a viral transport medium should be swirled vigorously and excess fluid removed by firmly pressing the swab against the side of the tube. Inoculate into susceptible cell cultures<sup>16</sup>. Do not accept dry swabs for viral isolation.

All aspirates, secretions and washes should be pipetted vigorously to obtain a homogeneous suspension. The suspension should then be treated with antibiotics and inoculated into susceptible cell culture<sup>16, 26, 27</sup>.

Stool specimens can either be used direct (see assay procedure) or prepared for cell culture isolation. For the latter, grind a few grams of fecal material with a mortar, pestle, sterile sand or other abrasive. Add that volume of diluent (e.g. culture medium) with antibiotics that results in a 10% suspension. Clarify by centrifugation and inoculate into Graham 293 cells<sup>16</sup>.

Ideally, specimen cell cultures should be monitored daily for the development of cytopathic effects (CPE) characteristic of adenoviruses. The cells of an adenovirus infected monolayer become enlarged, rounded and retracted, eventually aggregating into irregular grape-like clusters. These changes are generally associated with an increase in the acidity of the culture medium. The time required for the development of CPE depends on the strain of adenovirus and concentration of infectious virus in the specimen. Generally, CPE occurs within 2 to 7 days with adenovirus types 1, 2, 3, 4, 5, 6 and 7. Other adenovirus strains, especially those in subgroup D, can require up to 28 days before recognizable CPE develops or may even require a blind cell culture passage<sup>16</sup>. It is recommended that cultures are ready for confirmation when at least 25% of CPE becomes evident.

### SPECIMEN PREPARATION

#### Assay of Tissue Culture Specimen Fluid

When tissue cultures exhibit at least 25% to 50% CPE, add 2 drops (100 µL) of tissue culture fluid directly into the microwells, using the sample transfer pipette. **DO NOT DILUTE TISSUE CULTURE FLUID.**


#### Assay of Direct Stool Specimen

Add 1 mL of sample diluent to properly marked tube, using a sample transfer pipette. Add samples by one of the following:

1. Solid stool — press sample into sample transfer pipette to first mark (see figure below).
2. Liquid stool — aspirate sample into sample transfer pipette to first mark (see figure below).
3. Rectal swab — swirl swab in one ml of sample diluent to release fecal material. Firmly press swab against side of tube to remove liquid. Mix thoroughly. *Do not dilute further.*

Resuspend sample in 1 mL of sample diluent with sample transfer pipette by vigorously aspirating and eluting solution until a workable solution is obtained.

#### Sample Transfer Pipette

Using the sample transfer pipette, add 2 drops of either *undiluted* tissue culture fluid or *diluted* stool specimen to wells. See test procedures for details.

## TEST PROCEDURE

- Snap off a sufficient number of wells for samples and the controls and insert into the microwell holder. Record sample position.
- Add 2 drops (100 µL) of either diluted fecal sample, undiluted tissue culture fluid, negative control (sample diluent) or positive control to the bottom of separate microwells.
- Add 2 drops (100 µL) of enzyme conjugate to each microwell. Mix by gently swirling on table top.
- Incubate at room temperature (20-30 C) for 60 ± 5 minutes.
- Pour the liquid out of the wells into a discard vessel. Tap the microwell holder upside down vigorously against absorbent paper to ensure complete removal of liquid from the wells.
- Fill all the wells to overflowing with deionized water and pour the liquid out as in step 5.
- Repeat the washing procedure (steps 5 and 6) four more times for a total of 5 washes.
- Add 2 drops (100 µL) of Substrate A (urea peroxide) solution to each microwell.
- Add 2 drops (100 µL) of Substrate B (TMB) solution to each microwell.
- Incubate at room temperature for 10 minutes.
- After reading and recording results, discard the used wells and absorbent paper into biohazardous waste.

Visual determinations can be made after incubating for 10 minutes in step 10, and should be read within 10 minutes. Samples with blue color greater than the negative control are positive. Samples showing equal or less color than negative control are negative.

### Optional — Spectrophotometric Procedure

Spectrophotometric determinations can be made by adding 2 drops (100 µL) of stop solution (1 N sulfuric acid) to each microwell after the incubation in step 10. Read the absorbance of each microwell at 450 nm using a >600 nm reference filter (optional) against an air blank within 60 minutes.

## INTERPRETATION OF RESULTS

- Visual Determination. Any sample with blue color more intense than that of the negative control is considered positive. Any sample with color equal to or less intense than the negative control is considered negative.
- Spectrophotometric Determination. Specimens with absorbance units ( $A_{450}$ ) greater than 0.150 are considered positive. Specimens with absorbance equal to or less than 0.150 are considered negative. Occasionally, a discrepancy may occur between visual and spectrophotometric determinations with samples containing low amounts of antigen. Spectrophotometric determination, being an objective method, is slightly more accurate.  
NOTE: A precipitate may form in high positive samples. This will not affect the results.

## QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

**Premier Adenoclone** is supplied with a positive and negative (sample diluent) control which are to be used every time an assay is run. These reagents are used to insure that the kit reagents are working properly. The positive control well should be clearly blue with an absorbance at 450 nm of greater than 0.30. The negative control must have no visible color, and an absorbance at 450 nm of less than 0.150.

In addition to the positive and negative controls supplied, when performing tissue culture confirmation, two other types of controls may be used at the discretion of the laboratory. These controls insure that the tissue culture system employed is functioning properly:

- Uninfected cell supernatant from the same lot as those used to isolate the virus can be utilized as a negative control.
- Infected cell supernatant from a known adenovirus positive specimen (25% CPE or greater) can be utilized as a positive control.

The uninfected cell supernatant well should be colorless with an absorbance at 450 nm of less than 0.150. The infected cell supernatant from a known positive should be blue. The absorbance at 450 nm will depend on the degree of infection, but should be greater than 0.150.

If the positive and negative controls supplied with the kit do not perform properly, contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor. If the positive and negative controls supplied work properly but the tissue culture controls do not perform properly, an evaluation of the cell lines and procedures involved is required.

## EXPECTED VALUES

Rate of positivity may vary depending on geographic location, method of specimen collection, handling and transportation, cell culture system used or test employed and general health environment of patient population under study.

The frequency of adenovirus infection will vary with the clinical syndrome and the age of the individual. In children under 5 years, approximately 5% of acute respiratory disease is due to adenovirus<sup>5</sup>.

Approximately 10% of childhood pneumonia may be of adenovirus etiology<sup>7</sup>. Acute hemorrhagic cystitis in children may be caused by adenoviruses in 20 to 70% of cases<sup>8, 9</sup>. Enteric adenoviruses (types 40 and 41) have been implicated in approximately 10% of pediatric patients with gastroenteritis and appears most frequently in children under 2 years of age<sup>10</sup>. Performance characteristics have not yet been fully determined for direct testing of neonatal stool samples. Lack of false positive results should be determined in each individual laboratory before neonatal stools are tested directly.

In adults, adenoviruses have been implicated at times in cervicitis<sup>4</sup>, and in acute respiratory disease, especially in military recruits<sup>11</sup>. Ocular infections such as epidemic keratoconjunctivitis and so called "swimming pool conjunctivitis" due to adenoviruses may occur in any age group<sup>12, 13</sup> and those in any group who are immunosuppressed may become infected with adenoviruses<sup>14, 15</sup>.

## LIMITATIONS OF THE PROCEDURE

**Premier Adenoclone** assay is highly specific and sensitive for adenovirus antigen. The monoclonal antibody in this test reacts with the group specific hexon antigen. It will detect all 41 known serotypes, but cannot be used to differentiate types.

A negative result does not exclude the possibility of adenovirus infection in the patient. Failure to detect adenovirus may be a result of factors such as collection of specimen at an improper time in the disease when too few virions are present, improper sampling or handling of the specimen, failure of cell culture, etc. A negative result for a stool sample may not preclude the presence of an infection with non-enteric adenovirus in other body fluids. If respiratory or ophthalmic infection is suspected, specimens from the site of infection should be cultured.

All positive results must be interpreted with caution since adenovirus is capable of latency and recrudescence. Asymptomatic shedding may occur up to 18 months after infection<sup>16</sup>. Enteric adenoviruses may be found in the stools of asymptomatic children<sup>25</sup>. Test results should be interpreted in conjunction with information available from epidemiological studies or clinical evaluations of patient or other diagnostic procedures<sup>24</sup>.

Co-infection with bacterial pathogens is possible. Therefore, bacteriological tests should be performed in parallel with this test to rule out bacteriological etiology. The use of **Premier Adenoclone** for direct specimen testing other than stool specimens is not recommended as insufficient antigen, inadequate or improper specimen collection may cause negative results. This test is intended for use for culture confirmation or direct stool specimens. A positive result in feces, in association with diarrhea is highly suggestive of adenoviral gastroenteritis<sup>10</sup>. Adenovirus types 40 and 41 have been most often associated with viral gastroenteritis.

False positives could occur with high levels of *S. aureus* possessing protein A such as Cowan strain. *Staphylococcal enterocolitis* is an uncommon disease generally described in adult patients and is extremely rare in infants and children<sup>25, 26</sup>.

## SPECIFIC PERFORMANCE CHARACTERISTICS

**Premier Adenoclone** was tested in four laboratories in New York, California and Massachusetts for tissue culture confirmation or direct testing of stool samples. Results of **Premier Adenoclone** were compared to cytopathic effect (CPE) formation in tissue culture and Electron Microscopy (EM) for stool. **Premier Adenoclone** assay showed excellent correlation with both methods.

Method of Determination*	Site I		Site II		Site III		Site IV	
	S		S, V		V		S, V	
	Tissue Culture		Stool Suspension		Tissue Culture		Tissue Culture	
	+	-	+	-	+	-	+	-
<b>Premier Adenoclone®</b>	75	0	70	0	22	0	23	0
	0	115	1	160	0	51	1	49
Sensitivity	100%		99%		100%		96%	
Specificity	100%		100%		100%		100%	
Total Agreement	100%		99%		100%		99%	

\*S: Spectrophotometric determination; V: Visual determination

NOTE: There was 100% correlation between spectrophotometric and visual determination.

## LIMITS OF DETECTION

### STOOL SUSPENSION

Four adenovirus positive stool samples with known EM particle counts per ml were serially diluted and tested in **Premier Adenoclone** to show the limits of detection. These results show that adenovirus particle counts as low as  $10^7$  per mL can be detected by **Premier Adenoclone**.

Particles/mL	Non-enteric		Non-enteric		Ad 41		Ad 40	
	Stool 1	Stool 2	Stool 1	Stool 2	Stool 3	Stool 4	Stool 3	Stool 4
$10^8$	$A_{450}$	$A_{450}$	$A_{450}$	$A_{450}$	$A_{450}$	$A_{450}$	$A_{450}$	$A_{450}$
	1.78	-	-	-	-	-	-	-
$10^7$	.26*	>2.00	>2.00	>2.00	>2.00	>2.00	>2.00	>2.00
$10^6$	.04	1.06	1.06	1.06	1.06	1.06	1.06	1.06
$10^5$	.03	.17*	.17*	.17*	.17*	.17*	.17*	.17*
$10^4$	-	.04	.04	.04	.04	.04	.04	.04

## CELL CULTURE

Serial dilutions of different adenovirus supernatants, from 3+ CPE cultures, as tested by **Premier Adenoclone** Adenovirus diagnostic kit are shown in the following table. Asterisks indicate the relative sensitivity of **Premier Adenoclone** for the adenovirus serotype tested. The monoclonal antibody used in **Premier Adenoclone** has been shown to react with all human subgroups as well as with isolates that infect monkeys, dogs, cattle and swine<sup>2</sup>.

	TYPE A			TYPE B			TYPE C			
	AD12	AD18	AD31	AD3	AD7	AD21	AD1	AD2	AD5	AD6
1:10	>3.0	>3.0	>3.0	1.78	>3.0	>3.0	>3.0	>3.0	>3.0	>3.0
1:100	2.13	0.64*	2.09	0.49*	>3.0	2.36	>3.0	>3.0	>3.0	>3.0
1:1000	0.21**	0.05	0.15*	0.09	1.56	0.24	>3.0	1.67	0.57	2.52
1:2000	0.08	0.02	0.06	0.05	0.71*	0.08	1.66	0.70	0.23*	1.26
1:4000	0.04	0.01	0.03	0.04	0.15	0.04	0.65	0.19*	0.12	0.60
1:8000	0.02	0.02	0.02	0.03	0.07	0.02	0.22*	0.11	0.06	0.35
1:16000	0.02	0.02	0.02	0.02	0.07	0.02	0.11	0.07	0.03	0.19*

	TYPE D				TYPE E		ENTERIC	
	AD9	AD22	AD37	AD39	AD4	AD40	AD41	
1:10	>3.0	>3.0	>3.0	>3.0	>3.0	>3.0	>3.0	
1:100	2.75	2.44	2.62	2.28	>3.0	2.40	>3.0	
1:1000	0.24*	0.22*	0.25*	0.15*	1.60	0.26*	0.73	
1:2000	0.11	0.09	0.09	0.07	0.94	0.15	0.27*	
1:4000	0.05	0.04	0.06	0.03	0.46	0.08 -	0.12	
1:8000	0.05	0.03	0.02	0.02	0.25	0.06	0.06	
1:16000	0.02	0.02	0.02	0.02	0.15*	0.05	0.04	*limits of detection

## REPRODUCIBILITY

### Intra-assay

Each specimen\* was tested in a single assay twenty-one times and the mean and coefficient of variation determined.

	Neg Stool		Stool 1		Stool 2		Stool 3	
	N	21	21	21	21	21	21	21
$\bar{X}$	0.012	0.232	0.744	1.713				
CV	56%	10%	14%	8%				
	Neg Cult.		Cult.1		Cult.2		Cult.3	
	N	21	21	21	21	21	21	
$\bar{X}$	0.017	0.181	0.418	1.318				
CV	35%	6%	7%	7%				

\*A positive sample was diluted in a 5% stool suspension and in tissue culture fluid to represent a low (near the cut-off), medium and high sample.

### Inter-assay

Each specimen was run individually in twenty-one different assays performed by three technicians and the mean and coefficient of variation determined.

	Neg Stool		Stool 1		Stool 2		Stool 3	
	N	21	21	21	21	21	21	21
$\bar{X}$	0.026	0.129	0.474	0.990				
CV	48%	9%	12%	10%				
	Neg Cult.		Cult.1		Cult.2		Cult.3	
	N	21	21	21	21	21	21	
$\bar{X}$	0.017	0.198	0.508	1.547				
CV	59%	10%	10%	9%				

\*A positive sample was diluted in a 5% stool suspension and in tissue culture fluid to represent a low (near the cut-off), medium and high sample.

## CROSSREACTIVITY

The following microorganisms and viruses were tested in **Premier Adenoclone** and showed no cross reactivity in the test. All cross-reactivity tests utilized stock cultures containing an estimated  $10^7$ - $10^9$  infectious organisms per mL with the following exceptions: calicivirus, SRV, rotavirus, norwalk virus and astrovirus were tested as EM positive stools; *Clostridium difficile* was tested as a toxin positive. Stool and urine was tested and found to be negative.

Microorganisms	Staphylococcus aureus* (Cowan)	Mycoplasma arginini
<i>Pseudomonas aeruginosa</i>		
<i>Enterobacter aerogenes</i>	<i>Clostridium difficile</i>	<i>Mycoplasma hominis</i>
<i>Enterobacter cloacae</i>	<i>Candida albicans</i>	<i>Mycoplasma hyorhinis</i>
<i>Klebsiella pneumoniae</i>	<i>Proteus mirabilis</i>	<i>Mycoplasma orale</i>
<i>Serratia marcescens</i>	<i>Acholeplasma laidlawii</i>	<i>Mycoplasma pneumoniae</i>
<i>Escherichia coli</i>		

## Viruses

Coxsackievirus A-9	Rhinovirus 3	Echovirus 6
Coxsackievirus B-3	Rhinovirus 17	Echovirus 7
Coxsackievirus B-4	Rhinovirus 43	Echovirus 9
Coxsackievirus B-5	Varicella zoster	Echovirus 11
Poliovirus I	Chlamydia trachomatis	Echovirus 17
Poliovirus II	Simian foamyvirus 2	Echovirus 18
Poliovirus III	Simian foamyvirus 8	Echovirus 22
Human Coronavirus	Measles	Herpes simplex 1
Reovirus 1	Respiratory Syncytial virus	Herpes simplex 2
Reovirus 2	Cytomegalovirus	Small round virus (SRV)
Reovirus 3	Influenza A	Rotavirus
Parainfluenza-1	Influenza B	Nonwalk virus
Parainfluenza-2	Calicivirus	Astrovirus
Parainfluenza-3	Papovavirus SV-40	Myxoviruses-Paramyxoviruses SV-5

\*False positive results could occur with high levels of *S. aureus* possessing protein A such as Cowan strain. *Staphylococcal enterocolitis* is an uncommon disease generally described in adult patients and is extremely rare in infants and children<sup>25, 26</sup>.

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## REFERENCES

- Cepko CL, Whetstone CA, and Sharp PA. Adenovirus hexon monoclonal antibody that is group specific and potentially useful as a diagnostic reagent. J Clin Microbiol 1983;17:360-364.
- Herrmann JE, Perron-Henry DM, Stobbs-Walro D, and Blacklow NR. Preparation and characterization of monoclonal antibodies to enteric adenovirus types 40 and 41. Arch. Virology, 94:259-265, 1987.
- Horwitz MS. Adenoviral Diseases. In B.N. Fields, ed. Virology, Raven Press, New York, pp. 477-495.
- Laverty CR, Russell P, Black J, Kappagoda N, Benn RAV, and Booth N. Adenovirus infection of the cervix. Acta Cytol. 21:114-117, 1977.
- Takiff HR, Straus SE, and Garon CF. Propagation and *in vitro* studies of previously non-cultivable enteric adenoviruses in 293 cells, Lancet 11:832-834, 1981.
- Brandt CD, Kim HW, Vargosdo AJ, Jeffries BC, Arrobio JO, Rindge B, Parrott RH, and Chanock RM. Infections in 18,000 infants and children in a controlled study of respiratory tract disease. I. Adenovirus pathogenicity in relation to serologic type and illness syndrome. Am J Epidemiol 1969;90: 484-500.
- Mallet R, Riberre M, Bonenfant F, Labruno B, and Reyrole L. Les pneumopathies graves a adenovirus. Arch. Fr. Pediatr. 23:1057-1073,1966.
- Mufson MA, Belshe RB, Horrigan TJ, and Zollar LM. Cause of acute hemorrhagic cystitis in children. Am J Dis Child 1973;126:605-609.
- Numazaki Y, Kumasaka T, Yano N, Yamanaka M, Miyazawa T, Takai S, and Ishida N. Further study of acute hemorrhagic cystitis due to adenovirus type 11. N Engl J Med 1973;289:344-347.
- Uhnou I, Wadell G, Svensson L, Johansson ME. Importance of enteric adenoviruses 40 and 41 in acute gastroenteritis in infants and young children J Clin Microbiol 1984;20:365-372.
- McGagab WJ. *Mycoplasma pneumoniae* and adenovirus respiratory illnesses in military and university personnel, 1959-1966. Am. Rev. Respir. Dis. 97:345-358, 1968.
- Kemp MC, Hierholzer JC, Cabradilla CP, and Objesti JF. The changing etiology of epidemic keratoconjunctivitis: Antigenic and restriction enzyme analysis of adenovirus types 19 and 37 isolated over a 10 year period. J Infect Dis 1983;148:29-33.
- Foy HM, Cooney MK, and Hatlen JB. Adenovirus type 3 epidemic associated with intermittent chlorination of a swimming pool. Arch. Environ. Health. 17:795-802, 1968.
- De Jong PJ, Valderrama G, Spigland I, and Horwitz MS. Adenovirus isolates from the urines of patients with the acquired immunodeficiency syndrome. Lancet 1:1293-1296, 1983.
- Siegel RP, Kikman SH, Arayata RB, and Bottone EJ. Fatal disseminated adenovirus 11 pneumonia in an agammaglobulinemic patient. Am J Med 1981;71:1062-1067.
- Lennette EH, Schmidt NJ. Diagnostic Procedures for Viral Rickettsial and Chlamydial infections, Fifth edition, 1979, American Public Health Associates, Inc., Washington, D.C., pp 241-242, pp 229-255.
- Darougar S, Walpita P, Thaker U, Viswalingam N, and Wishart MS. Rapid culture test for adenovirus isolation. Br J Ophthalmol 1984;68:405-408.
- Gomel SA, Nascimento JP, Siquiera MM, Krawczuk MM, Pereira HG, and Russel WC. In situ hybridization with biotinylated DNA probes: a rapid diagnostic test for adenovirus upper respiratory infections. J Virol Methods 1985;12:105-110.
- Kidd AH, Harley EH, and Erasmus MJ. Specific detection and typing of adenovirus types 40 and 41 in stool specimens by dot-blot hybridization. J Clin Microbiol 1985;22:934-939.
- Lehtomaki K, Julkunen I, Sandelin K, Salonen J, Virtanen NL, Ranki M, and Hoyt T. Rapid diagnosis of respiratory infections in young adult men. J Clin Microbiol 1986;24:108-111.
- Shuster V, Matz B, Wiegand H, Traub B, and Neumann-Haefelin D. Detection of herpes simplex virus and adenovirus DNA by dot-blot hybridization using *in vitro* synthesized RNA transcripts. J Viral Methods 1986;13:291-299.
- Lennette EH, Schmidt NJ. Diagnostic procedures for viral rickettsial and chlamydial infections, Fifth edition, 1979, American Public Health Associates, Inc., Washington, D.C., pp 942-943.
- Cukor G, and Blacklow NR. Human viral gastroenteritis. Microbiological Reviews. 48:157-179,1984.
- Greenburg SS, and Krilov L. Laboratory diagnosis of viral respiratory disease. Cumitech No. 21, ASM Drew, W.L., Rubin, S.J. editors, March 1986.
- Gutman LT, Idriss ZH, Gehlbach S, Blackman L. Neonatal *staphylococcal enterocolitis*; Association with indwelling feeding catheters and *S.aureus* colonization. J Pediatr 1976;88:836-839.
- Hay P, and McKenzie P. Side effects of oxytetracycline therapy. Lancet 1:945-950. 1954.
- Lennette EH, Schmidt NJ. Diagnostic procedures for viral rickettsial and chlamydial infections, Fifth edition, 1979, American Public Health Associates, Inc., Washington, D.C., p. 81.

SN11246

REV. 01/15

	<b>Meridian Bioscience, Inc.</b> USA/Corporate Office 3471 River Hills Drive Cincinnati, Ohio 45244 Telephone: 513.271.3700 Orders/Customer Service: 800.543.1980 Technical Support Center: 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124
	<b>Meridian Bioscience Europe s.a./n.v.</b> 2 Avenue du Japon - 1420 Braine l'Alleud BELGIUM Tel: +32 (0) 67 89 59 59 Fax: +32 (0) 67 89 59 58 Email: info.bn@meridianbioscience.eu
	<b>Meridian Bioscience Europe France</b> 34 rue de Ponthieu - 75008 Paris FRANCE Tel: +33 (0) 1 42 56 04 40 Fax: +33 (0) 9 70 06 62 10 Email: info.fr@meridianbioscience.eu
	<b>Meridian Bioscience Europe b.v.</b> Postbus 301 - 5460 AH Veghel NETHERLANDS Tel: +31 (0) 411 62 11 66 Fax: +31 (0) 411 62 48 41 Email: info.bn@meridianbioscience.eu

**INTERNATIONAL SYMBOL USAGE**

You may see one or more of these symbols on the labeling/packaging of this product:

**Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de símbolos, Erläuterung der graphischen symbole)**

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	<b>CONTROL +</b>	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
<b>LOT</b>	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	<b>CONTROL -</b>	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
<b>IVD</b>	In vitro diagnostic medical device / Dispositivo médico-diagnóstico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	<b>EC REP</b>	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
<b>CE</b>	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG.	<b>SMP   PREP   DIL   SPE</b>	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdünnungspuffer befindet
<b>REF</b>	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	<b>BUF   RXN</b>	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contains sufficient for <n> tests / Contenuto sufficiente per <n> saggi / Contenu suffisant pour <n> tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	<b>SOLN   STOP</b>	Stopping Solution / Solutione di Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung	<b>CONJ   ENZ</b>	Enzyme Conjugate / Coniugato enzimatico / Conjugé enzymatique / Conjugado enzimático / Enzymkonjugat
<b>SN</b>	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	<b>CONTROL</b>	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
<b>TEST</b>	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät	<b>REAG</b>	Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	<b>BUF   WASH</b>	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
<b>BUF</b>	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
<b>CONJ</b>	Conjugate / Coniugato / Conjugé / Conjugado / Konjugat	<b>DIL   SPE</b>	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillon / Diluyente de muestra / Probenverdünnungspuffer
<b>SUBS</b>	Substrate / Substrato / Substrat / Substrato / Substrat	<b>BUF   WASH   20X</b>	Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
		<b>DET   REAG</b>	Detection Reagent / Reagente Diretta / Réactif de Détection / Reactivo de Detección / Nachweis Reagenz

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.