MSwab™

**Collection, Preservation and Transport System
Product Insert & How to Use Guide**

Copan MSwab™ Collection, Preservation and Transport System for molecular and culture applications - Product Insert & How To Use Guide

See symbol glossary at end of insert

INTENDED USE

MSwab™ is a Collection, Transport and Preservation System intended for the collection and transport of clinical specimens from the collection site to the testing laboratory. Once in the laboratory, it can be processed using standard operating procedure for:

- bacterial culture of gram positive aerobic and facultative anaerobic organisms
- viruses culture of HSV 1 and HSV 2
- nucleic acids detection of bacteria and viruses

SUMMARY AND PRINCIPLES

One of the routine procedures in the diagnosis of bacteriological infections involves the collection and safe transportation of swab samples. This can be accomplished using the Copan MSwab™ that is a Collection, Transport and Preservation System. The medium is designed to maintain the viability of Gram positive aerobic and facultative anaerobic bacteria and viruses during transit to the testing laboratory.

MSwab™ preserves nucleic acids of pathogenic agents to be detected with molecular amplification techniques (NAAT).

Copan MSwab™ Collection, Transport and Preservation System is supplied in three different formats:

- a) Collection kit format. Each collection kit consists of a package containing a plastic screw-cap tube with conical shaped bottom filled with 1 ml of MSwab™ transport and preservation medium and a small sterile peel pouch containing one specimen collection FLOQSwab™ that has a tip flocked with nylon fiber.
- b) Tube only format. A plastic screw-cap tube with conical shaped bottom filled with 1 ml or 2 ml of MSwab™ transport and preservation medium.
- c) Collection kit format with cleaning swab. Each collection kit consists of a package containing a plastic screw-cap tube with conical shaped bottom filled with 2 ml of MSwab™ transport and preservation medium and a small sterile peel pouch containing one specimen collection FLOQSwab™ that has a tip flocked with nylon fiber and a sterile peel pouch containing one cleaning swab with large tip in fiber wound for the removal of excess vaginal mucus.

Once a swab sample is collected, it should be placed immediately into the MSwab™ transport medium. Swab specimens for bacterial or viral investigations collected using MSwab™ should be transported directly to the laboratory, preferably within 2 hours of collection (1, 2, 7) to maintain optimum organism viability. If immediate delivery or processing is delayed, then specimens should be refrigerated at 4 – 8°C or stored at room temperature (20 – 25°C) and processed within 48 hours. Bacterial viability studies for *Staphylococcus aureus*, ATCC 29213 and ATCC 6538, and *Staphylococcus aureus* (Methicillin resistant) ATCC 43300 and ATCC 700698 show viability of the tested organisms up to 14 days at refrigerated temperature (4 – 8°C) or 72 hours at room temperature (20 – 25°C). Independent scientific studies on swab transport systems have shown that for certain bacteria viability is superior at refrigerated temperatures compared with room temperature (12 – 21). If viral samples must be frozen, they should be frozen at -70°C.

Swab specimens collected using MSwab™ for bacterial and viral nucleic acids investigations should be processed within 14 days when stored at room temperature (20 – 25°C), within 21 days when stored at 4°C and within 6 months when stored at -20°C.

REAGENTS

Formulation of MSwab™ Transport Medium

TRIS HCl
EDTA
TRIS Base
Dimethyl Sulfoxide (DMSO)
Bovine Serum Albumin
Distilled water

PRECAUTIONS

1. This product is For In Vitro Diagnostic Use.
2. Observe approved biohazard precautions and aseptic techniques. To be used only by adequately trained and qualified personnel.
3. All specimens and materials used to process them should be considered potentially infectious and handled in a manner which prevents infection of laboratory personnel. Sterilize all biohazard waste including specimens, containers and media after their use. Observe other CDC Biosafety Level 2 recommendations (31, 32, 33, 34).
4. Directions should be read and followed carefully.

STORAGE

This product is ready for use and no further preparation is necessary. The product should be stored in its original container at 5 – 25°C until used. Do not overheat. Do not incubate or freeze prior to use. Improper storage will result in a loss of efficacy. Do not use after expiration date, which is clearly printed on the outer box and on each individual collection unit and the specimen transport tube label.

PRODUCT DETERIORATION

Copan MSwab™ should not be used if (1) there is evidence of damage or contamination to the product, (2) there is evidence of leakage, (3) the expiration date has passed, (4) the swab package is open, or (5) there are other signs of deterioration.

SPECIMEN COLLECTION, STORAGE AND TRANSPORTATION

Specimens collected for microbiological investigations which comprise the isolation of bacteria or viruses should be collected and handled following published manuals and guidelines (7, 8, 4).

To maintain optimum organism viability and nucleic acids integrity, transport specimens collected using MSwab™ directly to the laboratory, preferably within 2 hours of collection (1, 2, 7). If immediate delivery or processing is delayed, then specimens should be refrigerated at 4 – 8°C or stored at room temperature (20 – 25°C) and processed within 48 hours. Bacterial viability studies for *Staphylococcus aureus*, ATCC 29213 and ATCC 6538, and *Staphylococcus aureus* (Methicillin resistant) ATCC 43300 and ATCC 700698 show viability of the tested organisms up to 14 days at refrigerated temperature (4 – 8°C) or 72 hours at room temperature (20 – 25°C). If viral samples must be frozen, they should be frozen at -70°C.

Specimens collected using MSwab™ to be analyzed with NAAT should be processed within 14 days when stored at room temperature (20 – 25°C), within 21 days when stored at 4°C and within 6 months when stored at -20°C.

Specific requirements for the shipment and handling of specimens should be in full compliance with state and federal regulations (34, 35, 36, 37). Shipping of spe-

cimens within medical institutions should comply with internal guidelines of the institution. All specimens should be processed as soon as they are received in the laboratory.

MATERIALS SUPPLIED

Fifty (50) MSwab™ collection units are contained in a shelf pack and 6 x 50 units are contained in a box.

Copan MSwab™ Collection, Transport and Preservation System is supplied in three different formats:

- Collection kit format. Each collection kit consists of a package containing a plastic screw-cap tube with conical shaped bottom filled with 1 ml of MSwab™ transport and preservation medium and a small sterile peel pouch containing one specimen collection FLOQSwab™ that has a tip flocked with nylon fiber.
- Tube only format. A plastic screw-cap tube with conical shaped bottom filled with 1 ml or 2 ml of MSwab™ transport and preservation medium.
- Collection kit format with cleaning swab. Each collection kit consists of a package containing a plastic screw-cap tube with conical shaped bottom filled with 2 ml of MSwab™ transport and preservation medium and a small sterile peel pouch containing one specimen collection FLOQSwab™ that has a tip flocked with nylon fiber and a sterile peel pouch containing one cleaning swab with large tip in fiber wound for the removal of excess vaginal mucus.

Two main types of collection applicators are available: regular size flocked nylon applicator intended for the collection of samples from nose, throat, vagina, rectum, feces or wounds; pernasal thin and flexible nylon flocked applicator for the collection of samples from the nasopharynx. All collection swab applicators provided with MSwab™ have a molded breakpoint in the shaft of the applicator which is highlighted with a colored indication line marked on the shaft of the applicator. After the sample is collected from the patient, the molded breakpoint facilitates easy breakage of the swab applicator into the MSwab™ tube of transport medium. MSwab™ tube caps have an internal molded design that is able to capture the swab shaft when it is broken off into the tube and the cap is closed (see Fig.2). The action of screwing the cap onto the tube moves the end of the broken swab shaft into a funnel shaped molded docking receptacle in the cap. This molded funnel shape effectively captures the end of the broken applicator shaft and secures it firmly in the dock by friction grip. In the testing laboratory when the swab cap is unscrewed and removed, the swab applicator is attached to the cap. This feature allows the operator to conveniently remove the swab from the transport tube and perform various microbiology analyses using the tube cap as a handle to hold the swab applicator.

Due to the flexibility of the shaft of pernasal swab (6E013N), the capture cap feature is not applicable, as the broken applicator may not firmly fit into the cap.

MATERIALS REQUIRED BUT NOT SUPPLIED

Appropriate materials for isolating and culturing aerobic and facultative anaerobic bacteria. These materials include culture media plates or tubes and incubation systems. Refer to laboratory reference manuals for recommended protocols for culture and identification techniques for aerobic and facultative anaerobic from clinical swab samples (2, 4).

Appropriate materials for isolating, differentiating and culturing viruses. These materials include tissue culture cell lines, tissue culture medium, incubation systems and reading equipment. Refer to appropriate references for recommended protocols for isolation and identification of viruses (1, 7).

Appropriate materials for nucleic acids extractions and amplifications for molecular assays. Heating block and centrifuge for rapid extraction method. Refer to laboratory reference manuals for nucleic acid detections and amplifications from clinical swab samples.

In case of rapid extraction method, it is suggested to transfer the appropriate aliquot of MSwab™ (see paragraph **Processing MSwab™ specimens for molecular testing in the laboratory point B**) in Copan tube with glass beads (code available separately 2E013S50)

DIRECTIONS FOR USE

Copan MSwab™ Collection, Transport and Preservation System is available in product configurations indicated in the table below.

Table 1

Catalog No.	Copan MSwab™ Product Descriptions	Pack Size
6E012N	Single use sample collection pack containing: - 12mm x 80 mm screw cap tube containing 1ml MSwab™ transport and preservation medium - One FLOQSwab™ applicator regular size with flocked nylon fiber tip, sterile and individually wrapped.	50 units/box 6x50units/case
6E013N	Single use sample collection pack containing: -12mm x 80 mm screw cap tube containing 1ml MSwab™ transport and preservation medium - One pernasal thin & flexible FLOQSwab™ applicator with flocked nylon fiber tip, sterile and individually wrapped.	50 units/box 6x50units/case
6E011N	Single transport and preservation tube: -1 ml MSwab™ transport and preservation medium in screw cap tube 12 x 80 mm with internal conical shape	50 units/box 6x50units/case
6U019N	Single transport and preservation tube: -2 ml MSwab™ transport and preservation medium in screw cap tube 12 x 80 mm with internal conical shape	50 units/box 6x50units/case
6E028N.MER	Single use sample collection pack containing: - 12mm x 80 mm screw cap tube containing 2 ml MSwab™ transport and preservation medium - One FLOQSwab™ applicator regular size with flocked nylon fiber tip, sterile and individually wrapped. -One swab with large tip in fiber wound, sterile and individually wrapped.	50 units/box 6x50units/case

Other product codes may be available. For updates please refer to our website: www.copanflock.com

Specimen Collection

Proper specimen collection from the patient is extremely critical for successful isolation and identification of infectious organisms. For specific guidance regarding specimen collection procedures, consult published reference manuals (7, 2).

MSwab™

HPC005 Rev.00 Date 2014.11

Do not use the MSwab™ medium for pre-moistening or pre-wetting the applicator swab prior to collecting the sample or for rinsing or irrigating the sampling sites.

For MSwab™ codes 6E012N e 6E013N:

1. Peel open the kit package and remove the tube of medium and inner pouch containing the sterile swab applicator (see Figure 1).
2. Remove the swab applicator from its peel pouch and use to collect the clinical specimen. The operator must touch the swab applicator only above the colored breakpoint line, as illustrated in Figure 1, which is the opposite end to the nylon fiber tip. At all times when handling the swab applicator, the operator must not touch the area below the colored breakpoint line (the area from the line to the tip of the nylon flocked swab) as this will lead to contamination of the applicator shaft and the subsequent culture.
3. Collect the sample from the patient.
4. Unscrew and remove the cap from MSwab™ tube making sure not to spill the medium.
5. Insert the swab into the tube until the red marked breaking point is at the level of the tube opening.
6. Bend and break the swab at the red marked breaking point holding the tube away from your face.
7. Discard the broken handle part of the swab shaft into an approved medical waste disposal container.
8. Replace cap on the tube and secure tightly.
9. Write patient information on the tube label or apply patient identification label. Send the sample to the test laboratory


For MSwab™ code 6E028N.MER:

1. Peel open the kit package and remove the tube of medium and the two inner pouches, one containing the sterile swab applicator and one containing the cleaning swab with large tip "Cleaning Swab".
2. The use of cleaning swab is recommended only for the endocervical specimen collection procedure. The cleaning swab, (made with large tip in fiber wound) – see Figure 1, is used only to remove the excess mucus present on the external cervical os and surrounding mucosa. **After use, discard the swab.** The large cleaning swab must not be used and must be discarded in case of specimen collection procedure other than endocervical.
3. For the collection of clinical specimens refer to the instructions reported for codes 6E012N and 6E013N from point 2 and subsequents.


For MSwab™ 6E011N e 6U019N:

1. After the swab sample is taken from the patient, unscrew and remove the cap from MSwab™ tube making sure not to spill the medium.
2. Insert the swab into the tube.
3. If the swab has a breaking point, bend and break the swab at the breaking point holding the tube away from your face.
4. Replace cap on the tube and secure tightly.
5. Write patient information on the tube label or apply patient identification label. Send the sample to the test laboratory


Fig 1. Collection swab showing breakpoint indication line and area for holding the applicator

COLLECTION SWAB

CLEANING SWAB

The operator must only handle the part of the swab applicator shaft above the breakpoint indication line as shown in Fig 1. After the swab sample is taken from the patient, the swab applicator shaft is broken off at the colored breakpoint indication line into the MSwab™ tube of transport medium. The operator then discards the handle part of the swab into an approved medical waste disposal container. The tube's screw cap is then replaced and secured tightly. The action of screwing the cap onto the tube moves the end of the broken swab shaft into a funnel shaped molded docking receptacle in the cap (see Fig 2). This molded funnel shape captures the end of the broken applicator shaft and secures it firmly in the dock by friction grip.

Fig 2. Capture of broken swab applicator stick by MSwab™ tube cap

Capture Cap Feature is guaranteed only when using Copan Regular Size Flock Swab.

In the testing laboratory when the MSwab™ cap is unscrewed and removed, the swab applicator stick is securely attached to the cap. This feature allows the operator to conveniently remove the swab and perform various microbiology analyses using the tube cap as a handle to hold and manipulate the swab.

NOTE ABOUT CAPTURE CAP FEATURE:

The tube cap is designed to have capturing feature when COPAN regular size flock swab is used and broken inside the tube. If pernasal flock swab is used and broken into the tube then the capturing feature is not guaranteed.

Pernasal thin and flexible swab (6E013N) must be carefully removed from the tube immediately before proceeding with sample analysis. This precaution is intended to avoid accidental falling of the swab from the cap

Processing MSwab™ Specimens in the Laboratory – Bacteriology

MSwab™ samples should be processed for bacteriological culture using recommended culture media and laboratory techniques which will depend on the specimen type and the organism under investigation. For recommended culture media and techniques for the isolation and identification of bacteria from clinical swab specimens refer to published microbiology manuals and guidelines (1-6).

Culture investigations of swab specimens for the presence of bacteria routinely involve the use of solid agar culture medium in Petri dish plates. The procedure for inoculation of MSwab™ samples onto solid agar in Petri dishes is as follows.

Note: Wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe other CDC Biosafety Level 2 recommendations (31, 32, 33, 34).

Vortex mix the MSwab™ tube containing the swab sample for 5 seconds to release the sample from the swab tip and evenly disperse and suspend the patient specimen in the medium

1. Unscrew the MSwab™ cap and remove the swab applicator.
2. Roll the tip of the MSwab™ applicator onto the surface of one quadrant of the culture media plate to provide the primary inoculum.
3. If it is necessary to culture the swab specimen onto a second culture media plate, return the MSwab™ applicator to the transport medium tube for two seconds to absorb and recharge the applicator tip with transport medium/patient sample suspension then repeat Step No. 3.
4. If it is necessary to inoculate additional culture media plates, return the MSwab™ applicator to the transport medium tube and recharge the swab applicator tip with the transport medium/patient sample suspension before inoculating each additional plate.

The procedure described above utilizes the MSwab™ applicator like an inoculation wand to transfer the suspension of patient sample in transport medium onto the surface of a culture plate creating the primary inoculum (see Fig 3).

Alternatively, the operator can vortex mix the MSwab™ tube with the swab inside for 5 seconds and then transfer 100µl volumes of the suspension onto each culture plate using a volumetric pipettor and sterile pipet tips. Standard laboratory techniques should then be used to streak the primary inoculum of patient sample across the surface of the culture plate (see Fig 4).

Fig 3. Procedures for inoculation of MSwab™ specimens onto solid agar in Petri dishes


1. Using swab to inoculate specimen

2. Using pipettor and sterile pipet tips to inoculate 100µl of specimen
This inoculation procedure should be considered preferable if the clinical specimen is also processed with NAAT assays.

Fig 4. Procedure for streaking MSwab™ specimens on agar Petri dishes for primary isolation (33)

Seed a primary inoculum of MSwab™ specimen onto the surface of an appropriate agar culture plate in the first quadrant.

Use a sterile bacteriology loop to streak the primary inoculum across the surface of the second, third and fourth quadrants of the agar culture plate.

Preparation of Gram Stain Smears of MSwab™ Specimens

Laboratory analysis of clinical swab samples collected from certain sites on the patient can routinely include microscopic examination of stained preparations ("direct smears") using the Gram stain procedure. This can provide valuable information to physicians who are managing patients with infectious diseases (22). There are many instances in which a Gram stain can assist in making a diagnosis (23, 27).

The Gram stain can also help to judge specimen quality and contribute to the selection of culture media especially with mixed flora.

Microscope slides of patient specimens transported in Copan MSwab™ transport system can be prepared for Gram stain analysis, as describe below, by sampling an aliquot of vortexed suspension of the swab (3, 4). Sample transported in MSwab™ elution medium represent a homogeneous suspension in liquid phase. It can be uniformly smeared allowing clear and easy reading.

Note: Wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe other CDC Biosafety Level 2

recommendations (31, 32, 33, 34).

1. Take a clean glass microscope slide, place it on a flat surface and inscribe an area using a diamond-tipped or similar glass marker to identify the location of the specimen inoculum. Note: a slide with a pre-marked 20 mm well can be used.

2. Vortex mix the MSwab™ tube containing the swab sample for 5 seconds to release the sample from the swab tip and evenly disperse and suspend the patient specimen in the medium.

3. Unscrew the MSwab™ cap and using a sterile pipet, transfer 1 – 2 drops of sample suspension to the inscribed area on the glass slide. Note: about 30µl would be a suitable amount of liquid for a pre-marked 20 mm diameter well slide.

In case of bloody or thicker specimens particular care should be taken to thinly spread the sample on the slide. Bacteria are difficult to detect if the sample shows many red cells and debris.

4. Allow the specimen on the slide to air dry at room temperature or place the slide in an electric slide warmer or incubator set at a temperature not exceeding 42°C.

5. Fix smears using methanol. Methanol fixation is recommended as it prevents lysis of Red Blood Cells, avoids damage to all host cells and results in a cleaner background (3, 4, 22).

6. Follow published laboratory reference manuals and guidelines for performing the Gram stain. If commercial Gram stain reagents are used, it is important to comply with instructions in the manufacturer's product insert for performance test procedure.

For further information or guidance on the preparation of specimen slides for microscopic analysis, for information on Gram staining procedures and the interpretation and reporting of microscopic analysis, consult published laboratory reference manuals (1 - 5, 22 - 27).

Processing MSwab™ Specimens in the Laboratory – Virology

The survival of HSV 1 and HSV 2 depends on many factors including the type and concentration of the microorganism, duration of transport and storage temperature. To maintain optimum viability, specimens should be transported directly to the laboratory preferably within 2 hours of collection (1, 2, 7, 29). If immediate delivery or processing is delayed, then specimens collected using Copan MSwab™ Collection, Transport and Preservation System should be refrigerated at 4-8°C or stored at room temperature (20-25°C) and processed within 48 hours. If samples must be frozen, they should be frozen placed at -70°C.

In simulated transport and storage studies, the Copan MSwab™ System was shown to be capable of maintaining the viability of HSV 1 and HSV 2 at refrigerated (4-8°C) and room temperature (20-25°C) conditions for up to 48 hours. Based on performance studies conducted by Copan and independent scientific publications, viability of certain microorganisms is superior at refrigerated temperatures compared with room temperature (12 – 21, 29).

MSwab™ specimens should be processed for virology culture using recommended cells line and laboratory techniques that will depend on the specimen type and the organism under investigation. For recommended shell vials and techniques for the isolation and identification of HSV 1 and HSV 2 from clinical swab specimens, refer to published virology manuals and guidelines (1 – 6, 29, 30).

Culture investigations of swab specimens for the presence of HSV 1 and HSV 2 routinely involve the use of cells culture in shell vials. The procedure for inoculation of MSwab™ specimens onto shell vials is described below.

NOTE: Wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe other BSL 2 recommendations.

1. Shake using a vortex mixer for 5 seconds the MSwab™ tube containing the swab specimen to release the specimen from the swab tip and evenly disperse and suspend the patient specimen in the liquid transport medium.
2. Unscrew the MSwab™ cap and remove the swab applicator.
3. Transfer 200 µl volumes of the suspension into the shell vial and proceed as per laboratory internal procedure.
4. **NOTE:** Patient specimen that may contain a high load of bacterial contaminants may require additional antibiotics into the refeed medium.

Processing MSwab™ specimens for molecular testing in the laboratory.

Specimens received in the laboratory for nucleic acid detection should be processed when received in the laboratory. In case of delay, please refer to the appropriate specimen storage conditions.

NOTE: Wear gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe other CDC Biosafety Level 2 recommendations (31, 32, 33, and 34).

When working with molecular methods care should be taken to prevent carry over contamination. Spatial separation of working areas and unidirectional workflow are essential to prevent amplicon carry-over. (35)

MSwab™ medium formulation has been developed to be compatible with PCR Master Mixes; this makes it suitable for direct NAAT assays without the need of a standard extraction and purification step.

Two different methods (A and B) can be applied to MSwab™ sample processing :

A) Standard extraction method:

1. Mix MSwab™ tube with a vortex mixer for 10 seconds, unscrew the cap and using it as a handle, spin the cap in between the thumb and index finger to drain most of the fluid from the tip.

NOTE: Due to the flexibility of the shaft of pernasal swabs (6E013N), the capture cap feature is not applicable, as the broken applicator may not firmly fit into the cap. In this case, use tweezers to extract the applicator from the tube after spinning the swab tip to extract most of the fluid from the tip.

2. Discard the swab, and transfer the appropriate amount of sample (e.g. 200ul) to an extraction tube as per laboratory SOP.
3. Continue as per extraction and amplification kits procedures

The MSwab™ has been tested with the following extraction methods: Silica-gel membrane (Qiagen, Macherey & Nagel) and Magnetic beads (EasyMAG, Qiasymphony, Magnapure). Other extraction methods may be also applicable prior validation.

B) Rapid extraction method:

A thermal shock step is suggested to lyse samples for viral detection.

1. Vortex original tube of MSwab™ sample for 10 seconds.

NOTE: Due to the flexibility of the shaft of pernasal swabs (code 6E013N), the capture cap feature is not applicable, as the broken applicator may not firmly fit into the cap. In this case, use tweezers to remove the swab from the tube.

2. Using a micropipette, transfer 200 µl to a sterile microtube containing glass beads (2E013S50: Copan Tube with glass beads not supplied with the kit) and vortex 10s. Store the original MSwab™ sample for culture or additional testing.
3. Heat in a thermostat set at 98°-100°C for at least:
 - 3 min FOR VIRUSES
 - 10 min FOR BACTERIA GRAM POSITIVE
 - 5 min FOR BACTERIA GRAM NEGATIVE
4. Cool the tube at Room Temperature for 5 min.
5. Centrifuge the microtube at 10000xg for 2 min to settle cells debris.
6. Transfer an aliquot of the heat treated MSwab™ sample to the master mix and amplify as per amplification kit procedure.

Rapid extraction method as described above has been tested with R-Biopharm RIDA®GENE Real Time PCR kits, with Seegene Anyplex FluA/B Typing Real Time detection kit and with Genesig Real Time PCR kits by Primer Design in comparison with standard extraction method. Other PCR amplification kit may be also applicable prior validation. Tested strains include : B. pertussis (ATCC 8467), S. aureus (methicillin resistant) (ATCC 43300), E. coli O157-H7 (ATCC 700728), S. tiphimurium (ATCC 14028), C. difficile (ATCC 9689), Influenza A virus (ATCC VR-822), Influenza B virus (ATCC VR 786), Neisseria gonorrhoeae (ATCC 43069) and Chlamydia trachomatis (ATCC VR880). Tests were not conducting using human matrixes. Results obtained will largely depend on proper and adequate specimen collection, as well as timely transport and processing in the laboratory.

Results obtained will largely depend on proper and adequate specimen collection, as well as timely transport and processing in the laboratory.

QUALITY CONTROL

MSwab™ applicators are tested to ensure they are non-toxic to bacteria. MSwab™ medium and applicators are tested to ensure they are non-toxic to cell lines used for HSV 1 and HSV 2 cultivation. MSwab™ transport medium is tested for pH stability and for absence of inhibitors with direct nucleic acid amplification (9). MSwab™ is quality control tested before release for ability to maintain viable Gram positive aerobic and facultative anaerobic bacteria and HSV viruses at room temperature (20 – 25°C) for specified time points. Procedures for quality control of microbiology transport devices should be conducted using testing methods described in Clinical and Laboratory Standards Institute M40-A and other publications (9). If aberrant quality control results are noted, patient results should not be reported.

LIMITATIONS

1. In the laboratory, wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe other CDC Biosafety Level 2 recommendations (31, 32, 33, 34) when handling or analyzing patient samples.
2. Condition, timing and volume of specimen collected for culture are significant variables in obtaining reliable culture results. Follow recommended guidelines for specimen collection (7, 8, 4).
3. MSwab™ is intended for use as a collection and transport medium for culture of gram positive aerobic and facultative anaerobic bacteria, HSV 1 and HSV 2 viruses and for nucleic acids detection of bacteria and viruses. MSwab™ cannot be used as enrichment, selective or differential medium.
4. MSwab™ is an antibiotics free medium. Patient specimen that may contain a high load of bacterial contaminants may require additional antibiotics into the refeed medium.
5. Performance testing with Copan MSwab™ was conducted using laboratory strains spiked onto a swab following the test protocols described in Clinical and Laboratory Standards Institute M40-A Approved Standard (9). Performance testing was not conducted using human specimens.
6. Performance testing with Copan MSwab™ was conducted using Copan flocked swabs.

WARNINGS

1. Not for diagnostic use.
2. Do not re-sterilize unused swabs.
3. ⚠ This product is for single use only; reuse may cause a risk of infection and/or inaccurate results.
4. Do not re-pack.
5. Not suitable for any other application than intended use.
6. The use of this product in association with a rapid diagnostic kit or with diagnostic instrumentation should be previously validated by the user.
7. Do not use if the swab is visibly damaged (i.e., if the swab tip or swab shaft is broken).
8. Do not use excessive force or pressure when collecting swab samples from patients as this may result in breakage of the swab shaft.
9. Do not use the same tube for more than one patient. This causes wrong or bad diagnosis.
10. Before transporting, make sure MSwab™ screw cap is tightly closed.
11. Dispose of unused reagents, waste and specimens in accordance with local regulations.
12. Avoid contact of MSwab™ medium with skin, mucous membranes. If contact does occurs, immediately wash with a large amounts of water.
13. Do not ingest the medium.
14. Directions for use must be followed carefully. The manufacturer cannot be held responsible for any unauthorized or unqualified use of the product.
15. To be handled by trained personnel only.
16. It must be assumed that all specimens contain infectious micro-organisms; therefore all specimens must be handled with appropriate precautions. After use, tubes and swabs must be disposed of according to laboratory regulations for infectious waste. Observe CDC Biosafety Level 2 recommendations (31, 32, 33, 34).
17. Do not use the MSwab™ medium for pre-moistening or pre-wetting the applicator swab prior to collecting the sample or for rinsing or irrigating the sampling sites.

RESULTS

Results obtained will largely depend on proper and adequate specimen collection, as well as timely transport and processing in the laboratory.

PERFORMANCE CHARACTERISTICS

Bacterial recovery

The test procedures employed for determining bacterial viability performance were based upon the quality control methods described in Clinical and Laboratory Standards Institute M40-A (9).

MSwab™ system has an intended use limited to Gram positive aerobic and facultative anaerobic bacteria and HSV 1 and HSV 2, therefore its field application is more restricted than some other devices. For this reason the bacterial recovery studies were conducted under the simulated transport and storage conditions as described and defined in CLSI M40-A, Quality Control of Microbiological Transport Systems: Approved Standard and included only the Gram positive aerobic and facultative anaerobic strains from the Group 1 of paragraph 7.11.1 of the CLSI M40-A document, in particular:

<i>Streptococcus pyogenes</i>	ATCC® 19615
<i>Streptococcus pneumoniae</i>	ATCC® 6305
<i>Pseudomonas aeruginosa</i>	ATCC® BAA-427

In addition, Copan included testing of additional Gram positive aerobic and facultative anaerobic microorganisms, clinically relevant, not required by CLSI M40-A.

The specific bacterial strains used in these studies are here listed:

<i>Enterococcus faecalis</i>	ATCC® 29212
<i>Staphylococcus epidermidis</i>	ATCC® 12228
<i>Staphylococcus aureus</i>	ATCC® 29213
<i>Streptococcus agalactiae</i> (Group B Strep)	ATCC® 13813
<i>Kocuria rhizophila</i>	ATCC® 9341
<i>Listeria monocytogenes</i>	ATCC® 19114
<i>Bacillus cereus</i>	ATCC® 10876
<i>Staphylococcus aureus</i> (Methicillin resistant)	ATCC® 43300
<i>Staphylococcus aureus</i>	ATCC® 6538
<i>Staphylococcus aureus</i> (Methicillin resistant)	ATCC® 700698
<i>Streptococcus pneumoniae</i>	ATCC® 49136

All cultures of bacteria were ATCC (American Type Culture Collection) and were obtained commercially.

The selection of these organisms also reflects those Gram positive aerobic and facultative anaerobic bacteria that would normally be encountered in specimens collected and analyzed in a typical clinical microbiology laboratory.

Bacterial viability studies were performed on the Copan MSwab™ at two different ranges of temperature, 4 – 8°C and 20 – 25°C, corresponding to refrigerator and room temperature, respectively. Swabs accompanying each transport system were inoculated in triplicate with 100µl of specific concentrations of organism suspension. Swabs were then placed in their respective transport medium tubes and were held for 0 hrs, 24 hrs and 48 hrs. At the appropriate time intervals, each swab was processed according to the Roll-Plate or Swab Elution Method.

Additional bacterial viability studies for *Staphylococcus aureus*, ATCC 29213 and ATCC 6538, and *Staphylococcus aureus* (Methicillin resistant) ATCC 43300 and ATCC 700698 were performed on the Copan MSwab™ at two different ranges of temperature, 4 – 8°C and 20 – 25°C, corresponding to refrigerator and room temperature, respectively.

Swab accompanying the transport system was inoculated in triplicate with 100µl of specific concentrations of organism suspension.

Swabs were then placed in their respective transport medium tubes and:

For studies performed at 4 – 8°C inoculated MSwab™ tubes were held for 0 hrs, 10 days and 14 days. At the appropriate time intervals, each MSwab™ was processed according to the Roll-Plate Method.

For studies performed at 20 – 25°C inoculated MSwab™ tubes were held for 0 hrs and 72 hrs. At the appropriate time intervals, each MSwab™ was processed according to the Roll-Plate Method.

Bacterial overgrowth studies were performed on the Copan MSwab™ at 4 – 8°C, corresponding to refrigerator temperature. Swabs accompanying each transport system were inoculated in triplicate with 100µl of specific concentrations of organism suspension. Swabs were then placed in their respective transport medium tubes and were held for 0 hrs and 48 hrs. At the appropriate time intervals, each swab was processed according to the Roll-Plate Method.

Bacterial overgrowth studies were performed using *Pseudomonas aeruginosa*.

Viral viability studies were performed using HSV 1 and HSV 2. Swabs accompanying each transport system were directly inoculated in triplicate with 100µl of organisms suspension. Swabs were then placed in their respective transport medium tubes and were held for 0, 24 and 48 hours at both 4°C and room temperature (20-25°C). At the appropriate time interval, each swab was vortexed, removed from its transport medium tube and then 200µl aliquots of this suspension was inoculated into shell vials. All cultures were processed by standard laboratory culture technique and examined after a specified incubation time. Organism viability was determined by fluorescing foci counts.

Organisms evaluated were:

Herpes Simplex Virus Type 1 (HSV 1)	ATCC VR-539
Herpes Simplex Virus Type 2 (HSV 2)	ATCC VR-734.

In accordance with Clinical and Laboratory Standards Institute M40-A, viability performance is measured for each test organism at the 48 hrs time point and compared with the acceptance criteria.

In both the Roll-Plate and Swab Elution viability performance studies, Copan MSwab™ System was able to maintain acceptable recovery of all organisms evaluated at both refrigerator (4 – 8°C) and room temperature (20 – 25°C). Acceptable recovery for the Roll-Plate Method is defined as ≥ 5 CFU following the specified holding time from the specific dilution that yielded zero-time plate counts closest to 300 CFU. Acceptable recovery for the Swab Elution Method is defined as no more than a $3 \log_{10}$ ($1 \times 10^3 \pm 10\%$) decline in CFU between the zero-time CFU count and the CFU of the swabs after the specified holding time.

Additional time-points were tested for *Staphylococcus aureus*, ATCC 29213 and ATCC 6538, and *Staphylococcus aureus* (Methicillin resistant) ATCC 43300 and ATCC 700698.

In the Roll-Plate viability performance studies, Copan MSwab™ System was able to maintain acceptable recovery of all organisms evaluated at both refrigerator (4 – 8°C) for 14 days and room temperature (20 – 25°C) for 72 hrs. Acceptable recovery for the Roll-Plate Method is defined as ≥ 5 CFU following the specified holding time from the specific dilution that yielded zero-time plate counts closest to 300 CFU.

Viability performance studies also include an assessment of bacterial overgrowth at refrigerated temperatures (4 – 8°C). For the Swab Elution Method, an overgrowth assessment is made on all bacteria species tested at the 48 hrs holding time point. Overgrowth assessment using the Swab Elution Method is defined as

MSwab™

HPC005 Rev.00 Date 2014.11

greater than 1 log₁₀ increase in CFU between the zero-time CFU count and the holding time point. For the Roll-Plate Method, an overgrowth assessment is made with a separate analysis in which swabs are dosed with 100µl containing 10² CFU of *Pseudomonas aeruginosa* culture. Overgrowth under these conditions is defined as greater than 1 log₁₀ increase in CFU between zero-time CFU and the 48 hrs holding time point.

Copan MSwab™ System demonstrated no overgrowth based on the acceptance criteria described in Clinical and Laboratory Standards Institute M40-A.

Copan MSwab™ System was able to maintain the viability of the following organisms for at least 48 hours at both room temperature (20-25°C) and in the refrigerator (2-8°C) under the test conditions described above: Herpes Simplex Virus Type 1, Herpes Simplex Virus Type 2.

Nucleic Acids preservation

Specimens collected using MSwab™ to be analyzed with NAAT should be processed within 14 days when stored at room temperature (20 – 25°C), within 21 days when stored at 4°C and within 6 months when stored at -20°C.

Performance testing with Copan MSwab™ was conducted using laboratory strains suspensions spiked onto the FLOQSwab associated to the medium. Tested strains include: *B. pertussis* ATCC 8467, *S. aureus* (methicillin resistant) ATCC 43300, *E. coli* O157-H7 ATCC 700728, *S. typhimurium* ATCC 14028, *C. difficile* ATCC 9689, Influenza A virus ATCC VR-822, Influenza B virus ATCC VR 786, *Neisseria gonorrhoeae* ATCC 43069 and *Chlamydia trachomatis* ATCC VR880. Performance testing was not conducted using human specimens or human matrixes.

Results obtained will largely depend on proper and adequate specimen collection, as well as timely transport and processing in the laboratory.

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Copan Flock Technologies Srl
Copan Diagnostics Inc.