Para-Pak® SpinCon®
Stool Concentration System
U.S. Patent No. 5,882,943

REF 972200

IVD

Rx Only

INTENDED USE
Para-Pak SpinCon is a unique, patented system for the concentration of parasite eggs, larvae, and protozoa from preserved fecal specimens. Specimens preserved in 10% Formalin, Sodium Acetate Formalin (SAF) and EcoFix may be used with the system. Kit systems are designed for easy use by individuals not trained in microbiologic procedures.

SUMMARY AND EXPLANATION OF THE TEST
Diagnosis of intestinal parasitic disease is confirmed by the recovery and identification of helminth eggs and larvae, or protozoan trophozoites and cysts. Concentration procedures facilitate this process and are of importance when small numbers of parasites are present in larger volumes of feces.¹⁻¹⁰

Traditional concentration procedures utilize gauze and funnel, or a similar preliminary filtration device followed by chemical treatment and centrifugation of the preserved stool. Often toxic chemicals such as diethyl ether, ethyl acetate and formalin are used in the concentration procedure.

The SpinCon procedure utilizes a device that concentrates parasites without the need for organic solvent treatment of the preserved stool specimen, thereby eliminating the hazard of working with and disposing of chemicals such as diethyl ether, ethyl acetate and formalin.

BIOLOGICAL PRINCIPLES
Para-Pak SpinCon involves passing a surfactant treated, preserved stool specimen through a preliminary screen by gravity flow. The surfactant helps to break down fecal aggregates, thus freeing parasites. The specimen is then forced by centrifugation through a series of two screens with successively smaller mesh. The series of screens, not present in other devices, trap stool debris, yet allow even the larger parasites to pass through. This second filtration step eliminates the need for organic solvent extraction of the stool specimen. The resulting pellet may be examined for the presence of parasites by standard wet mount procedures.

REAGENTS/MATERIALS PROVIDED
The maximum number of tests obtained from this test kit is listed on the outer box.
SpinCon Filtration Devices (item #CS3814A)
SpinCon Tubes with Caps (item #CS3813)
Para-Pak Surfactant (item #9601-02)

- MATERIALS NOT PROVIDED**
1. Microscope Slides and Cover Slips
 2. Physiological Saline
 3. Lugol's Iodine
 4. 10% Formalin, SAF or EcoFix
 5. Wooden Applicator Sticks
 6. Centrifuge
 7. Microscope

- PRECAUTIONS**
1. All reagents are for in vitro diagnostic use only.
 2. Patient specimens may contain HIV or other infectious agents and should be handled by properly trained personnel and disposed of as potential biohazards. Wear disposable gloves while handling specimens and performing test procedures.
 3. After centrifugation of surfactant treated stool in the SpinCon device, little to no liquid supernatant should be visible in the upper half of the device. Make sure proper quantity of surfactant is added to each stool, followed by adequate mixing prior to processing with the SpinCon device. Use only the surfactant provided with the kit.
 4. Any serious incident that has occurred in relation to the device should be reported to Meridian Bioscience, Inc., 3471 River Hills Drive, Cincinnati, Ohio 45244 USA or Technical Support Center 800-343-3858 and competent authority of the EU Member State in which the clinician and/or patient is established.
 5. **IMPORTANT:** See SDS for additional safety and hazard information.

HAZARD AND PRECAUTIONARY STATEMENTS

 Reagent A (MucoPenX), CON-Trate Reagent A and Macro-CON Surfactant	Signal word Warning Hazard Statements H302 - Harmful if swallowed Contains Polyethylene glycol 1,1,3,3-tetramethylbutylphenyl ether
---	---

SHELF LIFE AND STORAGE
Store all reagents and components at room temperature (15-30 C). Expiration dating is indicated on outer packaging.

- TEST PROCEDURE**
- A. Specimen Collection**
1. Para-Pak SpinCon is designed specifically for use with products from the Para-Pak line of parasitology products (10% Formalin, SAF, and EcoFix). Collect the specimens as directed by the appropriate Para-Pak package insert supplied with the collection system.
 2. The specimen must be well mixed in respective preservative and must be allowed to stand for a minimum of 30 minutes at room temperature prior to processing to allow adequate fixation.
 3. If a permanently stained smear is desired (Para-Pak SAF), it is recommended that a portion of preserved stool be handled separately prior to concentration with Para-Pak according to the respective package insert.
- B. Specimen Processing**
1. Add 2 drops Para-Pak Surfactant per 5 mL preserved stool specimen (8 to 10 drops per Para-Pak vial). Replace cap and shake vigorously for 60 seconds. Break-up of some formed or very mucoid specimens may be facilitated by vortexing.
 2. Make sure SpinCon device is assembled correctly with the upper and lower halves of the tube securely together and the primary funnel screen fitted securely in place in the upper half of the tube.
NOTE FOR Para-Pak ULTRA USERS: Remove tab from cap with key, "snap fit" SpinCon into top portion of ULTRA cap, invert and squeeze approximately 3 mL of specimen into the upper half portion of the tube. **PROCEED TO STEP #4.**
 3. Filter approximately 3 mL surfactant treated, preserved stool through the preliminary screen into the upper half of the SpinCon device (to the 3 mL mark on the tube). If necessary, use applicator sticks to help facilitate filtration of the specimens. Discard the preliminary screen into suitable biohazard receptacle.
 4. Add approximately 2 mL physiological saline to the upper half of the SpinCon device (to the 5 mL mark on the tube) for a total of approximately 5 mL diluted, filtered stool. Place cap into top of device securely and mix well by inverting the tube a few times.
 5. Centrifuge the device for ten minutes at 500 xg (1800-2200 rpm for most table top centrifuges).
 6. Discard top half of SpinCon device into a suitable biohazard receptacle. Decant supernatant from bottom half of device into a suitable biohazard receptacle.
 7. Resuspend final pellet in appropriate volume of 10% Formalin or SAF (EcoFix can also be used but is compatible for saline wet mounts only, not iodine).
 8. Perform a wet mount on the resuspended sediment in saline or Lugol's iodine and examine microscopically for the presence of parasites. (The proper ratio of sediment to saline or iodine is achieved when newsprint is legible through the prepared slide.)
- C. SpinCon Method for use only with EcoFix**
Follow the directions in Steps #1-6 under Section B, Specimen Processing. The following modification to the procedure should be followed after Step #6:
1. Prepare a smear (from the final pellet in Step #6) for permanent staining as follows:
 - a. Add a small amount of the suspended sediment to the slide.
 - b. Spread the mixture over the slide to produce a film which varies in thickness.
 - c. Allow the slide to dry overnight at room temperature, for several hours at 37 C, or under a laminar flow hood for 30-60 minutes.
 - d. Once slides prepared in EcoFix have dried completely, staining with EcoStain may begin. Refer to the package insert for Meridian catalog #801400.
 2. If using the MERIFLUOR C/G immunofluorescent product, use the same pellet from Step #6 above to prepare slides. Refer to the package insert for Meridian catalog #250050, Procedure Step #1.
 3. Resuspend the remaining pellet from Step #6 in appropriate volume of 10% formalin or SAF or saline. (EcoFix may be used, but is compatible for saline wet mounts only, not iodine.)
 4. Perform wet mount on the resuspended sediment in saline or Lugol's iodine and examine microscopically for the presence of parasites.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

1. SpinCon device should be free of cracks, with two screens fixed into top half of device. The two halves of the device should fit securely together. All screens should be uniform in appearance.
2. The surfactant should be free of observable bacterial or fungal contamination.

If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

SPECIFIC PERFORMANCE CHARACTERISTICS

1. The Para-Pak SpinCon, when used as directed (under Test Procedure) during clinical evaluation, was found to provide results equivalent to the standard Ritchie-Formalin ether (Ethyl Acetate) sedimentation concentration procedure and the Para-Pak CON-Trate System.
2. Listed below are the parasites recovered:

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REFERENCES

- ASMT, Recommended Procedures for the Examination of Clinical Specimens Submitted for the Diagnosis of Parasitic Infections. Am. J. Med. Technol., 1978;44:1101-1106.
- American Society of Parasitologist, Procedures Suggested for Use in Examination of Clinical Specimens to Parasitic Infection. J. Parasitol., 1977;63:959-960.
- Garcia, L.S. and Ash, L.R. "Diagnostic Parasitology", 2nd Edition, the C. B. Mosby Co. 1979;1-18.
- Garcia, L.S. and Shimizu, R. Comparison of Clinical Results for the use of Ethyl Acetate and Diethyl Ether in the Formalin-Ether Sedimentation Technique Performed on Polyvinyl Alcohol Preserved Specimens. J. Clin. Microbiol., 1981;13:709-713.
- Garcia, L.S. and Voge, M. Diagnostic Clinical Parasitology: I. Proper Specimen Collection and Processing. Am. J. Med. Technol., 1980;46:459-466.
- Maldonado, J.F., Acosta-Matienzo, J. and VelezHerrera, F. Comparative Value of Fecal Examination Procedures in the Diagnosis of Helminth Infections. Exptl. Parasitol., 1954;3:403-416.
- Melvin, D.M. and Brooke, M.M. "Laboratory Procedures for the Diagnosis of Intestinal Parasites", U.S.D.H.E.W., 1980;79:8282. CDC Atlanta, GA, 23-65.
- Ritchie, L.S., An Ether Sedimentation Technique for Routine Stool Examination. Bull. U.S. Army Med. Dept., 1948;8:326.
- Rohrbacher, G.H., The Recovery of Nematode Larvae by Baermann Apparatus as Affected by a Detergent. Proc. Helminth. Soc., 1957;24:24-25.
- Yang, J. and Scholten, Th. A Fixative for Intestinal Parasites Permitting the Use of Concentration and Permanent Staining Procedures. Am. J. Clin. Pathol., 1977;67:300-304.

SN10605

REV. 03/23

 Manufactured By	<p>Meridian Bioscience, Inc. 3471 River Hills Drive Cincinnati, OHIO - 45244 USA www.meridianbioscience.com</p> <p><u>Contacts:</u> Main Telephone (+1) 513.271.3700 Customer Service/Orders 800.543.1980 Technical Support Center 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124 E-mail: info@meridianbioscience.com</p>
 Authorized Representative	<p>Meridian Bioscience Europe, SRL Via Dell'Industria 7, 20035 Villa Cortese (Milano) ITALY www.meridianbioscience.com</p> <p><u>Contacts:</u> Main Telephone (+39) 0331.433636 E-mail: info@meridianbioscience.eu Technical Support: MBE-TechService@meridianbioscience.eu Customer Service/Orders:</p> <ul style="list-style-type: none"> For Italian Customers: ordini@meridianbioscience.com For Distributors / International Customers: Export.CustomerService@meridianbioscience.eu
<p>UK Authorised Representative</p>	<p>Launch Diagnostics Ash House Ash Road Longfield DA3 8JD UK</p>
 Authorized Representative	<p>MedEnvoy Switzerland Gotthardstrasse 28 6302 Zug Switzerland</p>

INTERNATIONAL SYMBOL CHART

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide es symboles, Guía de símbolos, Zeichenerklärung)

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For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.