


illumigene® Assay Workstation

Optimizing flexibility and efficiency for your molecular assay processing.



- **Optimized flexibility**

- Configurable to small or large batch sizes
- Adjustable to left- or right-handed operators


- **Enhanced efficiency**

- Simple tracking of samples throughout workflow
- *illumigene*® device transport rack matches *illumipro-10*™ configuration

- **Incredible ease of use**

- Clean/decontaminate through simple bleach submersion
- Quick assembly/disassembly/reconfiguration of all parts

illumigene[®] Molecular Diagnostic System



Flexible interlocking design allows for adaptation to run sizes as well as left- or right-handed operation.



Flow through wells for ease of cleaning and draining for optimized contamination control.



Interlocking stackable design to facilitate storage and maintain low footprint of **illumigene[®]** system.



Separate sample transport block matched with **illumipro-10[™]** design to enhance ease of use and tracking of patient samples.



The **illumigene[®]** workstation is designed to enhance the simplicity, efficiency, and workflow of the **illumigene[®]** assay system. Another innovation from the unparalleled global leader in medical diagnostics.

Ordering information

illumigene[®] Workstation.....610169



For more information, contact an **illumigene[®]** specialist at 1-888-763-6769 or visit us on the web at www.meridianbioscience.com.