




ImmunoCard®
MYCOPLASMA



The only rapid test that detects IgM antibody exclusively – for diagnosis of early/acute illness rather than convalescent disease.

“The Meridian ImmunoCard is the test of choice, especially for children, if only acute phase sera are available for testing. ¹”

“Immunocard Mycoplasma IgM assay is a valid and simple procedure which can reduce technologist time and turnaround time for labs analyzing small numbers of specimens. ²”


- Four simple steps to process
- No calculations necessary
- Rapid results in 9 minutes
- Sensitivity: 88% Specificity: 90%³

¹Talkington, D, Thaker, WL.. *Analysis of Complement Fixation and Commercial Enzyme Immunoassays for Detection of Antibodies to Mycoplasma pneumoniae in Human Serum*. Clinical and Diagnostic Laboratory Immunology, Sept 2000; 778-780.


²Alexander, et al. "Performance of Meridian ImmunoCard® Mycoplasma Test in a Multicenter Clinical Trial." *Journal of Clinical Microbiology*, May 1996; 1180-1183.

³ImmunoCard® Mycoplasma Package Insert.

Test Procedure



Interpretation of Results



*This chart does not contain complete instructions for use. For further information, please refer to the package insert.

Ordering Information

ImmunoCard® Mycoplasma

Test Size

30 Tests

Catalog #

709030

Other Meridian Products

ImmunoCard STAT!® Rotavirus

Test Size

30 Tests

Catalog #

750030

ImmunoCard STAT!® Strep A

30 Tests

750230

ImmunoCard STAT!® Mono

30 Tests

750330

USA/Corporate Office:
3471 River Hills Drive
Cincinnati, Ohio 45244
Telephone: 513.271.3700

Orders/Customer Service: 1.800.543.1980
Technical Support: 1.800.343.3858
Information Fax: 513.272.5432
Ordering Fax: 513.271.0124
www.meridianbioscience.com



Belgium / France / Italy / The Netherlands

