

breathID[®] Hp

H. pylori



Your trusted partner in diagnostic accuracy and improving patient outcomes

The state-of-the-art BreathID[®] Hp urea breath test system for in-office *H. pylori* detection is a CLIA non-regulated 15-minute test providing immediate results. The BreathID[®] Hp intuitively combines market-leading diagnostic accuracy with strong reimbursement potential via established CPT codes. Whether for initial active infection detection or confirmation of eradication, Meridian is your trusted partner in gastrointestinal diagnostics.

Market-leading accuracy. Right in your office.

Unparalleled diagnostic performance

- Market-leading non-invasive *H. pylori* diagnostic accuracy
 - 100% Sensitivity
 - 99.2% Specificity
- Clinically recommended for both initial active infection detection and confirmation of eradication post-treatment
- Testing for adult and pediatric (ages 3 — 17) patients using the same test kits
 - No pediatric calculation required
- Strong, steady reimbursements with established CPT codes




Intuitive Efficiency, Compact Footprint

- 15-minute in-office test with immediate point-of-care results
 - Continuous sampling of patient breath through nasal cannula
 - Test and treat in same office visit
- CLIA non-regulated device
 - Easy implementation and training
 - Begin testing patients immediately
 - Self-calibration, no validation necessary
- “Press and walk away” automation
 - User-friendly, one-button touchscreen interface
 - Step-by-step instructions on screen
 - Only 2-minutes hands-on time
 - Built-in automated monitoring manages testing analysis stage
- Data management options
 - Immediate printed results
 - Print or export results
 - EMR compatible



From start to finish in 3 easy steps

10-15 minutes to diagnose *H. pylori*

		
1 Baseline Breath Patient breathes normally via nasal cannula	2 Drink and Analyze Patient drinks ^{13}C -Urea substrate as instrument continuously analyzes breath sample	3 Results Diagnosis is complete in as little as 10 minutes

Reimbursement Advantages

BreathID® Systems Potential Reimbursements

Established CPT Codes	Description	2021 National Median Medicare Rate*
83013	<i>Helicobacter pylori</i> ; breath test analysis for urease activity, non-radioactive isotope (e.g., C-13)	\$ 67.36
83014	<i>Helicobacter pylori</i> ; drug administration and test collection.	\$ 7.86
Total Reimbursements:		\$ 75.22

- Clinical guidelines strongly recommend urea breath testing for detection AND testing again after treatment to confirm eradication. Both tests are reimbursable.
- Serology antibody testing is not recommended and is not reimbursed by most insurance companies.
- Over 30 different ICD-10 codes warrant testing for *H. pylori*.

* Reimbursement information utilize the national median Medicare rates for CPT codes as published in the 2021 Fee Schedule.

Ordering information:

Item	Catalog Number
BreathID® Hp	VS01000
IDkit Hp™ One Kit	AC00051

For more information on the Meridian Bioscience BreathID® Hp System, contact a specialist at **1-888-763-6769**.

WORLD HEADQUARTERS

3471 River Hills Drive

Cincinnati, Ohio 45244

customerservice@meridianbioscience.com



BreathID® Hp Physical Specifications

Height: 16" (406 mm) | **Width:** 15" (381 mm) |
Depth: 22" (559 mm) | **Weight:** 33 lbs (15 kg)

© 2021 Meridian Bioscience Corporation. Meridian and BreathID® are trademarks of Meridian Bioscience Corporation. All rights reserved.

MBI_BreathID_Hp_Platform_Brochure_2021-08_Digital