



# breathID<sup>®</sup> smart

*H. pylori*

---

## Your trusted partner in diagnostic accuracy and compact operational efficiencies

The state-of-the-art BreathID<sup>®</sup> Smart urea breath test system for *H. pylori* detection intuitively combines market-leading diagnostic accuracy with efficient automation, all in a compact footprint for maximizing mid-range testing throughput. Whether for initial active infection detection or confirmation of eradication, Meridian is your trusted partner in gastrointestinal diagnostics.

**meridian** BIOSCIENCE<sup>®</sup>  
LIFE DISCOVERED. LIFE DIAGNOSED.

# Market-leading performance. Compact efficiencies.

## Unparalleled diagnostic performance

---

- Market-leading non-invasive *H. pylori* diagnostic accuracy
  - 100% Sensitivity
  - 97.9% Specificity
- Clinically recommended for both initial active infection detection and confirmation of eradication
- Testing for adult and pediatric (ages 3 — 17) patients using the same test kits
  - No pediatric calculation required
- Strong, steady reimbursements with established CPT codes

## Automated efficiency, compact footprint

---

- Perfect platform for low to mid-volume laboratories and clinics
  - Test samples have a two-week window for analysis and require no special storage
  - 4 test samples per Auto Sampler
    - Under 10 minutes for full batch analysis, 2 minutes per test
- User-friendly interface
  - Touchscreen for intuitive process management
  - Scannable barcodes maximize speed of data entry
  - Automated QC selection and label generation
- Flexible data management options
  - Print, export, or transmit results
  - Secure, seamless integration with LIS
- Minimizes potential errors before they happen
  - Automatic identification and correction of reversed bags
  - Coordinated indicator signals guide batch assembly and bag placement
  - Scannable barcodes minimize human error



## From Scan to Results in 3 Steps

2 minutes to diagnose *H. pylori*



**1** Scan bag barcode



**2** Batch set up



**3** Press Start

# The BreathID® Platform Advantage

## Trusted *H. pylori* partner

- The BreathID® platform can be used in conjunction with other Meridian *H. pylori* platforms for a total testing solution across the healthcare system that offers:
  - Flexibility in testing options to support near-patient testing or centralized testing in a core lab
  - The ability for system administrators to be responsive to changes in testing requirements and demand fluctuations within the health network
  - An unparalleled *H. pylori* portfolio that accommodates any laboratory setting from low volume rapid to high volume automation
  - Physicians a choice in testing for active infection
  - Improved operational efficiencies and the capability to deliver standardization across the network by consolidating testing to one trusted partner

## Ordering information:

Item	Catalog Number
BreathID® Smart	VS07000W
IDkit Hp® Two	AC00063

**For more information** on the Meridian Bioscience BreathID® Smart, contact a specialist at **1-888-763-6769**.

## WORLD HEADQUARTERS

3471 River Hills Drive

Cincinnati, Ohio 45244

[customerservice@meridianbioscience.com](mailto:customerservice@meridianbioscience.com)



## BreathID® Smart Physical Specifications

**Height:** 16 " (406 mm) | **Width:** 15" (381 mm) |

**Depth:** 22" (559 mm) | **Weight:** 33 lbs (15 kg)