

Certificate of Analysis

COA No: CA_XBE-0005

Version: 10

Catalog No: MDX054 Lot No: EN070-B107820 Suitable for Research and further Manufacturing Use Storage Conditions: -20°C Component Lot No: UNG-122105A Expiry date: June2024

Quality Control Parameters

Enzyme that efficiently hydrolyzes uracil from ssDNA or dsDNA. Endonuclease, exonuclease, nickase and RNase-free

Analysis	Specification	Result
Functional	dUTP PCR products are degraded on the addition of UDGase in the presence of 10x UDGase reaction buffer. Single distinct bands were observed with agarose gel electrophoresis (ethidium stained).	Passed
DNA contamination	Quantitative PCR analysis with no template. Presence of <i>E. coli</i> and mouse genomic DNA checked. Test sample must amplify in line with a reference sample.	Passed
DNase contamination	Incubation of a 1Kb double stranded DNA fragment. Incubation for 4 hours at 37°C with dilution series of DNase I. Analysed by agarose gel electrophoresis. Test sample must show less degradation than the limit of detection 2.5 x 10 ⁻³ U DNase.	Passed

QA / QC Representative:

Alberta Newton

Date: 12th May 2022

United Kingdom

Tel: +44 (0)20 8830 5300 Fax: +44 (0)20 8452 2822 <u>USA</u> Tel: +1 901.382.8716 Fax: +1 901.382.0027 <u>Germany</u> Tel: +49 (0)3371 60222 00 Fax: +49 (0)3371 60222 01 <u>Australia</u>

Tel: +61 (0)2 9209 4180 Fax: +61 (0)2 9209 4763

Info@meridianlifescience.com; www.meridianlifescience.com