

Certificate of Analysis

COA No: CA_BSM-0105

Version: 01

Lyo-Ready™ LAMP Mix, 4x

For research or further manufacturing use only

Catalog No:	MDX097
Lot No:	B116280
Storage Conditions:	-20°C
Component Lot No:	123203A
Expiry date:	April 2025

Quality Control Parameters

Analysis	Specification	Result
Functional	TTR (time-to-result) values obtained by Loop mediated amplification reaction of Test Sample vs Reference sample must be within 1.5 minutes	Passed
DNA contamination	Quantitative PCR analysis with no template. Presence of <i>E. coli</i> and mouse genomic DNA checked. Test sample must amplify in concordance with control sample.	Passed
DNase contamination	Incubation of a 1 Kb double stranded DNA fragment. Incubation for 4 hours at 37° C with dilution series of DNase I. Analysed by agarose gel electrophoresis. Test sample must exhibit less degradation than the limit of detection 2.5 x 10 ⁻³ U DNase.	Passed
RNase contamination	Quantitative PCR analysis with high and low RNase standards. Test sample must show less RNase activity than the limit of detection 9.7 x 10^{-3} ng/ μ L RNase.	Passed

QA / QC Representative:



Alberta Newton

Date: 30th March 2023

United Kingdom

Tel: +44 (0)20 8830 5300 Fax: +44 (0)20 8452 2822 <u>USA</u>

Tel: +1 901.382.8716 Fax: +1 901.382.0027 <u>Germany</u>

Tel: +49 (0)3371 60222 00 Fax: +49 (0)3371 60222 01 <u>Australia</u>

Tel: +61 (0)2 9209 4180 Fax: +61 (0)2 9209 4763