

**Lyo-Ready™ Direct RNA/DNA qPCR  
Urine, 4x**

For R&D or further manufacturing use only

Cat No.: MDX153

Storage Conditions: -20 °C

Lot number: RD122H

Expiry Date: July 2024

**Test Parameters**

Analysis	Specification	Result
Functionality	Quantitative PCR analysis amplifying a multiplex reaction of 3 genes from mouse total RNA at different concentrations. Ct must be consistent for the test and reference sample with $\pm$ 0.5 Ct variance.	Passed

Authorised by Thomas Fouqueau  
Team Leader, R&D



Technical Support: [mbi.tech@meridianlifescience.com](mailto:mbi.tech@meridianlifescience.com)

For product handling guidelines please visit [meridianlifescience.com](http://meridianlifescience.com)

**Bioline Reagents Ltd**  
UNITED KINGDOM

**Meridian Life Science Inc**  
USA

**Bioline GmbH**  
GERMANY

**Bioline (Aus) Pty Ltd**  
AUSTRALIA

**Bioline France**  
FRANCE

Tel: +44 (0)20 8830 5300  
Fax: +44 (0)20 8452 2822

Tel: +1 901 382 8716  
Fax: +1 901 382 0027

Tel: +49 (0)337 168 1229  
Fax: +49 (0)337 168 1244

Tel: +61 (0)2 9209 4180  
Fax: +61 (0)2 9209 4763

Tel: +33 (0)1 42 56 04 40  
Fax: +33 (0)9 70 06 62 10