

**Glycerol-Free SP6 RNA  
Polymerase (HC), 500 U/μL**

Cat No.: MRNA106

Storage Conditions: -20 °C

Lot number: RD123K

**Test Parameters**

Analysis	Specification	Result
Unit Activity	<p>Activity is measured using a real-time <i>in vitro</i> transcription assay against a standard curve of known enzyme concentration.</p> <p><u>Pass Criteria:</u> Activity must be &gt; 500 U/μL</p>	545 U/μL

Authorised by Thomas Fouqueau  
R&D Team Leader

