

statID PRO™

COVID-19/Flu A&B



StatID PRO™ COVID-19/Flu A&B

Fast, reliable answers when every minute matters

StatID PRO COVID-19/Flu A&B combo rapid test delivers quick, reliable results for COVID-19, influenza A, and influenza B from a simple anterior nasal swab. FDA-cleared for use in patients 2 years and older, the test is designed for use within the first five days of symptoms. This easy-to-use rapid assay helps healthcare providers quickly differentiate between common respiratory infections with similar presentations.

meridian BIOSCIENCE®

StatID PRO COVID-19/Flu A&B

The Challenge

- COVID-19 and influenza present with overlapping symptoms — guessing delays care
- Seasonal surges overwhelm hospitals and clinics, slowing treatment decisions
- Running multiple tests means more swabs, more time, and higher costs
- Missed treatment windows reduce effectiveness of antivirals and isolation efforts

The Need

- A single swab test that detects COVID-19, Flu A, and Flu B at the point of care
- Rapid results in minutes
- Clear, easy-to-read answers with minimal steps and low error risk
- Simple, consistent workflow

StatID PRO Delivers

- One swab, three answers — COVID-19, Flu A, and Flu B antigen results
- FDA-cleared, not EUA — trusted regulatory approval with proven reliability
- Approved for use in patients 2 years and older
- 15-minute visual results for immediate decision-making
- Workflow flexibility - specimens may be held in extraction tube for up to 1 hour before testing

Specifications

CLIA Status

Waived

Turnaround Time

15 minutes

Sample Type

Anterior Nasal Swab

Sample Storage

Samples should be tested as soon as possible after collection. Based on data generated with the StatID PRO COVID-19/Flu A&B, nasal swabs are stable for up to 1 hour after sample collection

Kit Storage

36-86 F (2-30 C)

Performance

COVID-19	Flu A	Flu B
PPA: 92%	PPA: 92.5%	PPA: 90.5%
NPA: 99%	NPA: 99.9%	NPA: 99.9%

Catalog Number:

StatID PRO COVID-19 / Flu A&B 783025

StatID PRO COVID-19 / Flu A&B Controls 783010

CPT Code

87428-QW



For more information on the StatID PRO products, contact a specialist at meridianbioscience.com/contactus