

Cyfra-21-1 Antigen

Cat. No. A32340H

STORE AT -30 to -10°C

DESCRIPTION	Cyfra-21-1, Cytokeratin 19 Fragment, Low Cross-Reactivity (LCR) Calibrator Grade
TYPE	Antigen
SHIPPING & STORAGE INFORMATION	Store at -30°C to -10°C. Avoid multiple freeze/thaw cycles.
SAFETY PRECAUTIONS	Refer to the appropriate Safety Data Sheet (SDS) for additional information.
NOTES	This reagent has been manufactured under ISO13485 quality management system and is suitable for research or further manufacturing use only. Centrifuge before opening to ensure complete recovery of vial contents.

General Characteristics

SOURCE	MCF-7 cell supernatant
FORMAT	Calibrator, liquid
MOLECULAR WEIGHT	40 kDa
KNOWN CONTAMINANTS	As reported
PURIFICATION	Ion-exchange chromatography Product is 0.2 µm filtered.
CONCENTRATION	Cytokeratin 19 Ag: as reported Total protein: as reported (OD280)
INACTIVATION	Not applicable
BUFFER	PBS, pH 7.4 ± 0.1
PRESERVATIVES	None

Applications

Suitable for use in CLIA. Each laboratory should determine an optimum working titer for its particular assay format. Other assay formats have not been tested but should not necessarily be excluded.