

Bovine Gamma Globulin

Cat. No. A51300B

STORE AT 2 to 8°C

DESCRIPTION	Bovine Gamma Globulin
SHIPPING & STORAGE INFORMATION	Store at 2 to 8°C.
SAFETY PRECAUTIONS	Refer to the appropriate Safety Data Sheet (SDS) for additional information.
NOTES	This reagent has been manufactured under ISO13485 quality management system and is suitable for research or further manufacturing use only. Centrifuge before opening to ensure complete recovery of vial contents.

» For ordering and to learn about pack sizes, please [click here](#).

General Features/Characteristics

SOURCE	Normal bovine serum or plasma collected in a USDA certified facility within the contiguous 48 states of the United States.
FORMAT	Purified, Lyophilized. Off white, fine, flaky powder.
BUFFER	Not applicable
PRESERVATIVES	None
PURIFICATION	Purified using modified Cohn fractionation. Contains > 99% gamma globulin, which is less than 5% impurities as determined by serum protein electrophoresis (percent of the total protein present).
CONCENTRATION	Not applicable
RECONSTITUTION	Reconstitute in at least 0.2% Sodium Chloride.
INACTIVATION	Not applicable

Applications

Suitable for use in Immunoelectrophoresis (10 mg/mL). Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.