

Rabbit IgG - Differential PPT

Cat. No. A66105R

STORE AT -30 to -10°C

DESCRIPTION	Purified Rabbit Immunoglobulin G (IgG), Differential Precipitation
TYPE	Blocker
SHIPPING & STORAGE INFORMATION	Aliquot and store at -30 to -10°C.
SAFETY PRECAUTIONS	Refer to the appropriate Safety Data Sheet (SDS) for additional information.
NOTES	This reagent has been manufactured under ISO13485 quality management system and is suitable for research or further manufacturing use only. Centrifuge before opening to ensure complete recovery of vial contents.

» For ordering and to learn about pack sizes, please [click here](#).

General Features/Characteristics

HOST ANIMAL	Rabbit
COUNTRY OF ORIGIN	USA. This product was collected from normal healthy rabbits under veterinarian supervision.
FORMAT	Purified, Liquid
BUFFER	20 mM Potassium Phosphate Monobasic, 5 mM Sodium Phosphate Dibasic, 150 mM Sodium Chloride, pH 5.9 - 6.1
PRESERVATIVES	0.1% Sodium Azide
PURIFICATION	Differential Precipitation, plus ion-exchange chromatography if needed. Product is 0.2 µm filtered.
PURITY	SDS-PAGE (reduced and non-reduced) ≥ 90%
CONCENTRATION	OD280nm ($E^{0.1\%} = 1.4$) Acceptable Range: 50 - 60 mg/mL
BIOBURDEN	Direct Inoculation Acceptable Range: ≤ 5 CFU/mL

Applications

Specific methodologies have not been tested using this product.