

## TRU Block™

Cat. No. A66800H

STORE AT -30 to -10°C

|                                           |                                                                                                                                                                                                                     |
|-------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>DESCRIPTION</b>                        | TRU Block™ is a complex blocker that contains active ingredients and purified mouse IgG.                                                                                                                            |
| <b>TYPE</b>                               | Universal (Active & Passive) Blocker                                                                                                                                                                                |
| <b>SHIPPING &amp; STORAGE INFORMATION</b> | Store at -30 to -10°C.                                                                                                                                                                                              |
| <b>SAFETY PRECAUTIONS</b>                 | Refer to the appropriate Safety Data Sheet (SDS) for additional information.                                                                                                                                        |
| <b>NOTES</b>                              | This reagent has been manufactured under ISO13485 quality management system and is suitable for research or further manufacturing use only. Centrifuge before opening to ensure complete recovery of vial contents. |

» For ordering and to learn about pack sizes, please [click here](#).

### General Features/Characteristics

|                      |                                                                         |
|----------------------|-------------------------------------------------------------------------|
| <b>SOURCE</b>        | Mouse Serum<br>Normal mouse serum plus proprietary ingredients          |
| <b>FORMAT</b>        | Purified, Liquid                                                        |
| <b>BUFFER</b>        | 10 mM Sodium Phosphate, 150 mM Sodium Chloride. pH 7.0 - 7.4            |
| <b>PRESERVATIVES</b> | 0.05% Sodium Azide                                                      |
| <b>PURIFICATION</b>  | Protein A purification                                                  |
| <b>PURITY</b>        | SDS-PAGE<br>Acceptable Range: ≥ 95%                                     |
| <b>CONCENTRATION</b> | OD 280nm Assay ( $E^{0.1\%} = 1.4$ )<br>Acceptable Range: 24 – 26 mg/mL |

### Applications

Active blocking of heterophilic antibody interference in immunoassays.