

Human CA125 Calibrator Grade

Cat. No. A86125H

STORE AT -30 to -10°C

DESCRIPTION	Human Cancer Antigen 125 (CA125), (Ovarian Cancer) Calibrator Grade
TYPE	Antigen
SHIPPING & STORAGE INFORMATION	Store at -30 to -10°C. Avoid multiple freeze/thaw cycles.
SAFETY PRECAUTIONS	Refer to the appropriate Safety Data Sheet (SDS) for additional information.
NOTES	This reagent has been manufactured under ISO13485 quality management system and is suitable for research or further manufacturing use only. Centrifuge before opening to ensure complete recovery of vial contents.

» For ordering and to learn about pack sizes, please [click here](#).

General Features/Characteristics

SOURCE	Human Adenocarcinoma
FORMAT	Calibrator Grade, Liquid
BUFFER	Phosphate Buffered Saline, pH 7.4
PRESERVATIVES	0.05% Sodium Azide
INACTIVATION	Not Applicable

Applications

Calibrator in Immunoassays. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.