

Rabbit anti C. botulinum Neurotoxin Type A

Cat. No. B76750R

STORE AT -30 to -10°C

DESCRIPTION	Rabbit Antibody to Clostridium botulinum Neurotoxin Type A Complex
TYPE	Polyclonal Antibody
SHIPPING & STORAGE INFORMATION	Store at -30 to -10°C. Aliquot to avoid multiple freeze/thaw cycles
SAFETY PRECAUTIONS	Refer to the appropriate Safety Data Sheet (SDS) for additional information.
NOTES	This reagent has been manufactured under ISO13485 quality management system and is suitable for research or further manufacturing use only. Centrifuge before opening to ensure complete recovery of vial contents.

» For ordering and to learn about pack sizes, please [click here](#).

General Features/Characteristics

SPECIFICITY	C. botulinum Type A Neurotoxin Complex, neurotoxin or protein isolates
HOST ANIMAL	Rabbit
FORMAT	Purified, Liquid
IMMUNOGEN	Type A Neurotoxin and Complexing Proteins Toxoid
BUFFER	1X Phosphate Buffered Saline, pH 7.2 containing 20% Glycerol
PRESERVATIVES	None
PURIFICATION	Protein G Sepharose Chromatography
CONCENTRATION	OD 280nm Assay ($E^{0.1\%} = 1.4$)

Applications

Suitable for use in Dot Blot, ELISA (1:10,000) and Western Blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.