

Streptococcus pneumoniae Native Antigen

Cat. No. BN1188

STORE AT -20°C

DESCRIPTION*Streptococcus pneumoniae* Native Antigen**SHIPPING & STORAGE INFORMATION**

Store at -20°C. Avoid multiple freeze / thaw cycles.

SAFETY PRECAUTIONS

Refer to the appropriate Safety Data Sheet (SDS) for additional information.

NOTES

This reagent has been manufactured under ISO13485 quality management system and is suitable for research or further manufacturing use only. Centrifuge before opening to ensure complete recovery of vial contents.

» For ordering and to learn about pack sizes, please [click here](#).

General Features/Characteristics

SOURCE

Streptococcus pneumoniae (serotype 19A)

FORMAT

Liquid

PURIFICATION

None (crude lysate)

CONCENTRATION

CFU/mL as reported (measure of cell concentration in the lysate)

INACTIVATION

Ultrasonication. Confirmed by absence of colonies after 72 hours.

BUFFER

PBS

PRESERVATIVES

0.09% Sodium Azide

Applications

Suitable for use in Lateral Flow and ELISA. Each laboratory should determine an optimum working titer for its particular assay format. Other assay formats have not been tested but should not necessarily be excluded.