

Rubeola Nucleoprotein Antigen, Recombinant

Cat. No. EV9298

STORE AT $\leq -65^{\circ}\text{C}$

DESCRIPTION	Rubeola Nucleoprotein Antigen, Recombinant
SHIPPING & STORAGE INFORMATION	Store at $\leq -65^{\circ}\text{C}$. Avoid multiple freeze / thaw cycles.
SAFETY PRECAUTIONS	Refer to the appropriate Safety Data Sheet (SDS) for additional information.
NOTES	This reagent has been manufactured under ISO13485 quality management system and is suitable for research or further manufacturing use only. Centrifuge before opening to ensure complete recovery of vial contents.

» For ordering and to learn about pack sizes, please [click here](#).

General Features/Characteristics

SOURCE	Grown in <i>sf-9</i> cells. Product manufactured in the USA using only domestic animal materials processed in U.S. facilities.
VIRAL STRAIN	Edmonston
FORMAT	Purified, liquid
IDENTITY	SDS-PAGE, Western blot
CONCENTRATION	0.4 – 0.6 mg/mL (BCA assay)
BUFFER	50 mM HEPES, pH 7.4 – 7.6
PRESERVATIVES	None

Applications

Each laboratory should determine an optimum working titer for its particular assay format.