

Goat anti Human IgM

Cat. No. L66607G

STORE AT 2 to 8°C

DESCRIPTION	Goat Antibody to Human Immunoglobulin M (IgM). Turbidimetric Grade
TYPE	Blocker
SHIPPING & STORAGE INFORMATION	Store at 2 to 8°C. Slight precipitation may occur on storage which can be removed by centrifugation and should not affect performance characteristics.
SAFETY PRECAUTIONS	Refer to the appropriate Safety Data Sheet (SDS) for additional information.
NOTES	This reagent has been manufactured under ISO13485 quality management system and is suitable for research or further manufacturing use only. Centrifuge before opening to ensure complete recovery of vial contents.

» For ordering and to learn about pack sizes, please [click here](#).

General Features/Characteristics

SPECIFICITY	Monospecific by Immunoelectrophoresis (IEP) when tested against pooled human plasma and 2x concentrated pooled human serum.
HOST ANIMAL	Goat
IMMUNOGEN	Purified Human IgM
FORMAT	Monospecific Antiserum. Liquid
PURIFICATION	Delipidation followed by solid phase absorption. Product is 0.2 µm filtered.
BUFFER	0.05 M Tris, 0.5 M Sodium Chloride, pH 7.5.
PRESERVATIVES	0.05% Sodium Azide
CONCENTRATION	OD 280nm Units/ml

Applications

Specific methodologies have not been tested using this product. Each laboratory should determine an optimum working titer for use in its particular application.