

Chikungunya E1 Wild-Type Recomb.

Cat. No. R01653

STORE AT -30 to -10°C

DESCRIPTION	Chikungunya E1 Wild-Type Recombinant Antigen. Product contains histidinetag.
MOLECULAR WEIGHT	~50 kDa
TYPE	Recombinant Antigen
SHIPPING & STORAGE INFORMATION	Aliquot and store at -30 to -10°C.
SAFETY PRECAUTIONS	Refer to the appropriate Safety Data Sheet (SDS) for additional information.
NOTES	This reagent has been manufactured under ISO13485 quality management system and is suitable for research or further manufacturing use only. Centrifuge before opening to ensure complete recovery of vial contents.

» For ordering and to learn about pack sizes, please [click here](#).

General Features/Characteristics

SOURCE	Insect Cells
FORMAT	Purified
BUFFER	Phosphate Buffered Saline, pH 7.4
PRESERVATIVES	0.09% Sodium Azide
PURIFICATION	Affinity Chromatography
CONCENTRATION	Bradford Method
PURITY	SDS-PAGE > 85%

Applications

Suitable for use as a control and in ELISA, Lateral Flow application. Other applications have not been tested but use in such assays should not necessarily be excluded.