

Dengue Type 1 EP Recombinant

Cat. No. R01659

STORE AT -30 to -10°C

DESCRIPTION	Dengue Type 1 Envelope Protein (EP) (a.a. 281-675) Recombinant. Product contains histidine-tag.
MOLECULAR WEIGHT	~50 kDa
SHIPPING & STORAGE INFORMATION	Aliquot and store at -30 to -10°C.
SAFETY PRECAUTIONS	Refer to the appropriate Safety Data Sheet (SDS) for additional information.
NOTES	This reagent has been manufactured under ISO13485 quality management system and is suitable for research or further manufacturing use only. Centrifuge before opening to ensure complete recovery of vial contents.

» For ordering and to learn about pack sizes, please [click here](#).

General Features/Characteristics

SOURCE	Drosophila
FORMAT	Purified, Liquid
BUFFER	Phosphate Buffered Saline, pH 7.4
PRESERVATIVES	0.09% Sodium Azide
PURIFICATION	Affinity Chromatography
CONCENTRATION	Bradford Method
PURITY	SDS-PAGE > 85%
INACTIVATION	Not Applicable

Applications

Suitable for use in ELISA and Lateral Flow. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.